



碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301或800-8283301
 订货e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

细胞凋亡-DNA Ladder抽提试剂盒

产品编号	产品名称	包装
C0007	细胞凋亡-DNA Ladder抽提试剂盒	50次

产品简介:

- 碧云天生产的细胞凋亡-DNA Ladder抽提试剂盒, 是针对细胞凋亡过程中产生的核小体间DNA链断裂而设计的。可以非常有效地抽提最小片段为180-200bp的DNA ladder, 同时又可以抽提到50kb以上的基因组DNA。
- DNA ladder也称DNA fragmentation, 是细胞凋亡的一个重要指标。通常观察到DNA ladder, 就可以判定细胞发生了凋亡。
- 本试剂盒足够抽提50个细胞或组织样品。

包装清单:

产品编号	产品名称	包装
C0007-1	样品裂解液	30ml
C0007-2	蛋白酶K	130μl
C0007-3	10M 醋酸铵	6ml
C0007-4	TE	6ml
—	说明书	1份

保存条件:

-20°C保存, 一年有效。10M 醋酸铵和TE也可以室温保存。

注意事项:

- 需自备Tris平衡苯酚、氯仿和无水乙醇。
- 本产品仅限于专业人员的科学研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

使用说明:

1. 样品收集

- 对于组织样品:**
切下组织, 并剪切成小块, 置液氮中冻结, 研碎或捣碎。或直接冰浴上匀浆。
- 对于贴壁细胞:**
胰酶消化后, PBS或生理盐水洗一次, 1000-2000g离心1-2分钟, 弃上清, 收集细胞。
- 对于悬浮细胞:**
1000-2000g离心1-2分钟, 弃上清, 收集细胞。

2. DNA ladder抽提

- 每1毫升样品裂解液中加入5微升蛋白酶K, 混匀。
- 对于上述收集好的样品, 每5毫克组织或者 10^6 个细胞中加入500微升添加了蛋白酶K的样品裂解液, Vortex混匀, 充分裂解组织或细胞。
- 50°C水浴消化过夜(通常12-20小时皆可)。
- 加入500微升Tris平衡苯酚(pH8.0)。
- Vortex剧烈混匀, 使有机相和水相充分混合, 以达到抽提效果。4°C, 12,000g离心5分钟。
- 缓慢吸出上层水相至另一洁净离心管中。注意勿触及中间层, 中间层通常含有变性的蛋白等, 并用等体积Tris平衡酚再抽提一次(同步步骤e)。
- 缓慢吸出上层水相至另一洁净离心管中。注意勿触及中间层, 中间层通常含有变性的蛋白等, 并用等体积氯仿再抽提一次(同步步骤e)。
- 慢慢吸出约300微升上清液, 加入60微升10M醋酸铵和600微升无水乙醇, 颠倒数次混匀, 此时可见DNA沉淀产生。-20°C冻存1小时, 以充分沉淀小片段DNA。冻存过夜或-70°C冻存效果更佳。
- 4°C, 12,000g离心10分钟, 弃上清。
- 加入600微升70%乙醇, 轻轻颠倒约2次。4°C, 12,000g离心10分钟, 小心吸去上清。注意: 70%乙醇洗涤的时候, 千万注意避免损失一些细小的DNA沉淀, 这些沉淀中大部分是你所需的DNA ladder。

- k. 尽量吸除残余的乙醇，待看不到明显的液体时，立即加入50-100微升TE溶解DNA。注意：不可过分干燥基因组DNA沉淀，否则会极难溶解。如果发现DNA沉淀难以溶解，可以在4°C用摇床缓慢摇动过夜，以溶解DNA沉淀。
1. 取部分抽提得到的DNA，1%琼脂糖凝胶电泳分析。如果细胞发生凋亡，就可以观察到典型的DNA ladder。电泳时一定要注意换用新鲜配制的电泳液，DNA凝胶也要用新鲜配制的电泳液配制并新鲜配制后使用。电泳时为获取最佳电泳效果使ladder充分分开，电泳速度宜适当慢一些，凝胶宜适当长一些，而加样孔宜更加扁平一些。选取适当较薄的梳齿，往往会获得更好的ladder电泳效果。

使用本产品的文献：

1. Wang N, Hou L, Zhu ZY, Zhang JY, Ju CX, Ge YL, Yue W. Effect of 2, 3-dioxo indoline on apoptosis and human telomerase reverse transcriptase mRNA in SH2SY5Y cells. *Chinese Pharmacological Bulletin*. 2007 May;23(5):680-4.
2. Cai L, Wang H, Li Q, Qian Y, Yao W. Salidroside inhibits H₂O₂-induced apoptosis in PC12 cells by preventing cytochrome c release and inactivating of caspase cascade. *ACTA BIOCHIM BIOPH SIN*. 2008;40(9):796-802.
3. Qian YF, Wang H, Yao WB, Gao XD. Aqueous extract of the Chinese medicine, Danggui-Shaoyao-San, inhibits apoptosis in hydrogen peroxide-induced PC12 cells by preventing cytochrome c release and inactivating of caspase cascade. *Cell Biol Int*. 2008 Feb;32(2):304-11.
4. Bai Y, Li Q, Yang J, Zhou X, Yin X, Zhao D. p75(NTR) activation of NF- κ B is involved in PrP106-126-induced apoptosis in mouse neuroblastoma cells. *Neurosci Res*. 2008 Sep;62(1):9-14.
5. Wang H, Xu Y, Yan J, Zhao X, Sun X, Zhang Y, Guo J, Zhu C. Acteoside protects human neuroblastoma SH-SY5Y cells against beta-amyloid-induced cell injury. *Brain Res*. 2009 Aug 4;1283:139-47.
6. Chai J, Xiong Q, Zhang P, Zheng R, Peng J, Jiang S. Induction of Ca²⁺-signal mediated apoptosis and alteration of IP3R1 and SERCA1 expression levels by stress hormone in differentiating C2C12 myoblasts. *GEN COMP ENDOCR*. 2010;166(2):241-9.
7. Ou Y, Dong X, Liu XY, Cheng XC, Cheng YN, Yu LG, Guo XL. Mechanism of tetramethylpyrazine analogue CXCI95 inhibition of hydrogen peroxide-induced apoptosis in human endothelial cells. *Biol Pharm Bull*. 2010;33(3):432-8.
8. Cao H, Feng Q, Xu W, Li X, Kang Z, Ren Y, Du L. Genipin induced apoptosis associated with activation of the c-Jun NH2-terminal kinase and p53 protein in HeLa cells. *Biol Pharm Bull*. 2010;33(8):1343-8.
9. Li J, Xu Z, Tan M, Su W, Gong XG. 3-(4-(Benzo[d]thiazol-2-yl)-1-phenyl-1H-pyrazol-3-yl) phenyl acetate induced Hep G2 cell apoptosis through a ROS-mediated pathway. *CHEM-BIOL INTERACT*. 2010 Feb 12;183(3):341-8.
10. Li Y, Niu Y, Wu H, Sun Y, Li Q, Kong X, Liu L, Mei Q. Modified apple polysaccharides could induce apoptosis in colorectal cancer cells. *J Food Sci*. 2010 Oct;75(8):H224-9.
11. Wang C, Cao J, Qu J, Li Y, Peng B, Gu Y, He Z. Recombinant vascular basement membrane derived multifunctional peptide blocks endothelial cell angiogenesis and neovascularization. *J Cell Biochem*. 2010 Oct 1;111(2):453-60.
12. Gao LW, Zhang J, Yang WH, Wang B, Wang JW. Glaucocalyxin A induces apoptosis in human leukemia HL-60 cells through mitochondria-mediated death pathway. *Toxicol In Vitro*. 2011 Feb;25(1):51-63.
13. Feng Q, Cao HL, Xu W, Li XR, Ren YQ, Du LF. Apoptosis induced by genipin in human leukemia K562 cells: involvement of c-Jun N-terminal kinase in G2/M arrest. *Acta Pharmacol Sin*. 2011 Apr;32(4):519-27.
14. Huang H, Zhao X, Chen Y, Lu R, Wu Y. Apoptosis induced by ZnPcH1-based photodynamic therapy in Jurkat cells and HEL cells. *Int J Hematol*. 2011 Dec;94(6):539-44.
15. Wu G, Qian Z, Guo J, Hu D, Bao J, Xie J, Xu W, Lu J, Chen X, Wang Y. Ganoderma lucidum Extract Induces G1 Cell Cycle Arrest, and Apoptosis in Human Breast Cancer Cells. *AM J CHINESE MED*. 2012;40(3):631-42.
16. Chen HB, Zhou LZ, Mei L, Shi XJ, Wang XS, Li QL, Huang L. Gambogic acid-induced time- and dose-dependent growth inhibition and apoptosis involving Akt pathway inactivation in U251 glioblastoma cells. *J NAT MED-TOKYO*. 2012 Jan;66(1):62-9.
17. Shi ZY, Li YQ, Kang YH, Hu GQ, Huang-fu CS, Deng JB, Liu B. Piperonal ciprofloxacin hydrazone induces growth arrest and apoptosis of human hepatocarcinoma SMMC-7721 cells. *Acta Pharmacol Sin*. 2012 Feb;33(2):271-8.
18. Wu GS, Lu JJ, Guo JJ, Li YB, Tan W, Dang YY, Zhong ZF, Xu ZT, Chen XP, Wang YT. Ganoderic acid DM, a natural triterpenoid, induces DNA damage, G1 cell cycle arrest and apoptosis in human breast cancer cells. *Fitoterapia*. 2012 Mar;83(2):408-14.
19. Ye S, Jiang Y, Zhang H, Wang Y, Wu Y, Hou Z, Zhang Q. Multi-walled carbon nanotubes induce apoptosis in RAW 264.7 cell-derived osteoclasts through mitochondria-mediated death pathway. *J NANOSCI NANOTECHNO*. 2012 Mar;12(3):2101-12.
20. Tao LL, Cheng YY, Ding D, Mei S, Xu JW, Yu J, Ou-Yang Q, Deng L, Chen Q, Li QQ, Xu ZD, Liu XP. C/EBP- α ameliorates CCl₄-induced liver fibrosis in mice through promoting apoptosis of hepatic stellate cells with little apoptotic effect on hepatocytes in Apoptosis. 2012 May;17(5):492-502.
21. Gong K, Xie J, Yi H, Li W. CS055 (Chidamide/HBI-8000), a novel histone deacetylase inhibitor, induces G1 arrest, ROS-dependent apoptosis and differentiation in human leukaemia cells. *Biochem J*. 2012 May 1;443(3):735-46.
22. Jin M, Wang J, Xiao R, Liu X, Wu F, Pang Y, Feng B, Yang D, Li Q. Effects of the recombinant toxin protein rLj-RGD3 in multidrug-resistant human breast carcinoma cells. *ACTA BIOCHIM BIOPH SIN*. 2012 May;44(5):455-61.
23. Yuan Z, Long C, Junming T, Qihuan L, Youshun Z, Chan Z. Quercetin-induced apoptosis of HL-60 cells by reducing PI3K/Akt. *Mol Biol Rep*. 2012 Jul;39(7):7785-93.
24. Zhao P, Han T, Guo JJ, Zhu SL, Wang J, Ao F, Jing MZ, She YL, Wu ZH, Ye LB. HCV NS4B induces apoptosis through the mitochondrial death pathway. *Virus Res*. 2012 Oct;169(1):1-7.
25. Qiao Y, Xiang Q, Yuan L, Xu L, Liu Z, Liu X. Herbacetin induces apoptosis in HepG2 cells: Involvements of ROS and PI3K/Akt pathway. *Food Chem Toxicol*. 2013 Jan;51:426-33.
26. Li R, Cai L, Xie XF, Peng L, Wu TN, Li J. 7,3'-dimethoxy hesperetin inhibits inflammation by inducing synovial apoptosis in rats with adjuvant-induced arthritis. *IMMUNOPHARM IMMUNOT*. 2013 Feb;35(1):139-46.
27. Jin M, Xiao R, Wang J, Liu X, Liu Y, Xue Z, Lv L, Zheng Y, Li Q. Low concentrations of the recombinant toxin protein rLj-RGD3 suppress TNF- α -induced human renal carcinoma cell invasion. *ACTA BIOCHIM BIOPH SIN*. 2013 May;45(5):377-82.
28. Jiang Y, Zhang H, Wang Y, Chen M, Ye S, Hou Z, Ren L. Modulation of apoptotic pathways of macrophages by surface-functionalized multi-walled carbon nanotubes. *PLoS One*. 2013 Jun 6;8(6):e65756.
29. You P, Xing F, Huo J, Wang B, Di J, Zeng S, Liu J. In vitro and in vivo evaluation of anisomycin against Ehrlich ascites carcinoma. *Oncol Rep*. 2013 Jun;29(6):2227-36.
30. Han L, Tang C, Yin C. Effect of binding affinity for siRNA on the in vivo antitumor efficacy of polyplexes. *Biomaterials*. 2013 Jul;34(21):5317-27.
31. Sun JP, Shi ZY, Liu SM, Kang YH, Hu GQ, Huangfu CS, Deng JB, Liu B. Trimethoxy-benzaldehyde levofloxacin hydrazone inducing the growth arrest and apoptosis of human hepatocarcinoma cells. *Cancer Cell Int*. 2013 Jul 2;13(1):67.
32. Song Y, Li X, Li Y, Li N, Shi X, Ding H, Zhang Y, Li X, Liu G, Wang Z. Non-esterified fatty acids activate the ROS-p38-p53/Nrf2 signaling pathway to induce bovine hepatocyte apoptosis in vitro. *Apoptosis*. 2014 Jun;19(6):984-97.
33. Hu H, Tan C, Liu X, Luo F, Li K. Upregulation of the MCL1S protein variant following dihydroartemisinin treatment induces apoptosis in cholangiocarcinoma cells. *Oncology Letters*. 2015;10(6):3545-50.
34. Zhang N, Wang W, Li W, Liu C, Chen Y, Yang Q, Wang Y, Sun K. Inhibition of β -HSD2 expression by triclosan via induction of apoptosis in human placental syncytiotrophoblasts. *J CLIN ENDOCR METAB*. 2015 Apr;100(4):E542-9.
35. Zhang N, Wang WS, Li WJ, Liu C, Wang Y, Sun K. Reduction of progesterone, estradiol and hCG secretion by perfluorooctane sulfonate via induction of apoptosis in human placental

- syncytiotrophoblasts. *Placenta* . 2015 May;36(5):575-80
36. Li XM, Luo XG, Si CL, Wang N, Zhou H, He JF, Zhang TC. Antibacterial active compounds from *Hypericum ascyron* L. induce bacterial cell death through apoptosis pathway. *Eur J Med Chem* . 2015 May 26;96:436-44
 37. Chen P, Hu YF, Wang L, Xiao WF, Bao XY, Pan C, Yi HS, Chen XY, Pan MH, Lu C. Mitochondrial Apoptotic Pathway Is Activated by H₂O₂-Mediated Oxidative Stress in BmN-SWU1 Cells from Bombyx mori Ovary. *PLoS One* . 2015 Jul 30;10(7):e0134694
 38. Fan Y, Patima A, Chen Y, Zeng F, He W, Luo L, Jie Y, Zhu Y, Zhang L, Lei J, Xie X, Zhang H. Cytotoxic effects of β -carboline alkaloids on human gastric cancer SGC-7901 cells. *Int J Clin Exp Med* . 2015 Aug 15;8(8):12977-82
 39. Lu C, Liu X, Liu C, Wang J, Li C, Liu Q, Li Y, Li S, Sun S, Yan J, Shao J. Chlorpyrifos Induces MLL Translocations Through Caspase 3-Dependent Genomic Instability and Topoisomerase II Inhibition in Human Fetal Liver Hematopoietic Stem Cells. *Toxicol Sci* . 2015 Oct;147(2):588-606
 40. Yuan H, Wang J, Wang F, Zhang N, Li Q, Xie F, Chen T, Zhai R, Wang F, Guo Y, Ni W, Tai G. Mucin 1 gene silencing inhibits the growth of SMMC-7721 human hepatoma cells through Bax-mediated mitochondrial and caspase-8-mediated death receptor apoptotic pathways. *Mol Med Rep* . 2015 Nov;12(5):6782-8
 41. Shaochang Yao, Wenjing Huang, Chunliu Pan, Jie Zhan. Caspase-like proteases regulate aluminum-induced programmed cell death in peanut. *PLANT CELL TISS ORG* . 2016 Dec;127(3):691-703
 42. He L, Wang H, Luo L, Li Y, Huang R, Liao L, Zhu Z, Wang Y. Bid-deficient fish delay grass carp reovirus (GCRV) replication and attenuate GCRV-triggered apoptosis. *ONCOTARGET* . 2017 Jul 22;8(44):76408-76422
 43. Yang W, Zhao J, Wang Y, Xu H, Wu Z, Hu Y, Jiang K, Shen P, Ma C, Guan Z, Zhang Y, Ma J, Shang N, Yan G, Wang Z, Dai G. In vivo inhibitory activity of andrographolide derivative ADN-9 against liver cancer and its mechanisms involved in inhibition of tumor angiogenesis. *TOXICOL APPL PHARM* . 2017 Jul 15;327:1-12
 44. Shrestha S, Zhu J, Wang Q, Du X, Liu F, Jiang J, Song J, Xing J, Sun D, Hou Q, Peng Y, Zhao J, Sun X, Song X. Melatonin potentiates the antitumor effect of curcumin by inhibiting IKK β /NF- κ B/COX-2 signaling pathway. *Int J Oncol* . 2017 Oct;51(4):1249-1260
 45. Zhou C, Ma J, Su M, Shao D, Zhao J, Zhao T, Song Z, Meng Y, Jiao P. Down-regulation of STAT3 induces the apoptosis and G1 cell cycle arrest in esophageal carcinoma ECA109 cells. *Cancer Cell Int* . 2018 Apr 4;18:53.
 46. Shao J, Wang C, Li L, Liang H, Dai J, Ling X, Tang H. Luteoloside Inhibits Proliferation and Promotes Intrinsic and Extrinsic Pathway-Mediated Apoptosis Involving MAPK and mTOR Signaling Pathways in Human Cervical Cancer Cells. *Int J Mol Sci* . 2018 Jun 5;19(6): pii: E1664
 47. Que W, Liu H, Yang Q, Xu S. Antiproliferation activities of NK4 on multiple myeloma. *Exp Ther Med* . 2018 Oct;16(4):3668-3673
 48. Ma Y, Li B, Ke Y, Zhang Y. Effects of low doses Trichlorfon exposure on *Rana chensinensis* tadpoles. *Environ Toxicol* . 2019 Jan;34(1):30-36
 49. Lin F, Yang J, Muhammad U, Sun J, Huang Z, Li W, Lv F, Lu Z. Bacillomycin D-C16 triggers apoptosis of gastric cancer cells through the PI3K/Akt and FoxO3a signaling pathways. *ANTI-CANCER DRUG* . 2019 Jan;30(1):46-55
 50. He D, Wang N, Sai X, Li X, Xu Y. Camellia euphlebia protects against corticosterone-induced apoptosis in differentiated PC12 cells by regulating the mitochondrial apoptotic pathway and PKA/CREB/BDNF signaling pathway. *Food Chem Toxicol* . 2019 Apr;126:211-222
 51. Deng Z, Chai J, Zeng Q, Zhang B, Ye T, Chen X, Xu X. The anticancer properties and mechanism of action of tablysin-15, the RGD-containing disintegrin, in breast cancer cells. *Int J Biol Macromol* . 2019 May 15;129:1155-1167
 52. Zhang J, Ahmad S, Wang LY, Han Q, Zhang JC, Luo Y. Cell death induced by α -terthienyl via reactive oxygen species-mediated mitochondrial dysfunction and oxidative stress in the midgut of *Aedes aegypti* larvae. *FREE RADICAL BIO MED* . 2019 Jun;137:87-98
 53. Zhan H, Zhao H, Muhammad N, Li T, Liu Y, Wang J. Lytic peptide-grafted beta-cyclodextrin polymer based nano-scaled drug delivery system with enhanced camptothecin anti-cancer efficacy. *Nanotechnology* . 2020 Feb 7;31(7):075101
 54. Zhenzhou Xiao, Yan Chen, Zhaolei Cui. MicroRNA-21 depletion by CRISPR/Cas9 in CNE2 nasopharyngeal cells hinders proliferation and induces apoptosis by targeting the PI3K/AKT/MOTOR signaling pathway. *INT J CLIN EXP PATHO* . 2020 Apr 1;13(4):738-745
 55. Yi Yang, Xiaoyu Luo, Mayila Yasheng, Jun Zhao, Jinyu Li, Jinyao Li. Ergosterol peroxide from *Pleurotus ferulae* inhibits gastrointestinal tumor cell growth through induction of apoptosis via reactive oxygen species and endoplasmic reticulum stress. *Food Funct* . 2020 May 1;11(5):4171-4184.; doi: 10.1039/c9fo02454a

Version 2021.09.01