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1. Yongzhi Cui, Zhongying Li, Yuanyuan Guo, Xiangbei Qi, Yuehua Yang, Xiong Jia, Rui Li, Jingyu Shi, Weihang Gao, Zhengwei Ren, Guohui Liu, Qingsong Ye, Zhiping Zhang, Dehao Fu . Bioinspired Nanovesicles Convert the Skeletal Endothelium-Associated Secretory Phenotype to Treat Osteoporosis ACS Nano. 2022 Jul 26;16(7):11076-11091.
2. Jingzeng Cai, Jiaqiang Huang, Jie Yang, Xiaoming Chen, Haoran Zhang, Yue Zhu, Qi Liu, Ziwei Zhang . The protective effect of selenoprotein M on non-alcoholic fatty liver disease: the role of the AMPKα1-MFN2 pathway and Parkin mitophagy Cell Mol Life Sci. 2022 Jun 9;79(7):354.
3. Zhaosong Meng, Jiacheng Liu, Zhipeng Feng, Shuling Guo, Mingzhe Wang, Zheng Wang, Zhe Li, Hongjie Li, Lei Sui . N-acetylcysteine regulates dental follicle stem cell osteogenesis and alveolar bone repair via ROS scavenging Stem Cell Res Ther. 2022 Sep 8;13(1):466.
4. Qian Wang, Zhenzhen Sun, Shihan Cao, Xiuli Lin, Mengying Wu, Yuanyuan Li, Jie Yin, Wei Zhou, Songming Huang, Aihua Zhang, Yue Zhang, Weiwei Xia, Zhanjun Jia . Reduced Immunity Regulator MAVS Contributes to Non-Hypertrophic Cardiac Dysfunction by Disturbing Energy Metabolism and Mitochondrial Homeostasis Front Immunol. 2022 Jul 1:13:919038. doi: 10.3389/fimmu.2022.919038.
5. Lidan Cheng, Peng Yu, Fangfang Li, Xueling Jiang, Xiaojuan Jiao, Yunfeng Shen, Xiaoyang Lai . Human umbilical cord-derived mesenchymal stem cell-exosomal miR-627-5p ameliorates non-alcoholic fatty liver disease by repressing FTO expression Hum Cell. 2021 Nov;34(6):1697-1708.
6. Zhiwei Xie, Qiuping Xu, Lu Sun, Ruijing Li, Jizhou Shi, Qian Yang, Min Zong, Jianyong Qin . Effects of Y-27632 on the osteogenic and adipogenic potential of human dental pulp stem cells in vitro Hum Exp Toxicol. 2022 Jan-Dec;41:9603271221089003. doi: 10.1177/09603271221089003.
7. Songhan Li, Shuohui Dong, Qian Xu, Bowen Shi, Linchuan Li, Wenjie Zhang, Jiankang Zhu, Yugang Cheng, Guangyong Zhang, Mingwei Zhong . Sleeve Gastrectomy-Induced AMPK Activation Attenuates Diabetic Cardiomyopathy by Maintaining Mitochondrial Homeostasis via NR4A1 Suppression in Rats Front Physiol. 2022 Mar 10;13:837798. doi: 10.3389/fphys.2022.837798.
8. Helin Wang, Yuqian Shi, Feng He, Tao Ye, Shibin Yu, Hui Miao, Qian Liu, Mian Zhang . GDF11 inhibits abnormal adipogenesis of condylar chondrocytes in temporomandibular joint osteoarthritis Bone Joint Res. 2022 Jul;11(7):453-464.
9. Yue Wu, Yabin Qi, Yangqiu Bai, Haihui Zhang, Wenliang Zhu, Shengli Zhou, Yanrui Zhang . LncRNA 1700020I14Rik promotes AKR1B10 expression and activates Erk pathway to induce hepatocyte damage in alcoholic hepatitis Cell Death Discov. 2022 Aug 26;8(1):374.
10. Meng-Ke Liu, Li-Li Cheng, Hong-Mei Yi, Yang He, Xiao Li, Di Fu, Yu-Ting Dai, Hai Fang, Shu Cheng, Peng-Peng Xu, Ying Qian, Yan Feng, Qian Liu, Li Wang, Wei-Li Zhao . Enhanced lipid metabolism confers the immunosuppressive tumor microenvironment in CD5-positive non-MYC/BCL2 double expressor lymphoma Front Oncol. 2022 Oct 6:12:885011.
11. Jian Xiao, Kuan Shen, Kanghui Liu, Yuanhang Wang, Hao Fan, Quan Cheng, Xinyi Zhou, Li Hu, Gang Wang, Zekuan Xu, Li Yang . Obesity promotes lipid accumulation in lymph node metastasis of gastric cancer: a retrospective case-control study Lipids Health Dis. 2022 Nov 17;21(1):123.
12. Wanying Liu, Ke Li, Menglin Zheng, Ling He, Tong Chen . Genipin Attenuates Diabetic Cognitive Impairment by Reducing Lipid Accumulation and

Promoting Mitochondrial Fusion via FABP4/Mfn1 Signaling in Microglia Antioxidants (Basel). 2022 Dec 29;12(1):74.

13. Chang-Jing Huang, Chen-Yue Zhang, Ying-Ke Zhao, Dan Wang, Liping Zhuang, Ling Qian, Lin Xie, Ying Zhu, Zhi-Qiang Meng . Bufalin Inhibits Tumorigenesis and SREBP-1-Mediated Lipogenesis in Hepatocellular Carcinoma via Modulating the ATP1A1/CA2 Axis Am J Chin Med. 2023;51(2):461-485.
14. Si Qin, Qian Su, Xiang Li, Muqing Shao, Yindi Zhang, Fadong Yu, Yinxing Ni, Jian Zhong . Curcumin suppresses cell proliferation and reduces cholesterol absorption in Caco-2 cells by activating the TRPA1 channel Lipids Health Dis. 2023 Jan 14;22(1):6.
15. Ying Wang, Lingbin Che, Xi Chen, Zirui He, Dianwen Song, Yuan Yuan, Changsheng Liu . Repurpose dasatinib and quercetin: Targeting senescent cells ameliorates postmenopausal osteoporosis and rejuvenates bone regeneration Bioact Mater. 2023 Jan 20;25:13-28.
16. Li-Sha Ran, Ya-Zeng Wu, Yi-Wen Gan, Hong-Lian Wang, Li-Juan Wu, Chun-Mei Zheng, Yao Ming, Ran Xiong, Yong-Lin Li, Shi-Hang Lei, Xue Wang, Xiao-Qing Lao, Hong-Min Zhang, Li Wang, Chen Chen, Chang-Ying Zhao . Andrographolide ameliorates hepatic steatosis by suppressing FATP2-mediated fatty acid uptake in mice with nonalcoholic fatty liver disease J Nat Med. 2023 Jan;77(1):73-86.
17. Xue Zhang, Ying Zhi, Xinyan Zan, Kerui Fan, Kun Chen, Shuang Zhao, Xinyi Dai, Longjiang Li, Yongqiang Yang, Kai Hu, Xianqiong Gong, Li Zhang . Immune response gene 1 deficiency aggravates high fat diet-induced nonalcoholic fatty liver disease via promotion of redox-sensitive AKT suppression Biochim Biophys Acta Mol Basis Dis. 2023 Apr;1869(4):166656.
18. Yu Zhou, Huakan Zhao, Ran Ren, Mingyue Zhou, Jiangang Zhang, Zhijuan Wu, Yu Chen, Juan Lei, Yang Chen, Ying Yu, Yongsheng Li . GPAT3 is a potential therapeutic target to overcome sorafenib resistance in hepatocellular carcinoma Theranostics. 2024 Jun 1;14(9):3470-3485.
19. Ning Wang, Kelan Yuan, Shuo Yang, Xiuming Jin . 13- cis Retinoic Acid-Mediated Modulation of Human Meibomian Gland Epithelial Cells Development: Implications for In Vitro Modeling of Meibomian Gland Dysfunction J Ocul Pharmacol Ther. 2024 Dec;40(10):659-667.
20. Yu Cao, Daodao Wang, Dingzhou Zhou . MSC Promotes the Secretion of Exosomal lncRNA KLF3-AS1 to Regulate Sphk1 Through YY1-Musashi-1 Axis and Improve Cerebral Ischemia-Reperfusion Injury Mol Neurobiol. 2024 Dec;61(12):10462-10480.

Version 2022.07.26