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1. Xiaotong Wang,Shiming Xie,Zili Li,Zhen Ye,Xiuli Gu,Liquan Zhou,Honggang Li.Generation of an iPSC line (HUSTi002-A) from fibroblasts of a patient with Sertoli cell-only syndrome carrying c.731_732delAT in PIWIL2 gene Stem Cell Res. 2020 Jan;42:101703.;doi: 10.1016/j.scr.2020.101703

2. Gang Yang,Guangbing Xiong,Mengyu Feng,Fangyu Zhao,Jiangdong Qiu,Yueze Liu,Zhe Cao,Huanyu Wang,Jinshou Yang,Lei You,Lianfang Zheng,Taiping Zhang,Yupeai Zhao.OLR1 Promotes Pancreatic Cancer Metastasis via Increased c-Myc Expression and Transcription of HMGA2 Mol Cancer Res. 2020 May;18(5):685-697.;doi: 10.1158/1541-7786.MCR-19-0718

3. Wei Tang,Meng Li,Xin Qi,Jing Li. β 1,4-Galactosyltransferase V Modulates Breast Cancer Stem Cells through Wnt/ β -catenin Signaling Pathway Cancer Res Treat. 2020 Oct;52(4):1084-1102.;doi: 10.4143/crt.2020.093

4. Shibin Qu,Kunwei Niu,Jianlin Wang,Jimin Dai,Anutosh Ganguly,Chao Gao,Yuzi Tian,Zhibin Lin,Xisheng Yang,Xuan Zhang,Zhengcai Liu,Haimin Li.LINC00671 suppresses cell proliferation and metastasis in pancreatic cancer by inhibiting AKT and ERK signaling pathway Cancer Gene Ther. 2021 Apr;28(3-4):221-233.;doi: 10.1038/s41417-020-00213-4

5. Mu Zhang,Chunjie Yang,Meng Zhu,Li Qian,Yan Luo,Huimin Cheng,Rong Geng,Xiaojun Xu,Cheng Qian,Yu Liu.Saturated fatty acids entrap PDX1 in stress granules and impede islet beta cell function Diabetologia. 2021 May;64(5):1144-1157.;doi: 10.1007/s00125-021-05389-4

6. Zhiyong Xiong, Wei Xiong, Wen Xiao, Changfei Yuan, Jian Shi, Yu Huang, Cheng Wang, Xiangui Meng, Zhixian Chen, Hongmei Yang, Ke Chen, Xiaoping Zhang. NNT-induced tumor cell "slimming" reverses the pro-carcinogenesis effect of HIF2a in tumors. Clin Transl Med. 2021 Jan;11(1):e264.

7. Xiaonan Wang, Qianying Guo, Hao Wang, Xiaodong Yuan, Bijun Wang, Peter E Lobie, Tao Zhu, Sheng Tan, Zhengsheng Wu. PCBP2 Posttranscriptional Modifications Induce Breast Cancer Progression via Upregulation of UFD1 and NT5E. Mol Cancer Res. 2021 Jan;19(1):86-98.

8. Xin Zheng, Rui Liu, Chenchen Zhou, Haopeng Yu, Wanyi Luo, Jianhui Zhu, Jiaxin Liu, Zhe Zhang, Na Xie, Xian Peng, Xin Xu, Lei Cheng, Quan Yuan, Canhua Huang, Xuedong Zhou. ANGPTL4-Mediated Promotion of Glycolysis Facilitates the Colonization of Fusobacterium nucleatum in Colorectal Cancer. Cancer Res. 2021 Dec 15;81(24):6157-6170.

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10. Xianjie Lu, Na Song, Wei Wang, Yanming Liu, Hao Song, Li Xu, Yan Wang, Chuanfei Wei, Juanli Chen, Xiaofei Yang, Fabin Han.. Generation of integration-free human iPSC line LCPHi001-A from a Parkinson's disease patient carrying the RecNciI mutation in GBA gene. Stem Cell Res. 2021 Oct;56:102514.

11. Ping Liu, Yifei Bian, Jia Zhong, Yang Yang, Xiang Mu, Zhongjie Liu. Establishment and characterization of a rat intestinal microvascular endothelial cell line. Tissue Cell. 2021 Oct;72:101573.

12. Yue Wang, Weiyao Xiong, Shuangxia Zhao, Bin Li, Alex Chia Yu Chang. Generation of two induced pluripotent stem cell lines, SHIPMi001-A from a patient with hypertrophic cardiomyopathy caused by MYBPC3 gene mutation and SHIPMi002-A from a healthy male individual. Stem Cell Res. 2021 Nov 10;57:102594.

13. Ziwei Yuan, Bin Wang, Yilong Teng, William Ho, Bin Hu, Kofi Oti Boakye-Yiadom, Xiaoyang Xu, Xue-Qing Zhang. Rational design of engineered H-ferritin nanoparticles with improved siRNA delivery efficacy across an in vitro model of the mouse BBB. Nanoscale. 2022 May 5;14(17):6449-6464.

14. Xiaoyan Zhang, Ling Yang, Wanjun Lei, Qiang Hou, Ming Huang, Rongjing Zhou, Tariq Enver, Shixiu Wu. Single-cell sequencing reveals CD133+CD44--originating evolution and novel stemness related variants in human colorectal cancer. EBioMedicine. 2022 Aug;82:104125.

15. Xianjie Lu, Wei Wang, Yanming Liu, Na Song, Mengpeng Li, Xin Mu, Nan Zhang, Qingfa Chen, Licheng Jiang, Xianglin Kong, Peng Sun, Jiabei Tong, Yunping Zhang, Jingtao Li, Shengjun Ma, Fabin Han. Establishment and characterization of human induced pluripotent stem cell line from a Parkinson's disease patient harboring VPS13A gene mutation. Stem Cell Res. 2022 Apr;60:102685.

16. Maisumu Gulimiheranmu, Shuang Li, Junmei Zhou. Generation of a MIR5004 knockout cell line from human induced pluripotent stem cells by CRISPR/Cas9 gene editing. Stem Cell Res. 2022 Jul;62:102805.

17. Shuyu Fang, Anle Zeng, Qiling Xu, Lina Zhou, Zhiyong Zhang, Yunfei An, Xiaodong Zhao. Generation of human induced pluripotent stem cell line from peripheral blood mononuclear cells from an activated phosphoinositide 3-kinase δ syndrome patient. Stem Cell Res. 2022 Jul;62:102822.