碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

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- Zhang B, Peng X, Li G, Xu Y, Xia X, Wang Q. Oxidative stress is involved in Patulin induced apoptosis in HEK293 cells. *Toxicol.* 2015 Feb;94:1-7.
- Chen X, Li M, Li L, Xu S, Huang D, Ju M, Huang J, Chen K, Gu H. Trehalose, sucrose and raffinose are novel activators of autophagy in human keratinocytes through anmTOR-independent pathway. *Sci Rep.* 2016 Jun 22;6:28423.
- Xu S, Li L, Li M, Zhang M, Ju M, Chen X, Gu H. Impact on Autophagy and Ultraviolet B Induced Responses of Treatment withthe MTORInhibitors Rapamycin, Everolimus, Torin 1,and pp242 in Human Keratinocytes. *Oxid Med Cell Longev.* 2017;2017:5930639.
- Yang H,Wang S,Kang YJ,Wang C,Xu Y,Zhang Y,Jiang Z. Long non-coding RNA SNHG1 predicts a poor prognosis and promotes colon cancer tumorigenesis. *Oncol Rep.* 2018 Jul;40(1):261-271.
- Song N,Wang B,Feng G,Duan L,Yuan S,Jia W,Liu Y. Knockdown of high mobility group box 3 impairs cell viability and colony formation but increases apoptosis in A549 human non-small cell lung cancer cells. *Oncol Lett.* 2019 Mar;17(3):2937-2945.
- Ma Y,Bu D,Long J,Chai W,Dong J. LncRNA DSCAM-AS1 acts as a sponge of miR-137 to enhance Tamoxifen resistance in breast cancer. *J Cell Physiol.* 2019 Mar;234(3):2880-2894.
- Li K,Zhou F,Zhou Y,Zhang S,Li Q,Li Z,Liu L,Wu G,Meng R. Quinalizarin, a specific CK2 inhibitor, can reduce icotinib resistance in human lung adenocarcinoma cell lines. *Int J Mol Med.* 2019 Aug;44(2):437-446.
- Yan Qin, Peiling Ni, Qingye Zhang, Xiao Wang, Xiaoling Du, Zixi Yin, Lingling Wang, Lihong Ye, Lingyi Chen. Hbxip is essential for murine embryogenesis and regulates embryonic stem cell differentiation through activating mTORC1. *Development.* 2022 Jun 15;149(12):dev200527.
- Yang Yu, Yanwen Li, Xinjie Yang, Qiuyi Deng, Bin Xu, Hua Cao, Jianwen Mao. A Novel Imidazo[1,2-a]pyridine Compound Reduces Cell Viability and Induces Apoptosis of HeLa Cells by p53/Bax-Mediated Activation of Mitochondrial Pathway. *Anticancer Agents Med Chem.* 2022;22(6):1102-1110.

10. Jinsheng Liu, Li Li, Ping Xie, Xiaoyan Zhao, Dongjing Shi, Yan Zhang, Chuxiong Pan, Tianzuo Li. Sevoflurane induced neurotoxicity in neonatal mice links to a GSK3 β /Drp1-dependent mitochondrial fission and apoptosis. *Free Radic Biol Med.* 2022 Mar;181:72-81.
11. Haoya Chang, Jian Li, Chunhui Zhang, Weiping Qian. Octocrylene exposure impairs mouse oocyte quality by inducing spindle defects and mitochondria dysfunction. *Toxicology.* 2022 Sep;479:153306.
12. Zhaoyang Wang, Ziyi Li, Yaling Lei, Yuan Liu, Yuqing Feng, Derong Chen, Siying Ma, Ziyang Xiao, Meirong Hu, Jingxian Deng, Yuxin Wang, Qihao Zhang, Yadong Huang, Yan Yang. Recombinant Photolyase-Thymine Alleviated UVB-Induced Photodamage in Mice by Repairing CPD Photoproducts and Ameliorating Oxidative Stress. *Antioxidants (Basel).* 2022 Nov 22;11(12):2312.
13. Huifang Sun, Zhuoya Wang, Qi Zhang, Na Chen, Mibo Tang, Zihua Yang, Yuming Xu, Jiansheng Kang, Yanlin Wang. Golgi damage caused by dysfunction of PiT-2 in primary familial brain calcification. *Biochem Biophys Res Commun.* 2023 Jan 29;642:167-174.
14. Shun Yang, Yuanxin Ma, Xiaocong Lou, Zhewei Zhou, Huimin Zhang, Shunfa Yi, Yan Cheng, Shichao Qian, Mengmeng Huang, Hui Fei. The role of TNF- α in the phagocytosis of largemouth bass (*Micropterus salmoides*) leukocytes. *Fish Shellfish Immunol.* 2023 Jan;132:108488.

Version 2024.03.12