

一步法TUNEL细胞凋亡检测试剂盒

产品编号	产品名称	包装
C1088	一步法TUNEL细胞凋亡检测试剂盒(绿色荧光)	50次

产品简介:

- 碧云天生产的一步法TUNEL细胞凋亡检测试剂盒(One Step TUNEL Apoptosis Assay Kit)为您提供了一种高灵敏度又快速简便的细胞凋亡检测方法。对于经过固定和洗涤的细胞或组织, 只要经过一步染色反应, 洗涤后就可以通过荧光显微镜或流式细胞仪检测到呈现绿色荧光的凋亡细胞。
- 细胞在发生凋亡时, 会激活一些DNA内切酶, 这些内切酶会切断核小体间的基因组DNA。细胞凋亡时抽提DNA进行电泳检测, 可以发现180-200bp的DNA ladder。基因组DNA断裂时, 暴露的3'-OH可以在末端脱氧核苷酸转移酶(Terminal Deoxynucleotidyl Transferase, TdT)的催化下加上绿色荧光探针荧光素(FITC)标记的dUTP(fluorescein-dUTP), 从而可以通过荧光显微镜或流式细胞仪进行检测, 这就是TUNEL(TdT-mediated dUTP Nick-End Labeling)法检测细胞凋亡的原理。
注: FITC是fluorescein isothiocyanate的缩写, 实际上大多数情况下所谓的FITC即为fluorescein。
- 本试剂盒有如下优点。(1) 高灵敏度: 背景染色极低, 阳性染色明亮, 可以在单细胞水平检测到细胞凋亡, 同时由于凋亡早期就有DNA断裂, 可以检测到早期的细胞凋亡。(2) 特异性: TUNEL检测时通常更容易标记凋亡细胞, 而不容易标记坏死细胞。(3) 快速: 仅需约1-2个小时即可完成。(4) 方便: 只需一步染色反应, 洗涤后即可观察, 不必使用二抗等进行多步操作。(5) 应用范围广: 可以用于检测冷冻或石蜡切片中的细胞凋亡情况, 也可以检测培养的贴壁细胞或悬浮细胞的凋亡情况。(6) 实测效果好: 参考图1, 本试剂盒的实际检测效果和国外同类产品相当。




图1. 本产品的检测效果以及和Roche公司同类产品检测效果的比较。A. 碧云天的本试剂盒和Roche公司同类试剂盒分别检测正

常HeLa细胞和经DNase I处理的HeLa细胞的效果图，图中的绿色荧光为TUNEL染色阳性细胞。B. 碧云天的本试剂盒和Roche公司同类试剂盒分别检测正常HeLa细胞和经10μM喜树碱(Camptothecin)处理24小时的HeLa细胞的效果图。无论是经DNase I处理的阳性细胞样品，还是使用10μM喜树碱处理的阳性细胞样品，碧云天的本试剂盒和Roche公司的同类试剂盒的效果相当。本图中的染色实验均在12孔板中进行。本图仅作参考，不同的样品不同的检测条件，实际获得的结果可能有较明显的差别。

- TUNEL法特异性检测细胞凋亡时产生的DNA断裂，但不会检测出射线等诱导的DNA断裂(和细胞凋亡时的断裂方式不同)。这样一方面可以把凋亡和坏死区分开，另一方面也不会把射线等诱导发生DNA断裂的非凋亡细胞判断为凋亡细胞。
- 极少数细胞凋亡时没有DNA断裂，此时不适用TUNEL法检测。在个别类型的坏死细胞中也发现TUNEL检测呈阳性。在需要严格判断细胞凋亡的情况下，最好同时检测多个凋亡指标。
- 本试剂盒足够检测50个样品。

包装清单：

产品编号	产品名称	包装
C1088-1	TdT酶	250μl
C1088-2	荧光标记液	2×1.2ml
C1088-3	TdT酶稀释液(选用)	1ml
—	说明书	1份

保存条件：

-20℃保存，荧光标记液需避光保存。

注意事项：

- 需自备用于洗涤细胞的PBS或HBSS，用于封片的抗荧光淬灭封片液(P0126)，用于固定的4%多聚甲醛或向碧云天订购免疫染色固定液(P0098)，同时需自备含0.3% Triton X-100的PBS或向碧云天订购免疫染色强力通透液(P0097)。
- 如果用于石蜡切片的检测，需自备蛋白酶K，二甲苯。蛋白酶K(ST533)可以向碧云天订购。
- 本产品仅限于专业人员的科学研究用，不得用于临床诊断或治疗，不得用于食品或药品，不得存放于普通住宅内。
- 为了您的安全和健康，请穿实验服并戴一次性手套操作。

使用说明：

1. 对于贴壁细胞或细胞涂片：

- a. PBS或HBSS洗涤一次。
- b. 如果细胞贴得不牢，可以干燥样品使细胞贴得更牢。
- c. 用碧云天生产的免疫染色固定液(P0098) 或4%多聚甲醛固定细胞30分钟。
- d. 用PBS或HBSS洗涤一次。
- e. 加入碧云天生产的免疫染色强力通透液(P0097) 或含0.3% Triton X-100的PBS，室温孵育5分钟。
- f. 转步骤5。

2. 对于悬浮细胞或细胞悬液：

- a. 收集细胞(不超过200万细胞)，PBS或HBSS洗涤一次。
- b. 用碧云天生产的免疫染色固定液(P0098)或4%多聚甲醛固定细胞30分钟。为防止细胞聚集成团，宜在侧摆摇床或水平摇床上缓慢摇动的同时进行固定。
- c. 用PBS或HBSS洗涤一次。
- d. 用碧云天生产的免疫染色强力通透液(P0097)或含0.3% Triton X-100的PBS重悬细胞，室温孵育5分钟。
- e. 转步骤5。

3. 对于石蜡切片：

- a. 二甲苯中脱蜡5-10分钟。换用新鲜的二甲苯，再脱蜡5-10分钟。无水乙醇5分钟。90%乙醇2分钟。70%乙醇2分钟，蒸馏水2分钟。
- b. 滴加20μg/ml不含DNase的蛋白酶K(推荐使用碧云天的ST532/ST533 蛋白酶K(20mg/ml)，用P0106 免疫染色洗涤液或10mM Tris-HCl pH7.4-7.8稀释1000倍即为20μg/ml不含DNase的蛋白酶K)，20-37℃作用15-30分钟(不同组织的最佳作用温度和时间需自行摸索)。
- c. PBS或HBSS洗涤3次。注意：这一步必须把蛋白酶K洗涤干净，否则会严重干扰后续的标记反应。
- d. 转步骤5。

4. 对于冰冻切片：

- a. 用碧云天生产的免疫染色固定液(P0098)或4%多聚甲醛固定细胞30-60分钟。
- b. PBS或HBSS洗涤2次，每次10分钟。
- c. 加入碧云天生产的免疫染色强力通透液(P0097)或含0.5% Triton X-100的PBS，室温孵育5分钟。
- d. 转步骤5。

5. 配制TUNEL检测液：

参考下表配制适当量的TUNEL检测液，需充分混匀。注意：配制好的TUNEL检测液必须一次使用完毕，不能冻存。

	1个样品	5个样品	10个样品
TdT酶	5μl	25μl	50μl
荧光标记液	45μl	225μl	450μl
TUNEL检测液	50μl	250μl	500μl

6. 对于贴壁细胞、细胞涂片或组织切片：

- 用PBS或HBSS洗涤2次。
- 在样品上加50μl TUNEL检测液，37℃避光孵育60分钟。注意：50μl TUNEL检测液适合涂片、切片或96孔板、48孔板、24孔板或12孔板的一个孔，如果是6孔板中的一个孔TUNEL检测液宜使用100μl。如果待检测的样品为涂片、切片或在24孔板、12孔板或6孔板中，可以使用防蒸发膜，或自行尝试使用自封袋或者其它适当材料自行裁剪成比孔略小的圆形塑料片，滴加TUNEL检测液后覆盖在样品上，可以防止TUNEL检测液蒸发，并且使TUNEL检测液均匀覆盖样品。自行裁剪圆片时需要连着圆片突出一个角或连着一一条，并将圆形之外的突出部分折叠，方便染色结束后利用突出部分顺利取出圆片。也可以用PAP Pen圈出一个区域进行染色。孵育时需注意在多余的孔和多孔板的空隙中加入适量水以保持湿润，从而尽量减少TUNEL检测液的蒸发。
- PBS或HBSS洗涤3次。
- 用抗荧光淬灭封片液封片后荧光显微镜下观察。可以使用的激发波长范围为450-500nm，发射波长范围为515-565nm(绿色荧光)。染色效果可参见图1。

7. 对于悬浮细胞或细胞悬液：

- 用PBS或HBSS洗涤2次。
- 加入50μl TUNEL检测液，37℃避光孵育60分钟。
- PBS或HBSS洗涤2次。
- 用250-500μl PBS或HBSS悬浮。
- 此时可以用流式细胞仪进行检测或涂片后在荧光显微镜下观察。可以使用的激发波长范围为450-500nm，发射波长范围为515-565nm(绿色荧光)。

常见问题：

1. 出现非特异性荧光标记。

- 有些细胞或组织，例如平滑肌细胞或组织，nuclease或polymerase的酶活性水平较高，易导致出现非特异性的荧光标记。解决方法是，取细胞或组织后立即固定并且要充分固定，以阻止这些酶导致假阳性。
- 使用了不适当的固定液，例如一些酸性固定液，导致出现假阳性。建议采用推荐的固定液。
- TUNEL检测反应时间过长，或TUNEL检测反应过程中反应液渗漏，细胞或组织表面不能保持湿润，也可能出现非特异性荧光。注意控制反应时间，并确保TUNEL检测反应液能很好地覆盖样品。

2. 荧光背景很高。

- 支原体污染。请使用支原体染色检测试剂盒检测是否为支原体污染。支原体染色检测试剂盒(C0296)可以向碧云天订购。
- 高速分裂和增殖的细胞，有时也会出现细胞核中的DNA断裂。
- TUNEL反应过强。可以用试剂盒提供的TdT酶稀释液稀释TdT酶2-5倍后再按照说明书操作。稀释后的TdT酶需当日使用。
- 红细胞中血红蛋白导致的自发荧光产生严重干扰。此时宜选择其它细胞凋亡检测试剂盒。

3. 标记效率低。

- 使用乙醇或甲醇固定会导致标记的效率较低。
- 固定时间过长，导致交联程度过高。此时宜减少固定时间。
- 荧光淬灭。荧光在普通光照10分钟就会严重淬灭。解决方法是需注意避光操作。
- 碘化丙啶双染时，如果碘化丙啶染色过深会导致观察到的本试剂盒的TUNEL染色效果减弱。碘化丙啶可以接受fluorescein激发产生的荧光，从而起到淬灭作用。解决方法是用较低浓度的碘化丙啶染色，例如0.5μg/ml碘化丙啶。
- 贴壁细胞如果使用药物诱导凋亡，会使发生凋亡细胞的贴壁性会减弱，所以建议在凋亡诱导结束后，用可以对多孔板进行离心的离心机1000g离心5分钟，然后再吸除培养基并用PBS洗涤。如果没有适合的离心机，请注意操作轻缓，防止发生凋亡的细胞在洗涤时洗去。后续整个操作也需要轻缓。

相关产品：

产品编号	产品名称	包装
C1086	一步法TUNEL细胞凋亡检测试剂盒(绿色荧光)	20次
C1088	一步法TUNEL细胞凋亡检测试剂盒(绿色荧光)	50次
C1089	一步法TUNEL细胞凋亡检测试剂盒(红色荧光)	20次
C1090	一步法TUNEL细胞凋亡检测试剂盒(红色荧光)	50次
C1091	TUNEL细胞凋亡检测试剂盒(显色法)	20次
C1098	TUNEL细胞凋亡检测试剂盒(显色法)	50次
C1062	Annexin V-FITC细胞凋亡检测试剂盒	20次
C1063	Annexin V-FITC细胞凋亡检测试剂盒	50次

C1065	Annexin V-PE细胞凋亡检测试剂盒	20次
C1067	Annexin V-EGFP细胞凋亡检测试剂盒	20次
C1068	Annexin V-EGFP细胞凋亡检测试剂盒	50次
C1082	TUNEL检测阳性对照制备试剂盒	10次
P0098	免疫染色固定液	100ml
P0097-100ml	免疫染色强力通透液	100ml
P0097-500ml	免疫染色强力通透液	500ml
P0126	抗荧光淬灭封片液	5ml
ST533	Proteinase K (20mg/ml)	1ml

使用本产品的文献：

1. Li Y, Li Q, Wang Z, Liang D, Liang S, Tang X, Guo L, Zhang R, Zhu D. 15-HETE suppresses K(+) channel activity and inhibits apoptosis in pulmonary artery smooth muscle cells. Apoptosis. 2009 Jan;14(1):42-51.

2. Yin Z, Chen X, Chen JL, Shen WL, Hieu Nguyen TM, Gao L, Ouyang HW. The regulation of tendon stem cell differentiation by the alignment of nanofibers. Biomaterials. 2010;31(8):2163-75.

3. Zhang H, Wang W, Yin H, Zhao X, Du Y. Oligochitosan induces programmed cell death in tobacco suspension cells. Carbohydrate Polymers. 2012 Feb;87(3):2270–8.

4. Li XY, Jing CQ, Lei WL, Li J, Wang JJ. Apoptosis caused by imidazolium-based ionic liquids in PC12 cells. ECOTOX ENVIRON SAFE. 2012 Sep;83(1):102-7.

5. Li ZL, Hu J, Li YL, Xue F, Zhang L, Xie JQ, Liu ZH, Li H, Yi DH, Liu JC, Wang SW. The effect of hyperoside on the functional recovery of the ischemic/reperfused isolated rat heart: potential involvement of the extracellular signal-regulated kinase 1/2 signaling pathway. FREE RADICAL BIO MED. 2013 Apr;57:132-40.

6. Hu X, Sui X, Li L, Huang X, Rong R, Su X, Shi Q, Mo L, Shu X, Kuang Y, Tao Q, He C. Protocadherin 17 acts as a tumour suppressor inducing tumour cell apoptosis and autophagy, and is frequently methylated in gastric and colorectal cancers. J Pathol. 2013 Jan;229(1):62-73.

7. Zhang Y, Lin A, Zhang C, Tian Z, Zhang J. Phosphorothioate-modified CpG oligodeoxynucleotide (CpG ODN) induces apoptosis of human hepatocellular carcinoma cells independent of TLR9. CANCER IMMUNOL IMMUN. 2014 Apr;63(4):357-67.

8. Zhang Y, Lin A, Sui Q, Zhang C, Tian Z, Zhang J. Phosphorothioate modification of the TLR9 ligand CpG ODN inhibits poly(I:C)-induced apoptosis of hepatocellular carcinoma by entry blockade. Cancer Lett. 2014 Dec 1;355(1):76-84.

9. Guo Q, Lan P, Yu X, Han Q, Zhang J, Tian Z, Zhang C. Immunotherapy for hepatoma using a dual-function vector with both immunostimulatory and pim-3-silencing effects. Mol Cancer Ther. 2014 Jun;13(6):1503-13.

10. Zheng L, Fu Y, Zhuang L, Gai R, Ma J, Lou J, Zhu H, He Q, Yang B. Simultaneous NF-κB inhibition and E-cadherin upregulation mediate mutually synergistic anticancer activity of celastrol and SAHA in vitro and in vivo. Int J Cancer. 2014 Oct 1;135(7):1721-32.

11. Yu C, Cui S, Zong C, Gao W, Xu T, Gao P, Chen J, Qin D, Guan Q, Liu Y, Fu Y, Li X, Wang X. The Orphan Nuclear Receptor NR4A1 Protects Pancreatic β-Cells from Endoplasmic Reticulum (ER) Stress-mediated Apoptosis. J Biol Chem. 2015 Aug 21;290(34):20687-99.

12. Li PL, Zhang X, Wang LL, Du LT, Yang YM, Li J, Wang CX. MicroRNA-218 is a prognostic indicator in colorectal cancer and enhances 5-fluorouracil-induced apoptosis by targeting BIRC5. Carcinogenesis. 2015 Dec;36(12):1484-93.

13. Huang J, Zhang L, He C, Qu Y, Li J, Zhang J, Du T, Chen X, Yu Y, Liu B, Zhu Z. Claudin-1 enhances tumor proliferation and metastasis by regulating cell anoikis in gastric cancer. ONCOTARGET. 2015 Jan 30;6(3):1652-65.

14. Liu D, Wen Y, Tang TT, Lv LL, Tang RN, Liu H, Ma KL, Crowley SD, Liu BC. Megalin/Cubulin-Lysosome-mediated Albumin Reabsorption Is Involved in the Tubular Cell Activation of NLRP3 Inflammasome and Tubulointerstitial Inflammation. J Biol Chem. 2015 Jul 17;290(29):18018-28.

15. Wang L, Xu Y, Li H, Lei H, Guan R, Gao Z, Xin Z. Antioxidant icarisperide II combined with insulin restores erectile function in streptozotocin-induced type 1

diabetic rats. J Cell Mol Med. 2015 May;19(5):960-9.

16. Ding Y, Wang H, Niu J, Luo M, Gou Y, Miao L, Zou Z, Cheng Y. Induction of ROS Overload by Alantolactone Prompts Oxidative DNA Damage and Apoptosis in Colorectal Cancer Cells. Int J Mol Sci. 2016 Apr 14;17(4):558.

17. Du JK, Cong BH, Yu Q, Wang H, Wang L, Wang CN, Tang XL, Lu JQ, Zhu XY, Ni X. Upregulation of microRNA-22 contributes to myocardial ischemia-reperfusion injury by interfering with the mitochondrial function. FREE RADICAL BIO MED. 2016 Jul;96:406-17.

18. Guo R, Wu Z, Jiang J, Liu C, Wu B, Li X, Li T, Mo H, He S, Li S, Yan H, Huang R, You Q, Wu K. New mechanism of lipotoxicity in diabetic cardiomyopathy: Deficiency of Endogenous H2S Production and ER stress. Mech Ageing Dev. 2016 Nov 18. pii: S0047-6374(16)30222-6.

19. Dong R, Gong Y, Meng W, Yuan M, Zhu H, Ying M, He Q, Cao J, Yang B. The involvement of M2 macrophage polarization inhibition in fenretinide-mediated chemopreventive effects on colon cancer. Cancer Lett. 2016 Nov 30;388:43-53.

20. Guo W, Wang Y, Wang Z, Wang YP, Zheng H. Inhibiting autophagy increases epirubicin's cytotoxicity in breast cancer cells. Cancer Sci. 2016 Nov;107(11):1610-1621.

21. Yan N, Liu Y, Liu S, Cao S, Wang F, Wang Z, Xi S. Fluoride-Induced Neuron Apoptosis and Expressions of Inflammatory Factors by Activating Microglia in Rat Brain. Mol Neurobiol. 2016 Sep;53(7):4449-60.

22. Li T, Ding ZL, Zheng YL, Wang W. MiR-484 promotes non-small-cell lung cancer (NSCLC) progression through inhibiting Apaf-1 associated with the suppression of apoptosis. Biomed Pharmacother. 2017 Dec;96:153-164.

23. Song XL, Zhang YJ, Wang XF, Zhang WJ, Wang Z, Zhang F, Zhang YJ, Lu JH, Mei JW, Hu YP, Chen L, Li HF, Ye YY, Liu YB, Gu J. Casticin induces apoptosis and G0/G1 cell cycle arrest in gallbladder cancer cells. Cancer Cell Int. 2017 Jan 5;17:9.

24. Xu Y, Leng K, Li Z, Zhang F, Zhong X, Kang P, Jiang X, Cui Y. The prognostic potential and carcinogenesis of long non-coding RNA TUG1 in human cholangiocarcinoma. ONCOTARGET. 2017 Jul 22;8(39):65823-65835.

25. Zhang Z, Zhang L, Sun C, Kong F, Wang J, Xin Q, Jiang W, Li K, Chen O, Luan Y. Baicalin attenuates monocrotaline-induced pulmonary hypertension through bone morphogenetic protein signaling pathway. ONCOTARGET. 2017 Jun 28;8(38):63430-63441.

26. Pei N, Mao Y, Wan P, Chen X, Li A, Chen H, Li J, Wan R, Zhang Y, Du H, Chen B, Jiang G, Xia M, Sumners C, Hu G, Gu D, Li H. Angiotensin II type 2 receptor promotes apoptosis and inhibits angiogenesis in bladder cancer. J EXP CLIN CANC RES. 2017 Jun 9;36(1):77.

27. Guo R, Wu Z, Jiang J, Liu C, Wu B, Li X, Li T, Mo H, He S, Li S, Yan H, Huang R, You Q, Wu K. New mechanism of lipotoxicity in diabetic cardiomyopathy: Deficiency of Endogenous H2S Production and ER stress. Mech Ageing Dev. 2017 Mar;162:46-52.

28. Liao T, Xue Y, Zhao D, Li S, Liu M, Chen J, Brand DD, Zheng H, Zhang Y, Zheng SG, Sun Q. In Vivo Attenuation of Antibody-Mediated Acute Renal Allograft Rejection by Ex Vivo TGF-β-Induced CD4+Foxp3+ Regulatory T Cells. Front Immunol. 2017 Oct 16;8:1334.

29. Pu Y, Meng K, Gu C, Wang L, Zhang X. Thrombospondin-1 modified bone marrow mesenchymal stem cells (BMSCs) promote neurite outgrowth and functional

- p>recovery in rats with spinal cord injury.
- ONCOTARGET*
- . 2017 Oct 24;8(56):96276-96289.
30. Xue X, Wang L, Meng X, Jiao J, Dang N. Regulator of G protein signaling 4 inhibits human melanoma cells proliferation and invasion through the PI3K/AKT signaling pathway. *ONCOTARGET*. 2017 Sep 11;8(45):78530-78544.
31. Li Z, Xu K, Zhang N, Amador G, Wang Y, Zhao S, Li L, Qiu Y, Wang Z. Overexpressed SIRT6 attenuates cisplatin-induced acute kidney injury by inhibiting ERK1/2 signaling. *Kidney Int*. 2018 Apr;93(4):881-892.
32. Fu Y, Li C, Luo Y, Li L, Liu J, Gui R. Silencing of Long Non-coding RNA MIAT Sensitizes Lung Cancer Cells to Gefitinib by Epigenetically Regulating miR-34a. *Front Pharmacol*. 2018 Feb 13;9:82.
33. Yin J, Ni B, Liao WG, Gao YQ. Hypoxia-induced apoptosis of mouse spermatocytes is mediated by HIF-1 α through a death receptor pathway and a mitochondrial pathway. *J Cell Physiol*. 2018 Feb;233(2):1146-1155.
34. Pan Y, Wang N, Xia P, Wang E, Guo Q, Ye Z. Inhibition of Rac1 ameliorates neuronal oxidative stress damage via reducing Bcl-2/Rac1 complex formation in mitochondria through PI3K/Akt/mTOR pathway. *Exp Neurol*. 2018 Feb;300:149-166.
35. Zhang LN, Tian H, Zhou XL, Tian SC, Zhang XH, Wu TJ. Upregulation of microRNA-351 exerts protective effects during sepsis by ameliorating skeletal muscle wasting through the Tead-4-mediated blockade of the Hippo signaling pathway. *FASEB J*. 2018 Jul 24:fj201800151RR.
36. Deng X, Zhang X, Li W, Feng RX, Li L, Yi GR, Zhang XN, Yin C, Yu HY, Zhang JP, Lu B, Hui L, Xie WF. Chronic Liver Injury Induces Conversion of Biliary Epithelial Cells into Hepatocytes. *Cell Stem Cell*. 2018 Jul 5;23(1):114-122.e3.
37. He W, Ye S, Zeng C, Xue S, Hu X, Zhang X, Gao S, Xiong Y, He X, Vivalda S, Li L, Wang Y, Ye Q. Hypothermic oxygenated perfusion (HOPE) attenuates ischemia/reperfusion injury in the liver through inhibition of the TXNIP/NLRP3 inflammasome pathway in a rat model of donation after cardiac death. *FASEB J*. 2018 Jun 5:fj201800028RR.
38. Wu D, Han R, Deng S, Liu T, Zhang T, Xie H, Xu Y. Protective Effects of Flagellin A N/C Against Radiation-Induced NLR Pyrin Domain Containing 3 Inflammasome-Dependent Pyroptosis in Intestinal Cells. *INT J RADIAT ONCOL*. 2018 May 1;101(1):107-117.
39. Gong Y, Zhang Z, Chang Z, Zhou H, Zhao R, He B. Inactivation of glycogen synthase kinase-3 α is required for mitochondria-mediated apoptotic germ cell phagocytosis in Sertoli cells. *AGING-US*. 2018 Nov 6;10(11):3104-3116.
40. Zhang YQ, Wen ZH, Wan K, Yuan D, Zeng X, Liang G, Zhu J, Xu B, Luo H. A novel synthesized 3', 5'-diprenylated chalcone mediates the proliferation of human leukemia cells by regulating apoptosis and autophagy pathways. *Biomed Pharmacother*. 2018 Oct;106:794-804.
41. Song Y, Li S, Geng W, Luo R, Liu W, Tu J, Wang K, Kang L, Yin H, Wu X, Gao Y, Zhang Y, Yang C. Sirtuin 3-dependent mitochondrial redox homeostasis protects against AGEs-induced intervertebral disc degeneration. *Redox Biol*. 2018 Oct;19:339-353.
42. Li R, Xin T, Li D, Wang C, Zhu H, Zhou H. Therapeutic effect of Sirtuin 3 on ameliorating nonalcoholic fatty liver disease: The role of the ERK-CREB pathway and Bnip3-mediated mitophagy. *Redox Biol*. 2018 Sep;18:229-243.
43. Zhou L, Yang L, Li YJ, Mei R, Yu HL, Gong Y, Du MY, Wang F. MicroRNA-128 Protects Dopamine Neurons from Apoptosis and Upregulates the Expression of Excitatory Amino Acid Transporter 4 in Parkinson's Disease by Binding to AXIN1. *CELL PHYSIOL BIOCHEM*. 2018;51(5):2275-2289.
44. Zhang J, An H, Ni K, Chen B, Li H, Li Y, Sheng G, Zhou C, Xie M, Chen S, Zhou T, Yang G, Chen X, Wu G, Jin S, Li M. Glutathione prevents chronic oscillating glucose intake-induced β -cell dedifferentiation and failure. *Cell Death Dis*. 2019 Apr 11;10(4):321.
45. Xu X, Wang D, Zheng C, Gao B, Fan J, Cheng P, Liu B, Yang L, Luo Z. Progerin accumulation in nucleus pulposus cells impairs mitochondrial function and induces intervertebral disc degeneration and therapeutic effects of sulforaphane. *Theranostics*. 2019 Apr 12;9(8):2252-2267.
46. Xu P, Shi X, Zhang X, Liu Q, Xie Y, Hong Y, Li J, Peng M, Liu X, Xu G. Overexpression of BRCA1 in Neural Stem Cells Enhances Cell Survival and Functional Recovery after Transplantation into Experimental Ischemic Stroke. *Oxid Med Cell Longev*. 2019 Apr 3;2019:8739730.
47. Li W, Lin J, Shi Z, Wen J, Li Y, Liu Z, Diao A. Clomiphene citrate induces nuclear translocation of the TFEB transcription factor and triggers apoptosis by enhancing lysosomal membrane permeabilization. *Biochem Pharmacol*. 2019 Apr;162:191-201.
48. Tao S, Sun Q, Cai L, Geng Y, Hua C, Ni Y, Zhao R. Caspase-1-dependent mechanism mediating the harmful impacts of the quorum-sensing molecule N-(3-oxo-dodecanoyl)-l-homoserine lactone on the intestinal cells. *J Cell Physiol*. 2019 Apr;234(4):3621-3633.
49. Liu Y, Xu J, Liu L, Tan J, Gao L, Wang J, Wang Y, Yan Z, Yu L. Amorphous Manganese Dioxide Coated Polydopamine Nanoparticles for Acid-Sensitive Magnetic Resonance Imaging-Guided Tumor Photothermal Therapy. *J BIOMED NANOTECHNOL*. 2019 Aug 1;15(8):1771-1780.
50. Duan PY, Ma Y, Li XN, Qu FZ, Ji L, Guo XY, Zhang WJ, Xiao F, Li L, Hu JS, Sun B, Wang G. Inhibition of RIPK1-dependent regulated acinar cell necrosis provides protection against acute pancreatitis via the RIPK1/NF- κ B/AQP8 pathway. *Exp Mol Med*. 2019 Aug 2;51(8):1-17.
51. Hua H, Cheng JW, Bu WB, Liu J, Ma WW, Ni N, Shi J, Zhou BR, Luo D. 5-aminolaevulinic acid-based photodynamic therapy inhibits ultraviolet B-induced skin photodamage. *Int J Biol Sci*. 2019 Aug 7;15(10):2100-2109.
52. Xiang Q, He X, Mu J, Mu H, Zhou D, Tang J, Xiao Q, Jiang Y, Ren G, Xiang T, Peng W. The phosphoinositide hydrolase phospholipase C delta1 inhibits epithelial-mesenchymal transition and is silenced in colorectal cancer. *J Cell Physiol*. 2019 Aug;234(8):13906-13916.
53. Huang X, Ni B, Xi Y, Chu X, Zhang R, You H. Protease-activated receptor 2 (PAR-2) antagonist AZ3451 as a novel therapeutic agent for osteoarthritis. *AGING-US*. 2019 Dec 16;11(24):12532-12545.
54. Dong G, Huang X, Jiang S, Ni L, Chen S. Simvastatin Mitigates Apoptosis and Transforming Growth Factor-Beta Upregulation in Stretch-Induced Endothelial Cells. *Oxid Med Cell Longev*. 2019 Dec 17;2019:6026051.
55. Hao SL, Yang WX. KIFC1 is essential for normal spermatogenesis and its depletion results in early germ cell apoptosis in the Kuruma shrimp, *Penaeus (Marsupenaeus) japonicus*. *AGING-US*. 2019 Dec 29;11(24):12773-12792.
56. Dang J, Ye H, Li Y, Liang Q, Li X, Yin L. Multivalency-assisted membrane-penetrating siRNA delivery sensitizes photothermal ablation via inhibition of tumor glycolysis metabolism. *Biomaterials*. 2019 Dec;223:119463.
57. Zhou WY, Zhang MM, Liu C, Kang Y, Wang JO, Yang XH. Long noncoding RNA LINC00473 drives the progression of pancreatic cancer via upregulating programmed death-ligand 1 by sponging microRNA-195-5p. *J Cell Physiol*. 2019 Dec;234(12):23176-23189.
58. Guo Y, Feng Y, Cui X, Wang Q, Pan X. Autophagy inhibition induces the repolarisation of tumour-associated macrophages and enhances chemosensitivity of laryngeal cancer cells to cisplatin in mice. *CANCER IMMUNOL IMMUN*. 2019 Dec;68(12):1909-1920.
59. Li X, Wu M, An D, Yuan H, Li Z, Song Y, Liu Z. Suppression of Tafazzin promotes thyroid cancer apoptosis via activating the JNK signaling pathway and enhancing INF2-mediated mitochondrial fission. *J Cell Physiol*. 2019 Feb 11
60. Xu K, Han B, Bai Y, Ma XY, Ji ZN, Xiong Y, Miao SK, Zhang YY, Zhou LM. MiR-451a suppressing BAP31 can inhibit proliferation and increase apoptosis through inducing ER stress in colorectal cancer. *Cell Death Dis*. 2019 Feb 15;10(3):152.
61. Wang L, Liu Y, Zhou Y, Wang J, Tu L, Sun Z, Wang X, Luo F. Zoledronic acid inhibits the growth of cancer stem cell derived from cervical cancer cell by attenuating their stemness phenotype and inducing apoptosis and cell cycle arrest through the Erk1/2 and Ak. *J EXP CLIN CANC RES*. 2019 Feb 21;38(1):93.
62. Zhang L, Li S, Wang R, Chen C, Ma W, Cai H. Cytokine augments the sorafenib-induced apoptosis in Huh7 liver cancer cell by inducing mitochondrial fragmentation and activating MAPK-JNK signalling pathway. *Biomed Pharmacother*. 2019 Feb;110:213-223.

63. Xu B, Zang SC, Li SZ, Guo JR, Wang JF, Wang D, Zhang LP, Yang HM, Lian S. HMGB1-mediated differential response on hippocampal neurotransmitter disorder and neuroinflammation in adolescent male and female mice following cold exposure. *Brain Behav Immun*. 2019 Feb;76:223-235.
64. Jin G, He R, Liu Q, Lin M, Dong Y, Li K, Tang BZ, Liu B, Xu F. Near-infrared light-regulated cancer theranostic nanoplatfrom based on aggregation-induced emission luminogen encapsulated upconversion nanoparticles. *Theranostics*. 2019 Jan 1;9(1):246-264.
65. Meng Y, Sun J, Wang X, Ma Y, Kong C, Zhang G, Dou H, Nan N, Shi M, Yu T, Piao H. The biological macromolecule *Nocardia rubra* cell-wall skeleton as an avenue for cell-based immunotherapy. *J Cell Physiol*. 2019 Jan 29
66. Deng J, Liu AD, Hou GQ, Zhang X, Ren K, Chen XZ, Li SSC, Wu YS, Cao X. N-acetylcysteine decreases malignant characteristics of glioblastoma cells by inhibiting Notch2 signaling. *J EXP CLIN CANC RES*. 2019 Jan 3;38(1):2.
67. Sun L, Li Q, Zhang L, Chai H, Yu L, Xu Z, Kang Y, Xue P. Stimuli responsive PEGylated bismuth selenide hollow nanocapsules for fluorescence/CT imaging and light-driven multimodal tumor therapy. *BIOMATER SCI-UK*. 2019 Jul 1;7(7):3025-3040.
68. Xu F, Yang J, Shang J, Lan F, Li M, Shi L, Shen L, Wang Y, Ge J. MicroRNA-302d promotes the proliferation of human pluripotent stem cell-derived cardiomyocytes by inhibiting LATS2 in the Hippo pathway. *CLIN SCI*. 2019 Jul 2;133(13):1387-1399.
69. Lin S, Yang L, Yao Y, Xu L, Xiang Y, Zhao H, Wang L, Zuo Z, Huang X, Zhao C. Flubendazole demonstrates valid antitumor effects by inhibiting STAT3 and activating autophagy. *J EXP CLIN CANC RES*. 2019 Jul 8;38(1):293.
70. Weng J, Tu M, Wang P, Zhou X, Wang C, Wan X, Zhou Z, Wang L, Zheng X, Li J, Wang Z, Wang Z, Chen C. Amiodarone induces cell proliferation and myofibroblast differentiation via ERK1/2 and p38 MAPK signaling in fibroblasts. *Biomed Pharmacother*. 2019 Jul;115:108889.
71. Li Y, Dang J, Liang Q, Yin L. Carbon monoxide (CO)-Strengthened cooperative bio-reductive anti-tumor therapy via mitochondrial exhaustion and hypoxia induction. *Biomaterials*. 2019 Jul;209:138-151.
72. Li H, Wang C, He T, Zhao T, Chen YY, Shen YL, Zhang X, Wang LL. Mitochondrial Transfer from Bone Marrow Mesenchymal Stem Cells to Motor Neurons in Spinal Cord Injury Rats via Gap Junction. *Theranostics*. 2019 Mar 17;9(7):2017-2035.
73. Yang R, Hou M, Gao Y, Zhang L, Xu Z, Kang Y, Xue P. Indocyanine green-modified hollow mesoporous Prussian blue nanoparticles loading doxorubicin for fluorescence-guided tri-modal combination therapy of cancer. *Nanoscale*. 2019 Mar 21;11(12):5717-5731.
74. Cao Y, Shi H, Sun Z, Wu J, Xia Y, Wang Y, Wu Y, Li X, Chen W, Wang A, Lu Y. Protective Effects of Magnesium Glycyrhizinate on Methotrexate-Induced Hepatotoxicity and Intestinal Toxicity May Be by Reducing COX-2. *Front Pharmacol*. 2019 Mar 25;10:119.
75. Yuzhu Cao, Hang Shi, Zhiguang Sun, Jiawei Wu, Yawen Xia, Yufei Wang, Yuanyuan Wu, Xiaoman Li, Wenxing Chen, Aiyun Wang, Yin Lu. Protective Effects of Magnesium Glycyrhizinate on Methotrexate-Induced Hepatotoxicity and Intestinal Toxicity May Be by Reducing COX-2. *Front Pharmacol*. 2019 Mar 25;10:119.
76. Du S, Li C, Lu Y, Lei X, Zhang Y, Li S, Liu F, Chen Y, Weng D, Chen J. Dioscin Alleviates Crystalline Silica-Induced Pulmonary Inflammation and Fibrosis through Promoting Alveolar Macrophage Autophagy. *Theranostics*. 2019 Mar 7;9(7):1878-1892.
77. Zhao Z, Qu W, Wang K, Chen S, Zhang L, Wu D, Chen Z. Bisphenol A inhibits mucin 2 secretion in intestinal goblet cells through mitochondrial dysfunction and oxidative stress. *Biomed Pharmacother*. 2019 Mar;111:901-908.
78. Jin M, Ji X, Zhang B, Sheng W, Wang R, Liu K. Synergistic effects of Pb and repeated heat pulse on developmental neurotoxicity in zebrafish. *ECOTOX ENVIRON SAFE*. 2019 May 15;172:460-470.
79. Zhao X, Fang J, Li S, Gaur U, Xing X, Wang H, Zheng W. Artemisinin Attenuated Hydrogen Peroxide (H₂O₂)-Induced Oxidative Injury in SH-SY5Y and Hippocampal Neurons via the Activation of AMPK Pathway. *Int J Mol Sci*. 2019 May 31;20(11). pii: E2680.
80. Liu S, Zhou H, Liu H, Ji H, Fei W, Luo E. Fluorine-contained hydroxyapatite suppresses bone resorption through inhibiting osteoclasts differentiation and function in vitro and in vivo. *CELL PROLIFERAT*. 2019 May;52(3):e12613.
81. Ge M, Liu C, Li L, Lan M, Yu Y, Gu L, Su Y, Zhang K, Zhang Y, Wang T, Liu C, Liu F, Li M, Xiong L, Wang K, He T, Dai Y, Zhao Y, Li N, Yu Z, Meng Q. miR-29a/b1 Inhibits Hair Follicle Stem Cell Lineage Progression by Spatiotemporally Suppressing WNT and BMP Signaling. *Cell Rep*. 2019 Nov 19;29(8):2489-2504.e4.
82. Li S, Zhang S, Chen J. c-Myc induced upregulation of long non-coding RNA SNHG16 enhances progression and carcinogenesis in oral squamous cell carcinoma. *Cancer Gene Ther*. 2019 Nov;26(11-12):400-410.
83. Yu M, Huang H, Dong S, Sha H, Wei W, Liu C. High mobility group box-1 mediates hippocampal inflammation and contributes to cognitive deficits in high-fat high-fructose diet-induced obese rats. *Brain Behav Immun*. 2019 Nov;82:167-177.
84. Zhang X, Hou G, Liu A, Xu H, Guan Y, Wu Y, Deng J, Cao X. Matrine inhibits the development and progression of ovarian cancer by repressing cancer associated phosphorylation signaling pathways. *Cell Death Dis*. 2019 Oct 10;10(10):770.
85. Ma JY, Feng X, Tian XY, Chen LN, Fan XY, Guo L, Li S, Yin S, Luo SM, Ou XH. The repair of endo/exogenous DNA double-strand breaks and its effects on meiotic chromosome segregation in oocytes. *Hum Mol Genet*. 2019 Oct 15;28(20):3422-3430.
86. Jian R, Yang M, Xu F. Lentiviral-mediated silencing of mast cell-expressed membrane protein 1 promotes angiogenesis of rats with cerebral ischemic stroke. *J Cell Biochem*. 2019 Oct;120(10):16786-16797.
87. Gong Y, Dong R, Gao X, Li J, Jiang L, Zheng J, Cui S, Ying M, Yang B, Cao J, He Q. Neohesperidin prevents colorectal tumorigenesis by altering the gut microbiota. *Pharmacol Res*. 2019 Oct;148:104460.
88. Huang J, Zhao Q, Li M, Duan Q, Zhao Y, Zhang H. The effects of endothelium-specific CYP2J2 overexpression on the attenuation of retinal ganglion cell apoptosis in a glaucoma rat model. *FASEB J*. 2019 Oct;33(10):11194-11209.
89. Xia Z, Wu S, Wei X, Liao Y, Yi P, Liu Y, Liu J, Liu J. Hypoxic ER stress suppresses β -catenin expression and promotes cooperation between the transcription factors XBP1 and HIF1 α for cell survival. *J Biol Chem*. 2019 Sep 13;294(37):13811-13821.
90. Li Z, Zhang X, Xue W, Zhang Y, Li C, Song Y, Mei M, Lu L, Wang Y, Zhou Z, Jin M, Bian Y, Zhang L, Wang X, Li L, Li X, Fu X, Sun Z, Wu J, Nan F, Chang Y, Yan J, Yu H, Feng X, Wang G, Zhang D, Fu X, Zhang Y, Young KH, Li W, Zhang M. Recurrent GNAQ mutation encoding T96S in natural killer/T cell lymphoma. *Nat Commun*. 2019 Sep 16;10(1):4209.
91. Shi X, Jiang X, Yuan B, Liu T, Tang Y, Che Y, Shi Y, Ai Q. LINC01093 Upregulation Protects against Alcoholic Hepatitis through Inhibition of NF- κ B Signaling Pathway. *MOL THER-NUCL ACIDS*. 2019 Sep 6;17:791-803.
92. Tan Z, Guo F, Huang Z, Xia Z, Liu J, Tao S, Li L, Feng Y, Du X, Ma L, Fu P. Pharmacological and genetic inhibition of fatty acid-binding protein 4 alleviated cisplatin-induced acute kidney injury. *J Cell Mol Med*. 2019 Sep;23(9):6260-6270.
93. Wang J, Zhu Y, Chen J, Yang Y, Zhu L, Zhao J, Yang Y, Cai X, Hu C, Rosell R, Sun X, Cao P. Identification of a novel PAK1 inhibitor to treat pancreatic cancer. *Acta Pharm Sin B*. 2020 Apr;10(4):603-614.
94. Kai-You Song, Xian-Zhao Zhang, Feng Li, Qing-Rong Ji. Silencing of ATP2B1-AS1 contributes to protection against myocardial infarction in mouse via blocking NFKBIA-mediated NF- κ B signalling pathway. *J Cell Mol Med*. 2020 Apr;24(8):4466-4479.
95. Kai Yuan, Jingtian Mei, Dandan Shao, Feng Zhou, Han Qiao, Yakun Liang, Kai Li, Tingting Tang. Cerium Oxide Nanoparticles Regulate Osteoclast Differentiation Bidirectionally by Modulating the Cellular Production of Reactive Oxygen Species. *INT J NANOMED*. 2020 Aug 25;15:6355-6372.
96. Tao Shu-Ya, Zhang Qiu-Yang, Li Jing-Jing, Yao Jin, Yan Biao. Suppression of pathological ocular neovascularization by a small molecule, SU1498. *Biomed Pharmacother*. 2020 Aug;128:110248.
97. Qin Shi, Baolan Sun, Di Wang, Yi Zhu, Xinxin Zhao, Xiaoqing Yang, Yuquan Zhang.

- Circ6401, a novel circular RNA, is implicated in repair of the damaged endometrium by Wharton's jelly-derived mesenchymal stem cells through regulation of the miR-29b-1-5p/RAP1B axis. *Stem Cell Res Ther.* 2020 Dec 1;11(1):520.
98. Jia Liu, Yang Zhan, Jiefu Wang, Junfeng Wang, Jiansheng Guo, Dalu Kong. Long noncoding RNA LINC01578 drives colon cancer metastasis through a positive feedback loop with the NF- κ B/YY1 axis. *Mol Oncol.* 2020 Dec;14(12):3211-3233.
 99. Dahong Yao, Jin Zhang, Jinhui Wang, Dabo Pan, Zhendan He. Discovery of novel ATAD2 bromodomain inhibitors that trigger apoptosis and autophagy in breast cells by structure-based virtual screening. *J ENZYM INHIB MED CH.* 2020 Dec;35(1):713-725.
 100. Pingping Xia, Fan Zhang, Yajing Yuan, Cheng Chen, Yan Huang, Longyan Li, E Wang, Qulian Guo, Zhi Ye. ALDH 2 conferred neuroprotection on cerebral ischemic injury by alleviating mitochondria-related apoptosis through JNK/caspase-3 signaling pathway. *Int J Biol Sci.* 2020 Feb 17;16(8):1303-1323.
 101. Yao D, Pan D, Zhen Y, Huang J, Wang J, Zhang J, He Z. Ferulin C triggers potent PAK1 and p21-mediated anti-tumor effects in breast cancer by inhibiting Tubulin polymerization in vitro and in vivo. *Pharmacol Res.* 2020 Feb;152:104605.
 102. Xiu He, Shi Chen, Chao Li, Jiaqi Ban, Yungeng Wei, Yangyang He, Fangwei Liu, Ying Chen, Jie Chen. Trehalose Alleviates Crystalline Silica-Induced Pulmonary Fibrosis via Activation of the TFEB-Mediated Autophagy-Lysosomal System in Alveolar Macrophages. *Cells.* 2020 Jan 4;9(1):122.
 103. Wei Li, Dan Li, Hongying Kuang, Xiaoling Feng, Wenxia Ai, Ying Wang, Simao Shi, Jing Chen, Ruifeng Fan. Berberine increases glucose uptake and intracellular ROS levels by promoting Sirtuin 3 ubiquitination. *Biomed Pharmacother.* 2020 Jan;121:109563.
 104. Deng Y, Wei Z, Huang M, Xu G, Wei W, Peng B, Nong S, Qin H. Long non-coding RNA F11-AS1 inhibits HBV-related hepatocellular carcinoma progression by regulating NR1I3 via binding to microRNA-211-5p. *J Cell Mol Med.* 2020 Jan;24(2):1848-1865.
 105. Yibin Deng, Zhongheng Wei, Meijin Huang, Guidan Xu, Wujun Wei, Bin Peng, Shunqiang Nong, Houji Qin. Long non-coding RNA F11-AS1 inhibits HBV-related hepatocellular carcinoma progression by regulating NR1I3 via binding to microRNA-211-5p. *J Cell Mol Med.* 2020 Jan;24(2):1848-1865.
 106. Chenxing Fu, Jie Zhan, Junqi Huai, Shaodan Ma, Minghui Li, Guoqin Chen, Minsheng Chen, Yanbin Cai, Caiwen Ou. Furin-instructed molecular self-assembly actuates endoplasmic reticulum stress-mediated apoptosis for cancer therapy. *Nanoscale.* 2020 Jun 11;12(22):12126-12132.
 107. Qing-Zhou Chen, Xin Wang, Fan Luo, Ning Li, Ni Zhu, Shuang Lu, Yu-Xing Zan, Chao-Jie Zhong, Mei-Rong Wang, Hai-Tao Hu, Yong-Zhen Zhang, Hai-Rong Xiong, Wei Hou. HTNV Sensitizes Host Toward TRAIL-Mediated Apoptosis-A Pivotal Anti-hantaviral Role of TRAIL. *Front Immunol.* 2020 Jun 19;11:1072.
 108. Xin Dang, Nai-Kei Wong, Yongli Xie, Vengatesen Thiagarajan, Fan Mao, Xiangyu Zhang, Yue Lin, Zhiming Xiang, Jun Li, Shu Xiao, Zohaib Noor, Yuanqiu He, Yang Zhang, Ziniu Yu. Autophagy Dually Induced by AMP Surplus and Oxidative Stress Enhances Hemocyte Survival and Bactericidal Capacity via AMPK Pathway in *Crassostrea hongkongensis*. *Front Cell Dev Biol.* 2020 Jun 3;8:411.
 109. Yuping Chen, Ni Jiao, Ming Jiang, Liping Liu, Yihui Zhu, Hongyan Wu, Jing Chen, Yingxue Fu, Qiu Du, Huiqin Xu, Jihu Sun. Loganin alleviates testicular damage and germ cell apoptosis induced by AGEs upon diabetes mellitus by suppressing the RAGE/p38MAPK/NF- κ B pathway. *J Cell Mol Med.* 2020 Jun;24(11):6083-6095.
 110. Jing Wang, Hong-Yan Li, Zhen-Bo Su. Stabilization of the histone acetyltransferase Tip60 by deubiquitinating enzyme USP7 stimulates the release of pro-inflammatory mediators in acute lung injury. *J MOL MED.* 2020 Jun;98(6):907-921.
 111. Wang H, Wei H, Wang J, Li L, Chen A, Li Z. MicroRNA-181d-5p-Containing Exosomes Derived from CAFs Promote EMT by Regulating CDX2/HOXA5 in Breast Cancer. *MOL THER-NUCL ACIDS.* 2020 Mar 6;19:654-667.
 112. Li B, Luan S, Chen J, Zhou Y, Wang T, Li Z, Fu Y, Zhai A, Bi C. The MSC-Derived Exosomal lncRNA H19 Promotes Wound Healing in Diabetic Foot Ulcers by Upregulating PTEN via MicroRNA-152-3p. *MOL THER-NUCL ACIDS.* 2020 Mar 6;19:814-826.
 113. Qian Xiang, Liang Kang, Juntan Wang, Zhiwei Liao, Yu Song, Kangcheng Zhao, Kun Wang, Cao Yang, Yukun Zhang. CircRNA-CIDN mitigated compression loading-induced damage in human nucleus pulposus cells via miR-34a-5p/SIRT1 axis. *EBioMedicine.* 2020 Mar;53:102679.
 114. Shaochen Pang, Yue Gao, Fengbang Wang, Yuanyuan Wang, Mengxi Cao, Wenjuan Zhang, Yong Liang, Maoyong Song, Guibin Jiang. Toxicity of silver nanoparticles on wound healing: A case study of zebrafish fin regeneration model. *Sci Total Environ.* 2020 May 15;717:137178.
 115. Xiaotong Yang, Xiaoying Ma, Ocholi Don, Yuanlin Song, Xiaoyan Chen, Jianwen Liu, Jieming Qu, Yun Feng. Mesenchymal stem cells combined with liraglutide relieve acute lung injury through apoptotic signaling restrained by PKA/ β -catenin. *Stem Cell Res Ther.* 2020 May 19;11(1):182.
 116. Wanxia Fang, Xiaofang Che, Guohui Li, Anhui Wang, Yizhe Wang, Xiaonan Shi, Kezuo Hou, Xiaojie Zhang, Xiujuan Qu, Yunpeng Liu. Sur-X, a novel peptide, kills colorectal cancer cells by targeting survivin-XIAP complex. *J EXP CLIN CANC RES.* 2020 May 7;39(1):82.
 117. Shuli Man, Li Yang, Hao Xiang, Guanyu Lu, Yijia Wang, Changxiao Liu, Wenyuan Gao. Antihypertensive and renal protective effect of Shunaoxin pill combined with captopril on spontaneous hypertension rats. *Biomed Pharmacother.* 2020 May;125:109977.
 118. Song Z, Zhang Y, Zhang H, Rajendran RS, Wang R, Hsiao CD, Li J, Xia Q, Liu K. Isoliquiritigenin triggers developmental toxicity and oxidative stress-mediated apoptosis in zebrafish embryos/larvae via Nrf2-HO1/JNK-ERK/mitochondrion pathway. *Chemosphere.* 2020 May;246:125727.
 119. Zhenzhen Song, Yun Zhang, Huazheng Zhang, R Samuel Rajendran, Rongchun Wang, Chung-Der Hsiao, Jianheng Li, Qing Xia, Kechun Liu. Isoliquiritigenin triggers developmental toxicity and oxidative stress-mediated apoptosis in zebrafish embryos/larvae via Nrf2-HO1/JNK-ERK/mitochondrion pathway. *Chemosphere.* 2020 May;246:125727.
 120. Wenlong Yang, Shuning Zhang, Tiantong Ou, Hao Jiang, Daile Jia, Zhiyong Qi, Yunzeng Zou, Juying Qian, Aijun Sun, Junbo Ge. Interleukin-11 regulates the fate of adipose-derived mesenchymal stem cells via STAT3 signalling pathways. *CELL PROLIFERAT.* 2020 May;53(5):e12771.
 121. Wen-Ning Xu, Run-Ze Yang, Huo-Liang Zheng, Lei-Sheng Jiang, Sheng-Dan Jiang. NDUFA4L2 Regulated by HIF-1 α Promotes Metastasis and Epithelial-Mesenchymal Transition of Osteosarcoma Cells Through Inhibiting ROS Production. *Front Cell Dev Biol.* 2020 Nov 20;8:515051.
 122. Feng Ding, Jinhua Liu, Xiaofei Zhang. microRNA-375 released from extracellular vesicles of bone marrow mesenchymal stem cells exerts anti-oncogenic effects against cervical cancer. *Stem Cell Res Ther.* 2020 Oct 27;11(1):455.
 123. Hongyan Li, Hong Zhang, Guomin Huang, Zhihui Dou, Yi Xie, Jing Si, Cuixia Di. Heavy ion radiation-induced DNA damage mediates apoptosis via the Rpl27a-Rpl5-MDM2-p53/E2F1 signaling pathway in mouse spermatogonia. *ECOTOX ENVIRON SAFE.* 2020 Sep 15;201:110831.
 124. Changzhen Ren, Kaiqiang Sun, Yanda Zhang, Yangxi Hu, Bowen Hu, Jian Zhao, Zhiqing He, Ru Ding, Weizhong Wang, Chun Liang. Sodium-Glucose CoTransporter-2 Inhibitor Empagliflozin Ameliorates Sunitinib-Induced Cardiac Dysfunction via Regulation of AMPK-mTOR Signaling Pathway-Mediated Autophagy. *Front Pharmacol.* 2021 Apr 29;12:664181.
 125. Jingyu Xu, Zonghao Tang, Youwu He, Shufang Cai, Beini Wang, Susu Zhang, Man Wu, Kai Qian, Kailun Zhang, Bo Chai, Guorong Chen, Ke Xu, Hao Ji, Jian Xiao, Yanqing Wu. D1-3-n-Butylphthalide Ameliorates Diabetic Nephropathy by Ameliorating Excessive Fibrosis and Podocyte Apoptosis. *Front Pharmacol.* 2021 Aug 23;12:628950.
 126. Huiqing Qu, Ye Zhu. SMPDL3B Predicts Poor Prognosis and Contributes to Development of Acute Myeloid Leukemia. *Front Mol Biosci.* 2021 Aug 24;8:695601.
 127. Yun Cao, Zhaowei Chen, Jijia Hu, Jun Feng, Zijiang Zhu, Yanqin Fan, Qiaoxuan Lin, Guohua Ding. Mfn2 Regulates High Glucose-Induced MAMs Dysfunction and

- Apoptosis in Podocytes via PERK Pathway. *Front Cell Dev Biol.* 2021 Dec 20;9:769213.
128. Changwei Yang, Wuqiong Wang, Pengxi Deng, Chen Li, Liangcai Zhao, Hongchang Gao. Fibroblast Growth Factor 21 Modulates Microglial Polarization That Attenuates Neurodegeneration in Mice and Cellular Models of Parkinson's Disease. *Front Aging Neurosci.* 2021 Dec 22;13:778527.
129. Yuxin Luo, Jinbo Guo, Wenxiu Jia, Mengyao Wu, Fengrong Yin, Guochao Niu, David Q Shih, Stephan R Targan, Xiaolan Zhang. TNF-Like Ligand 1 Aberrance Aggravates Nonalcoholic Steatohepatitis via M1 Macrophage Polarization. *Oxid Med Cell Longev.* 2021 Dec 31;2021:3877617.
130. Wenhao Ge, Tao Wang, Yang Zhao, Yunxia Yang, Qi Sun, Xiao Yang, Yan Gao, Xi Xu, Jianfa Zhang. Period1 mediates rhythmic metabolism of toxins by interacting with CYP2E1. *Cell Death Dis.* 2021 Jan 12;12(1):76.
131. Linlin Chen, Yuming Ma, Mengya Peng, Wenbo Chen, Huiqing Xia, Jingya Zhao, Yake Zhang, Zhuo Fan, Xiaoping Xing, Honglian Li. Analysis of Apoptosis-Related Genes Reveals that Apoptosis Functions in Conidiation and Pathogenesis of *Fusarium pseudograminearum*. *mSphere.* 2021 Jan 6;6(1):e01140-20.
132. Duo Zhang, Ran Zhang, Xiaoyuan Song, Karen Chang Yan, Haiyi Liang. Uniaxial Cyclic Stretching Promotes Chromatin Accessibility of Gene Loci Associated With Mesenchymal Stem Cells Morphogenesis and Osteogenesis. *Front Cell Dev Biol.* 2021 Jul 7;9:664545.
133. Xiaoyong Sun, Bin Wu, Liguogeng, Jinfeng Zhang, Guozheng Qin. Xiaokang Liuwei Dihuang decoction ameliorates the immune infertility of male rats induced by lipopolysaccharide through regulating the levels of sex hormones, reactive oxygen species, pro-apoptotic an. *Biomed Pharmacother.* 2021 Jul;139:111514.
134. Lei Mao, Qingcui Zeng, Wenjie Su, Menglong Song, Jiachen Li, Min Xie. Elevation of miR-146a Inhibits BTG2/BAX Expression to Ameliorate Postoperative Cognitive Dysfunction Following Probiotics (VSL#3) Treatment. *Mol Neurobiol.* 2021 Jul;58(7):3457-3470.
135. Shuo Wang, Chen Li, Panpan Sun, Jianli Shi, Xiaoyan Wu, Chang Liu, Zhe Peng, Hong Han, Shaojian Xu, Ying Yang, Yao Tian, Jiaxin Li, Hongbin He, Jun Li, Zhao Wang. PCV2 Triggers PK-15 Cell Apoptosis Through the PLC-IP3R-Ca²⁺ Signaling Pathway. *Front Microbiol.* 2021 Jun 15;12:674907.
136. Weilin Zheng, Jie Wang, Jiayi Wu, Tao Wang, Yangxue Huang, Xuefang Liang, Lixing Cao. Exploration of the Modulatory Property Mechanism of ELeng Capsule in the Treatment of Endometriosis Using Transcriptomics Combined With Systems Network Pharmacology. *Front Pharmacol.* 2021 Jun 18;12:674874.
137. Yuping Chen, Jing Chen, Anmei Shu, Liping Liu, Qin Wu, Juansong Wu, Siyuan Song, Weiping Fan, Yihui Zhu, Huiqin Xu, Jihu Sun, Liucai Yang. Combination of the Herbs *Radix Rehmanniae* and *Cornus Officinalis* Mitigated Testicular Damage From Diabetes Mellitus by Enhancing Glycolysis via the AGEs/RAGE/HIF-1 α Axis. *Front Pharmacol.* 2021 Jun 28;12:678300.
138. Huimin Huang, Mingxing Li, Yi Wang, Xiaoxiao Wu, Jing Shen, Zhangang Xiao, Yueshui Zhao, Fukuan Du, Yu Chen, Zhigui Wu, Huijiao Ji, Chunyuan Zhang, Jing Li, Qinglian Wen, Parham Jabbarzadeh Kaboli, Chi Hin Cho, Shengpeng Wang, Yitao Wang, Yisheng He, Xu Wu. Excessive Intake of Longan *Arillus* Alters gut Homeostasis and Aggravates Colitis in Mice. *Front Pharmacol.* 2021 Mar 26;12:640417.
139. Guangchan Jing, Huanyuan Wang, Fengwei Nan, Yuqin Liu, Mengren Zhang. Naofucong Ameliorates High Glucose Induced Hippocampal Neuron Injury Through Suppressing P2X7/NLRP1/Caspase-1 Pathway. *Front Pharmacol.* 2021 May 20;12:647116.
140. Lei Yu, Yuxi Chen, Shaohui Yuan, Yang Cao, Zhenggang Bi. Peiminine Induces G0/G1-Phase Arrest, Apoptosis, and Autophagy via the ROS/JNK Signaling Pathway in Human Osteosarcoma Cells in Vitro and in Vivo. *Front Pharmacol.* 2021 Nov 12;12:770846.
141. Zhi-Chao Sun, Zeyin Jiang, Xiaowen Xu, Meifeng Li, Qing Zeng, Ying Zhu, Shanghong Wang, Yuanyuan Li, Xiao-Li Tian, Chengyu Hu. Fish Paralog Proteins RNASEH-A and -B Enhance Type I Interferon Secretion and Promote Apoptosis. *Front Immunol.* 2021 Nov 22;12:762162.
142. Xuemei Yang, Jialing Sun, Bin Wen, Yu Wang, Mingjia Zhang, Weicong Chen, Wenting Zhao, Chunyu He, Xiaodan Zhong, Yang Liu, Tong Li, Haitao Sun, Songqi He. Biejiajian Pill Promotes the Infiltration of CD8⁺ T Cells in Hepatocellular Carcinoma by Regulating the Expression of CCL5. *Front Pharmacol.* 2021 Nov 26;12:771046.
143. Min Chen, Chungren, Bangqi Ren, Yajie Fang, Qiyu Li, Yuxiao Zeng, Yijian Li, Fang Chen, Baishijiao Bian, Yong Liu. Human Retinal Progenitor Cells Derived Small Extracellular Vesicles Delay Retinal Degeneration: A Paradigm for Cell-free Therapy. *Front Pharmacol.* 2021 Nov 29;12:748956.
144. Dingyi Lv, Minghao Luo, Zhe Cheng, Ruiyu Wang, Xiyang Yang, Yongzheng Guo, Longxiang Huang, Xiang Li, Bi Huang, Jian Shen, Suxin Luo, Jianghong Yan. Tubeimoside I Ameliorates Myocardial Ischemia-Reperfusion Injury through SIRT3-Dependent Regulation of Oxidative Stress and Apoptosis. *Oxid Med Cell Longev.* 2021 Nov 9;2021:5577019.
145. Changzheng Li, Jie Chen, Yishan Li, Binghuo Wu, Zhitao Ye, Xiaobin Tian, Yan Wei, Zechen Hao, Yuan Pan, Hongli Zhou, Keyue Yang, Zhiqiang Fu, Jingbo Xu, Yanan Lu. 6-Phosphogluconolactonase Promotes Hepatocellular Carcinogenesis by Activating Pentose Phosphate Pathway. *Front Cell Dev Biol.* 2021 Oct 26;9:753196.
146. Jin Liu, Lihong Pan, Wenxuan Hong, Siqin Chen, Peiyuan Bai, Wei Luo, Xiaolei Sun, Furong He, Xinlin Jia, Jialiang Cai, Yingjie Chen, Kai Hu, Zhenju Song, Junbo Ge, Aijun Sun. GPR174 knockdown enhances blood flow recovery in hindlimb ischemia mice model by upregulating AREG expression. *Nat Commun.* 2022 Dec 6;13(1):7519.
147. Lu Wang, Zihui Deng, Yan Sun, Yan Zhao, Yun Li, Mengmeng Yang, Rui Yuan, Yuyan Liu, Zhirong Qian, Feihu Zhou, Hongjun Kang. The Study on the Regulation of Th Cells by Mesenchymal Stem Cells Through the JAK-STAT Signaling Pathway to Protect Naturally Aged Sepsis Model Rats. *Front Immunol.* 2022 Feb 7;13:820685.
148. Wujiang Gao, Lu Chen, Li Lin, Meiling Yang, Taoqiong Li, Hong Wei, Chunli Sha, Jie Xing, Mengxue Zhang, Shijie Zhao, Qi Chen, Wenlin Xu, Yuefeng Li, Xiaolan Zhu. SIAH1 reverses chemoresistance in epithelial ovarian cancer via ubiquitination of YBX-1. *Oncogenesis.* 2022 Mar 10;11(1):13.
149. Hongjuan Yao, Wenping Song, Rui Cao, Cheng Ye, Li Zhang, Hebing Chen, Junting Wang, Yuchen Shi, Rui Li, Yi Li, Xiujun Liu, Xiaofei Zhou, Rongguang Shao, Liang Li. An EGFR/HER2-targeted conjugate sensitizes gemcitabine-sensitive and resistant pancreatic cancer through different SMAD4-mediated mechanisms. *Nat Commun.* 2022 Sep 20;13(1):5506.
150. Shengqiang Chen, Yanxia Zhu, Qingqing Xu, Qi Jiang, Danyang Chen, Ting Chen, Xishen Xu, Zhaokui Jin, Qianjun He. Photocatalytic glucose depletion and hydrogen generation for diabetic wound healing. *Nat Commun.* 2022 Sep 27;13(1):5684.
151. Guoqiang Sun, Yandong Zheng, Xiaolong Fu, Weiqi Zhang, Jie Ren, Shuai Ma, Shuhui Sun, Xiaojuan He, Qiaoran Wang, Zhejun Ji, Fang Cheng, Kaowen Yan, Ziyi Liu, Juan Carlos Izpisua Belmonte, Jing Qu, Si Wang, Renjie Chai, Guang-Hui Liu. Single-cell transcriptomic atlas of mouse cochlear aging. *Protein Cell.* 2023 Apr 13;14(3):180-201.

注：更多使用本产品的文献请参考产品网页。

Version 2024.03.12