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碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

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- Li J, Xia X, Ke Y, Nie H, Smith MA, Zhu X. Trichosanthin induced apoptosis in HL-60 cells via mitochondrial and endoplasmic reticulum stress signaling pathways. Biochim Biophys Acta. 2007 Aug; 1770(8):1169-80.
- Chen R, Jin R, Wu L, Ye X, Yang Y, Luo K, Wang W, Wu D, Ye X, Huang L, Huang T, Xiao G. Reticulon 3 attenuates the clearance of cytosolic prion aggregates via inhibiting autophagy. Autophagy. 2011 Feb; 7(2):205-16.
- Chen X, Ran ZH, Tong JL, Nie F, Zhu MM, Xu XT, Xiao SD. RNA interference (RNAi) of Ufd1 protein can sensitize a hydroxycamptothecin-resistant coloncancer cell line SW1116/HCPT to hydroxycamptothecin. J Dig Dis. 2011 Apr; 12(2):110-6.
- Ning B, Bai M, Shen W. Reduced glutathione protects human hepatocytes from palmitate-mediated injury by suppressingendoplasmic reticulum stress response. Hepatogastroenterology. 2011 Sep-Oct; 58(110-111):1670-9.
- Chen H, Wang DL, Liu YL. Poly (I:C) transfection induces mitochondrial-mediated apoptosis in cervical cancer. Mol Med Rep. 2016 Mar;13(3):2689-95.
- Xiao G,Zhou J,Lu X,Huang R,Zhang H. Excessive UDPG resulting from the mutation of UAP1 causes programmed cell death by triggering reactive oxygen species accumulation and caspase-like activity in rice. New Phytol. 2018 Jan;217(1):332-343.
- Huazhong Xie, Pengfei Qiang, Yao Wang, Fan Xia, Peiqing Liu, Min Li. Discovery and mechanism studies of a novel ATG4B inhibitor Ebselen by drug repurposing and its anti-colorectal cancer effects in mice. Cell Biosci. 2022 Dec 21;12(1):206.