碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301或800-8283301
 订货e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

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- Yang F, Sun X, Shen J, Yu LP, Liang JY, Zheng HQ, Wu ZD. A recombinant protein (rSj16) derived from Schistosoma japonicum induces cell cycle arrest and apoptosis of murine myeloid leukemia cells. Parasitol Res. 2013 Mar; 112(3):1261-72.
- Yi S, Wen L, He J, Wang Y, Zhao F, Zhao J, Zhao Z, Cui G, Chen Y. Deguelin, a selective silencer of the NPM1 mutant, potentiates apoptosis and induces differentiation in AML cells carrying the NPM1 mutation. Ann Hematol. 2014 Sep 23; 94(2):201-10.
- Xiao-Xiao Ni, Jing Nie, Qiu-You Xie, Rong-Hao Yu, Lei Su, Zhi-Feng Liu. Protective Effects of Hyperbaric Oxygen Therapy on Brain Injury by Regulating the Phosphorylation of Drp1 Through ROS/PKC Pathway in Heatstroke Rats Cell Mol Neurobiol. 2020 Nov; 40(8):1253-1269.
- Yanting Zhu, Jintuo Zhou, Peiguang Niu, Huajiao Chen, Daohua Shi. Cardamonin inhibits cell proliferation by caspase-mediated cleavage of Raptor N-S ARCH PHARMACOL. 2021 Apr; 394(4):809-817.
- Huazhong Xie, Pengfei Qiang, Yao Wang, Fan Xia, Peiqing Liu, Min Li. Discovery and mechanism studies of a novel ATG4B inhibitor Ebselen by drug repurposing and its anti-colorectal cancer effects in mice. Cell Biosci. 2022 Dec 21; 12(1):206.
- Baigang Zhang, Chenghui Huang, Qikun Lu, Hairong Liang, Jinliang Li, Dongmei Xu. Involvement of caspase in patulin-induced hepatotoxicity in vitro and in vivo. Toxicon. 2022 Jan 30; 206:64-73.
- Xin Liang, Ruyi Qian, Yiqun Ou, Dan Wang, Xianrong Lin, Chengliang Sun. Lipid peroxide-derived short-chain aldehydes promote programmed cell death in wheat roots under aluminum stress. J Hazard Mater. 2023 Feb 5; 443(Pt A):130142.

Version 2024.03.12