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1. Qingji Lyu, Qianying Yang, Jia Hao, Yuan Yue, Xiaodong Wang, Jianhui Tian, Lei An. A small proportion of X-linked genes contribute to X chromosome upregulation in early embryos via BRD4-mediated transcriptional activation. *Curr Biol.* 2022 Oct 24;32(20):4397-4410.e5.

2. Zong Miao, Wei Tian, Yangfan Ye, Wei Gu, Zhongyuan Bao, Lei Xu, Guangchi Sun, Chong Li, Yiming Tu, Honglu Chao, Sin Man Lam, Ning Liu, Jing Ji. Hsp90 induces Aclsl4-dependent glioma ferroptosis via dephosphorylating Ser637 at Drp1. *Cell Death Dis.* 2022 Jun 13;13(6):548.

3. Yun Cao, Zhaowei Chen, Jijia Hu, Jun Feng, Zijong Zhu, Yanqin Fan, Qiaoxuan Lin, Guohua Ding. Mfn2 Regulates High Glucose-Induced MAMs Dysfunction and Apoptosis in Podocytes via PERK Pathway. *Front Cell Dev Biol.* 2021 Dec 20;9:769213.

4. Jianmin Ling, Yanqing Wu, Xiaojing Zou, Yanmin Chang, Gang Li, Minghao Fang. (-)-Epicatechin Reduces Neuroinflammation, Protects Mitochondria Function, and Prevents Cognitive Impairment in Sepsis-Associated Encephalopathy. *Oxid Med Cell Longev.* 2022 May 24;2022:2657713.

5. Yong-Jie Chen, Yi Li, Xian Guo, Bo Huo, Yue Chen, Yi He, Rui Xiao, Xue-Hai Zhu, Ding-Sheng Jiang, Xiang Wei. Upregulation of IRF9 Contributes to Pulmonary Artery Smooth Muscle Cell Proliferation During Pulmonary Arterial Hypertension. *Front Pharmacol.* 2021 Dec 1;12:773235.

6. Xiaoxia Kong, Daopeng Lin, Liling Lu, Lidan Lin, Hongyu Zhang, Hailin Zhang. Apelin-13-Mediated AMPK ameliorates endothelial barrier dysfunction in acute lung injury mice via improvement of mitochondrial function and autophagy. *Int Immunopharmacol.* 2021 Dec;101(Pt B):108230.

7. Hulinyue Peng, Longtai You, Chunjing Yang, Kaixin Wang, Manting Liu, Dongge Yin, Yuchen Xu, Xiaoxv Dong, Xingbin Yin, Jian Ni. Ginsenoside Rb1 Attenuates Triptolide-Induced Cytotoxicity in HL-7702 Cells via the Activation of Keap1/Nrf2/ARE Pathway. *Front Pharmacol.* 2022 Jan 3;12:723784.

8. Hui Tian, Hao Yan, Yong Zhang, Qiaofen Fu, Chunyan Li, Juan He, Hui Li, Yong Zhou, Youguang Huang, Rongqing Li. Knockdown of mitochondrial threonyl-tRNA synthetase 2 inhibits lung adenocarcinoma cell proliferation and induces apoptosis. *Bioengineered.* 2022 Mar;13(3):5190-5204.

9. Pengpeng Duan, Jin Tan, Yuyang Miao, Qiang Zhang. PINK1/Parkin-Mediated Mitophagy Plays a Protective Role in Albumin Overload-Induced Renal Tubular Cell Injury. *Front Biosci (Landmark Ed).* 2022 Jun 7;27(6):184.

10. Ziyu Meng, Binchao Wang, Yiqiang Liu, Yejian Wan, Qianshi Liu, Huasheng Xu, Renchuan Liang, Ying Shi, Peng Tu, Hong Wu, Chuan Xu. Mitochondria-targeting polydopamine-coated nanodrugs for effective photothermal- and chemo-synergistic therapies against lung cancer. *Regen*

Biomater. 2022 Aug 1;9:rbac051.

11. Jiantao Guo, Yiping Yang. Parecoxib sodium alleviates ischemia reperfusion-induced pulmonary injury via inhibiting ERK/NF- $\kappa$ B and further activating the HIF-1 $\alpha$  pathway. *Immun Inflamm Dis*. 2022 Sep;10(9):e684.
12. Pengzhan Zhao, Yutian Wei, Guangchi Sun, Lei Xu, Tian Wang, Yufei Tian, Honglu Chao, Yiming Tu, Jing Ji. Fetuin-A alleviates neuroinflammation against traumatic brain injury-induced microglial necroptosis by regulating Nrf-2/HO-1 pathway. *J Neuroinflammation*. 2022 Nov 4;19(1):269.
13. Ya-Ling Yin, Huan-Huan Wang, Zi-Chen Gui, Shan Mi, Shuang Guo, Yue Wang, Qian-Qian Wang, Rui-Zhu Yue, Lai-Biao Lin, Jia-Xin Fan, Xue Zhang, Bing-Yan Mao, Tian-Heng Liu, Guang-Rui Wan, He-Qin Zhan, Mo-Li Zhu, Lin-Hua Jiang, Peng Li. Citronellal Attenuates Oxidative Stress-Induced Mitochondrial Damage through TRPM2/NHE1 Pathway and Effectively Inhibits Endothelial Dysfunction in Type 2 Diabetes Mellitus. *Antioxidants (Basel)*. 2022 Nov 14;11(11):2241.
14. Heng Sun, Ziying Sun, Xingquan Xu, Zhongyang Lv, Jiawei Li, Rui Wu, Yuxiang Fei, Guihua Tan, Zizheng Liu, Yuan Liu, Dongquan Shi. Blocking TRPV4 Ameliorates Osteoarthritis by Inhibiting M1 Macrophage Polarization via the ROS/NLRP3 Signaling Pathway. *Antioxidants (Basel)*. 2022 Nov 23;11(12):2315.
15. Xiuying Liang, Jingwen Sun, Haijing Guan, Qingyu Zhu, Wenjuan Yao. Angiotensin II Inhibits Adipogenic Differentiation and Promotes Mature Adipocyte Browning through the Corepressor CtBP1. *Biomedicines*. 2022 Dec 4;10(12):3131.
16. Qin Li, Xiangjian Meng, Qiang Hua. Circ ASAP2 decreased inflammation and ferroptosis in diabetic nephropathy through SOX2/SLC7A11 by miR-770-5p. *Acta Diabetol*. 2023 Jan;60(1):29-42.
17. Li-Wei Liu, Yu Xie, Guan-Qun Li, Tao Zhang, Yu-Hang Sui, Zhong-Jie Zhao, Yang-Yang Zhang, Wen-Bo Yang, Xing-Long Geng, Dong-Bo Xue, Hua Chen, Yong-Wei Wang, Tian-Qi Lu, Li-Ren Shang, Zhi-Bo Li, Le Li, Bei Sun. Gut microbiota-derived nicotinamide mononucleotide alleviates acute pancreatitis by activating pancreatic SIRT3 signalling. *Br J Pharmacol*. 2023 Mar;180(5):647-666.
18. Yu Kang, Rui Zhu, Shuang Li, Kun-Peng Qin, Hao Tang, Wen-Shan Shan, Zong-Sheng Yin. Erythropoietin inhibits ferroptosis and ameliorates neurological function after spinal cord injury. *Neural Regen Res*. 2023 Apr;18(4):881-888.

Version 2024.03.12