碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

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附表1. EndA- and EndA+ strains of E. coli.

EndA ⁻	EndA ⁺
BJ5183	BL21 (DE3)
DH1	CJ236
DH20	HB101
DH21	JM83
DH5α	JM101
JM103	JM110
JM105	LE392
JM106	MC1061
JM107	NM522 (all NM series are EndA ⁺)
JM108	NM554
JM109	P2392
MM294	PR700 (all PR series are EndA ⁺)
SK1590	Q358
SK1592	RR1
SK2267	TB1
SRB	TG1
TOP10	Y1088 (all Y10 series are EndA ⁺)
XL1-Blue	BMH 71-18
XLO	ES1301

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2. Jin Ding, Jun Liu, Cai-Fang Xue, Ying-Hui Li, Ya Zhao, Jun Chen, Yu-Xiao Huang, Zhong-Xiang Liu. TatPTD can introduce HBV targeted ribonuclease into hepatocytes. World Chin J Digestol. 2005 Apr 15;13(8):958-962.

3. Wu Y, Wang XJ, Liu XM, Zeng CY, Liang XQ. Recombinant PIG11 Gene Retroviral Vector Construction and Its High Expression in HepG2 Cells. Journal of Nanhua University(Medical Edition). 2007 Sep;35(5):653-59.

4. Wang F, Huang K, Yang L, Gong J, Tao Q, Li H, Zhao Y, Zeng S, Wu X, Stöckigt J, Li X, Qu J. Preparation of C-23 esterified silybin derivatives and evaluation of their lipid peroxidation inhibitory and DNA protective properties. Bioorg Med Chem. 2009 Sep 1;17(17):6380-9.

5. Jing RR, Cui M, Sun BL, Yu J, Wang HM. Tissue-specific expression profiling of receptor for advanced glycation end products and its soluble forms in esophageal and lung cancer. Genet Test Mol Biomarkers. 2010 Jun;14(3):355-61.

6. Liang QL, Wang BR, Li ZY, Chen GQ, Zhou Y. Basic research Construction of eukaryotic expression vector of TSLC1 gene. Arch Med Sci. 2011 Aug;7(4):579-85.

7. Zhou Y, Wang L, Yang F, Lin X, Zhang S, Zhao ZK. Determining the extremes of the cellular NAD(H) level by using an Escherichia coli NAD(+)-auxotrophic mutant. APPL ENVIRON MICROB. 2011 Sep;77(17):6133-40.

8. Yongjin J. Zhou, Fan Yang, Sufang Zhang, Haidong Tan and Zongbao K. Zhao. Efficient gene disruption in Saccharomyces cerevisiae using marker cassettes with long homologous arms prepared by the restriction-free cloning strategy. WORLD J MICROB BIOT. 2011 Dec; 27(12):2999-3003.

9. Yang F, Zhang S, Zhou YJ, Zhu Z, Lin X, Zhao ZK. Characterization of the mitochondrial NAD+ -dependent isocitrate dehydrogenase of the oleaginous yeast Rhodosporidium toruloides. APPL MICROBIOL BIOT. 2012 May;94(4):1095-105.

10. Liang QL, Wang BR, Li ZY, Chen GQ, Zhou Y. Effect of TSLC1 gene on growth and apoptosis in human esophageal carcinoma Eca109 cells. Arch Med Sci . 2012 Dec 20;8(6):987-92.

11. Yang W, Zhou Y, Ke Z. Production of dihydroxyacetone from glycerol by engineered Escherichia coli cells co-expressing gldA and nox genes. AFR J BIOTECHNOL . 2013 Jul;12(27): 4387-4392.

12. Wang X, Jiang Q, Wang W, Su L, Han Y, Wang C. Molecular mechanism of polypeptides from Chlamys farreri (PCF)'s anti-apoptotic effect in UVA-exposed HaCaT cells involves HSF1/HSP70, JNK, XO, iNOS and NO/ROS. J PHOTOCH PHOTOBIO B. 2014 Jan 5;130:47-56.

13. Tian W, Yin X, Wang L, Wang J, Zhu W, Cao J, Yang H. The key role of miR-21-regulated SOD2 in the medium-mediated bystander responses in human fibroblasts induced by α -irradiated keratinocytes. MUTAT RES-REV MUTAT. 2015 Oct;780:77-85.

14. Lu B, Chen L, Zhang Y, Shi Y, Zhou N. Quantitative analysis of G-protein-coupled receptor internalization using DnaE intein-based assay. METHOD CELL BIOL. 2016;132:293-318.

15. Xie Y, Wu B, Zhang XX, Yin J, Mao L, Hu M. Influences of graphene on microbial community and antibiotic resistance genes in mouse gut as determined by high-throughput sequencing. Chemosphere. 2016 Feb;144:1306-12.

16. Zhang B, Wang Y, Tan Z, Li Z, Jiao Z, Huang Q. Screening of Probiotic Activities of Lactobacilli Strains Isolated from Traditional Tibetan Qula, A Raw Yak Milk Cheese. ASIAN AUSTRAL J ANIM. 2016 Oct;29(10):1490-9.

17. Zeng K, Zhong B, Shen XL, Fang M, Lin BT, Ma DH. RNAi targeting Nogo Receptor enhanced survival and proliferation of murine retinal ganglion cells during N-methyl-D-aspartate-induced optic nerve crush. ONCOTARGET. 2017 Apr 21;8(39):65009-65021.

18. Xi LC, Ji YX, Yin D, Zhao ZX, Huang SC, Yu SL, Liu BY, Li HY. Effects of Dermatopontin gene silencing on apoptosis and proliferation of osteosarcoma MG-63 cells. Mol Med Rep. 2018 Jan;17(1):422-427.

19. Wang Y, Liu M, Yang P, Peng H. Peroxiredoxin 1 (PRDX1) Suppresses Progression and Metastasis of Osteosarcoma and Fibrosarcoma of Bone. MED SCI MONITOR. 2018 Jun 16;24:4113-4120.

20. Wang Z, Yu X, Niu X, Yang J, Tang Y, Hu J, Diao Y. Screening of the proteins interacting with NS1 of TMUV by yeast two-hybrid system and the identification of the function of the interacted protein. Infect Genet Evol. 2018 Sep;63:277-284.

21. Rong-Rong Ma, Jing Sun, Wen-Hong Fang, Ya-Ping Dong, Ji-Ming Ruan, Xian-Le Yang, Kun Hu. Identification of Carassius Auratus Gibelio Liver Cell Proteins Interacting With the GABA A Receptor γ 2 Subunit Using a Yeast Two-Hybrid System. Fish Physiol Biochem. 2019 Feb;45(1):199-208.

22. Min-Ke Shi, Yu-Long Xuan, Xiao-Feng He. FHL1 Overexpression as a Inhibitor of Lung Cancer Cell Invasion via Increasing RhoGDI β mRNA Expression. Cell J. 2022 May;24(5):239-244.

23. Changchang Meng, Shiyu Chen, Qi He, Junyi Tan, Jingxian Wu, Jing Zhao. IKZF3 modulates cerebral ischemia/reperfusion injury by inhibiting neuroinflammation. Int Immunopharmacol. 2023 Jan;114:109480.

24. Ju Yu, Junwei Zou, Xuan Liu, Ying Pan, Yuanyuan Mu, Shuyan Li, Juhua Wang, Fazhi Xu, Yong Wang. TaqMan-probe-based multiplex real-time RT-qPCR for simultaneous detection of GoAstV, GPV, and GoCV. Poult Sci. 2023 Feb;102(2):102396.

25. Yejuan Qiu, Zhongyang Qiu, Jun Xia, Xiaoyan Liu, Hanwen Zhang, Yuxiang Yang, Wenyi Hou, Xiangqian Li, Jianlong He. Co-expression of Xylose Transporter and Fructose-Bisphosphate Aldolase Enhances the Utilization of Xylose by Lactococcus lactis IO-1. Appl Biochem Biotechnol. 2023 Feb;195(2):816-831.

26. Haixia Du, Yanpeng Ma, Xiqiang Wang, Yong Zhang, Ling Zhu, Shuang Shi, Shuo Pan, Zhongwei Liu. Advanced glycation end products induce skeletal muscle atrophy and insulin resistance via activating ROS-mediated ER stress PERK/FOXO1 signaling. Am J Physiol Endocrinol Metab. 2023 Mar 1;324(3):E279-E287.

27. Xiaohan Jiang, Guoxun Li, Benzhi Zhu, Jingnan Zang, Tian Lan, Rui Jiang, Bing Wang. p20BAP31 induces cell apoptosis via both AIF caspase-independent and the ROS/JNK mitochondrial pathway in colorectal cancer. Cell Mol Biol Lett. 2023 Mar 28;28(1):25.

28. Yuliang Ren, Hui Wu, Miao Tan, Junjie Chen, Zhongqi Duan, Bingxin

Zhu, Xuzhi Ruan, Qingqing Yu, Shuzhen Li, Xuewen Liu, Ying Liu, Yuan Si. Acetylation of MOB1 mediates polyphyllin II-reduced lysosome biogenesis in breast cancer by promoting the cytoplasmic retention of the YAP/TFEB coactivator complex. *Phytomedicine*. 2024 Jan;122:155152.

Version 2024.03.12