碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

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1. Liu Z, Li T, Liu Y, Jia Z, Li Y, Zhang C, Chen P, Ma K, Affara N, Zhou C. WNT signaling promotes Nkx2.5 expression and early cardiomyogenesis via downregulation of Hdac1. *Biochim Biophys Acta*. 2009 Feb;1793(2):300-11.
2. Zhan X, Hu H, Ke C, Hu S, Wang D, Chen F. Isolation and characterization of eleven microsatellite loci in small abalone, *Haliotis diversicolor* Reeve. *Conserv Genet*. (2009) 10:1185–1187.
3. Shi G, Zhang Z, Feng D, Xu Y, Lu Y, Wang J, Jiang J, Zhang Z, Li X, Ning G. Selection of reference genes for quantitative real-time reverse transcription-polymerase chain reaction in concanavalin A-induced hepatitis model. *Anal Biochem*. 2010 Jun 1;401(1):81-90.
4. Yang W, Qin W, Hu Z, Suo Y, Zhao R, Ma X, Ma W, Wang T, Liang J, Tian J, Wang J. Comparison of Cerenkov Luminescence Imaging (CLI) and gamma camera imaging for visualization of let-7 expression in lung adenocarcinoma A549 Cells. *Nucl Med Biol*. 2012 Oct;39(7):948–953.
5. Li L, Cao W, Zheng W, Fan C, Chen T. Ruthenium complexes containing 2,6-bis(benzimidazolyl)pyridine derivatives induce cancer cell apoptosis by triggering DNA damage-mediated p53 phosphorylation. *Dalton Trans*. 2012 Nov 7;41(41):12766-72.
6. Sun Y, Mu F, Li C, Wang W, Luo M, Fu Y, Zu Y. Aspidin BB, a phloroglucinol derivative, induces cell cycle arrest and apoptosis in human ovarian HO-8910 cells. *Chem Biol Interact*. 2013 Jul 5;204(2):88-97.
7. Wang J, Yuan L, Xiao H, Xiao C, Wang Y, Liu X. Momordin Ic induces HepG2 cell apoptosis through MAPK and PI3K/Akt-mediated mitochondrial pathways. *Apoptosis*. 2013 Jun;18(6):751-65.
8. Luo Z, Yu L, Yang F, Zhao Z, Yu B, Lai H, Wong KH, Ngai SM, Zheng W, Chen T. Ruthenium polypyridyl complexes as inducer of ROS-mediated apoptosis in cancer cells by targeting thioredoxin reductase. *Metallomics*. 2014 Aug;6(8):1480-90.
9. Bo YX, Song XH, Wu K, Hu B, Sun BY, Liu ZJ, Fu JG. Characterization of interleukin-1 β as a proinflammatory cytokine in grass carp (*Ctenopharyngodon idella*). *Fish Shellfish Immunol*. 2015 Oct;46(2):584-95.
10. Qu W, Wang M, Wu Y, Lv Y, Wang Q, Xu R. Calcium-ion-modulated ceramic hydroxyapatite resin for the scalable purification of recombinant Adeno-Associated Virus serotype 9. *J Chromatogr B Analyt Technol Biomed Life Sci*. 2015 May 15;990:15-22.
11. Yang Q, Cai J, Sun S, Kang X, Guo J, Zhu Y, Yan L, Jing X, Wang Z. Polymer nanoparticle delivery of dichloroacetate and DACH-Pt to enhance antitumor efficacy and lower systemic toxicity. *Biomater Sci*. 2016 Apr;4(4):661-9.
12. Hou X, Lin L, Xing W, Yang Y, Duan X, Li Q, Gao X, Lin Y. Spleen tyrosine kinase regulates mammary epithelial cell proliferation in mammary glands of dairy cows. *J Dairy Sci*. 2016 May;99(5):3858-68.
13. Li Y, Shang Q, Li P, Cao J, Zhu L, Jager MJ, Fan X, Ge S, Jia R. Characterization of a conjunctival melanoma cell line CM-AS16, newly-established from a metastatic Han Chinese. *Exp Eye Res*. 2018 Apr 10. pii: S0014-4835(18)30054-X
14. Li Y, Shang Q, Li P, Cao J, Zhu L, Jager MJ, Fan X, Ge S, Jia R. Characterization of a conjunctival melanoma cell line CM-AS16, newly-established from a metastatic Han Chinese patient. *Exp Eye Res*. 2018 Aug;173:51-63.
15. Haiqin Song, Feifei Quan, Zhiqiang Yu, Minhua Zheng, Yan Ma, Haihua Xiao, Fang Ding. Carboplatin Prodrug Conjugated Fe₃O₄ Nanoparticles for Magnetically Targeted Drug Delivery in Ovarian Cancer Cells. *J Mater Chem B*. 2019 Jan 21;7(3):433-442.;doi: 10.1039/c8tb02574f
16. Ha S, Zhou H, Gautam M, Song Y, Wang C. Reduced ribosomal RNA expression and unchanged ribosomal DNA promoter methylation in oral squamous cell carcinoma. *Mol Genet Genomic Med*. 2019 Jul;7(7):e00783.
17. Xin H, Wang Z, Wu S, Wang P, Tao X, Xu C, You L. Calcified decellularized arterial scaffolds impact vascular smooth muscle cell transformation via downregulating α -SMA expression and upregulating OPN expression. *Exp Ther Med*. 2019 Jul;18(1):705-710
18. Li Yin, Jie Chen, Chengxian Ma, Shuai Pei, Mingyu Du, Yufeng Zhang, Yong Feng, Rong Yin, Xiuhua Bian, Xia He, Jifeng Feng. Hsa_circ_0046263 functions as a ceRNA to promote nasopharyngeal carcinoma progression by upregulating IGFBP3. *Cell Death Dis*. 2020 Jul 23;11(7):562.;doi: 10.1038/s41419-020-02785-3
19. Ying Chen, Jinjun Qian, Pinggang Ding, Wang Wang, Xinying Li, Xiaozhu Tang, Chao Tang, Ye Yang, Chunyan Gu. Elevated SFXN2 limits mitochondrial autophagy and increases iron-mediated energy production to promote multiple myeloma cell proliferation. *Cell Death Dis*. 2022 Sep 26;13(9):822.
20. Minghui Zhang, Xueying Li, Xu Cui, Rui Li, Zonghua Ma, Xuejun Gao. Selenomethionine promotes ANXA2 phosphorylation for proliferation and protein synthesis of myoblasts and skeletal muscle growth. *J Nutr Biochem*. 2023 May;115:109277.

Version 2024.03.12