

QuickBlock™ Western封闭液

产品编号	产品名称	包装
P0252-100ml	QuickBlock™ Western封闭液	100ml
P0252-500ml	QuickBlock™ Western封闭液	500ml

产品简介:

- 碧云天生产的QuickBlock™ Western封闭液(QuickBlock™ Blocking Buffer for Western Blot)是最新一代的快速高效Western封闭液, 总体效果显著优于传统的基于BSA(牛血清白蛋白)、脱脂奶粉、酪蛋白(Casein)等的封闭液及国外同类产品, 主要用于Western blot (WB)实验中PVDF膜或硝酸纤维素膜(nitrocellulose membrane, NC膜)的封闭, 也可以用于Western实验中一抗或二抗的稀释。
- QuickBlock™ Western封闭液**快速高效**。封闭时间通常仅需5-15分钟, 并且和BSA、脱脂奶粉、酪蛋白等传统的Western封闭液以及国外同类的快速封闭液相比, 显示出更强的信噪比(参考图1)。
- QuickBlock™ Western封闭液封闭后**背景极低**。本封闭液不含血清和白蛋白, 确保极高的信噪比。
- QuickBlock™ Western封闭液**兼容性好**, 兼容辣根过氧化物酶(Horseradish peroxidase, HRP)、碱性磷酸酶(Alkaline phosphatase, AP)和生物素标记的二抗。本产品中添加了不影响HRP和AP活性的防腐剂, 不会干扰HRP或AP标记二抗的检测。同时本产品不含生物素, 不会干扰基于生物素的检测。
- QuickBlock™ Western封闭液**使用便捷**。本产品无需添加任何额外的试剂, 可以直接用于印迹膜的封闭。
- 本产品与BSA及国外同类产品的封闭效果对比参见图1。在相同样品和实验条件下, 仅封闭液及封闭时间存在如下图所示的差异时, 碧云天的QuickBlock™ Western封闭液封闭后的整体背景明显低于BSA封闭后的背景, 而且目的条带亮度明显高于国外同类品牌产品。




图1. QuickBlock™ Western封闭液与BSA及国外同类产品的封闭效果比较。每组实验从左到右依次为: 5 μ l蛋白Marker, 2.5 μ g蛋白量的HeLa细胞裂解液, 5 μ g蛋白量的HeLa细胞裂解液。请注意对于Antibody A, QuickBlock™封闭液比BSA和国外同类产品具有更低的背景, 并且和国外同类产品相比信号明显更强。对于Antibody B, QuickBlock™封闭液比BSA具有更低的背景, 和国外同类产品相比信号明显更强。实际实验结果会因样品、抗体、实验条件等的不同而存在差异, 图中数据仅供参考。

- 关于不同封闭液的选择, 请参考碧云天的相关网页: <http://www.beyotime.com/support/blocking-buffer.htm>。
- 按照每张膜封闭需要5-10ml QuickBlock™ Western封闭液计算, 一个100ml包装的本产品可以封闭10-20张膜, 一个500ml包装的本产品可以封闭50-100张膜。

包装清单:

产品编号	产品名称	包装
P0252-100ml	QuickBlock™ Western封闭液	100ml
P0252-500ml	QuickBlock™ Western封闭液	500ml
—	说明书	1份

保存条件:

4°C保存, 一年有效。长期不使用可以-20°C保存。

注意事项:

- 通常本产品用于PVDF膜及NC膜时的封闭时间为5-15分钟。对于一些背景非常高的抗体，可以尝试将封闭时间延长为30-60分钟。此外，如有特殊需要，也完全可以4°C封闭过夜。
- 由于没有任何一种封闭液是适用于所有实验体系的，因此对于一些特殊的实验，可能需要根据具体情况考虑使用其它更合适的封闭液。
- 取放PVDF膜和NC膜应使用平头镊子，并仅轻轻夹取其边角，操作过程须避免膜表面产生划痕、折痕或压痕等痕迹。
- PVDF膜一经浸润和活化，需一直保持湿润，根据Western进行到的具体步骤可放置于western转膜液或洗涤液等适当溶液中，否则可能会产生难以封闭的异常背景。
- 为进一步提高信噪比，推荐同时使用QuickBlock™ Western一抗稀释液(P0256)和QuickBlock™ Western二抗稀释液(P0258)进行一抗及二抗的稀释。
- 本产品仅限于专业人员的科学研究用，不得用于临床诊断或治疗，不得用于食品或药品，不得存放于普通住宅内。
- 为了您的安全和健康，请穿实验服并戴一次性手套操作。

使用说明：

1. 完成转膜后，用Western洗涤液洗涤印迹膜1-2分钟。
2. 根据膜的大小，在平皿或者其它适当容器中倒入一定体积的QuickBlock™ Western封闭液，确保封闭液后续能充分覆盖膜即可。对于常规的western，一张约6.6×8.5cm的膜推荐使用约10ml左右的封闭液。
3. 用平头镊夹住膜的一角，将膜放置在QuickBlock™ Western封闭液中，使封闭液完全浸没膜，置于水平摇床上封闭约10分钟(通常5-15分钟均可；经多种抗体的测试封闭10分钟的效果很多时候会显著优于常规的BSA封闭1小时的效果)。
4. 封闭后的膜即可用于一抗孵育等后续实验。详细的Western操作可以参考如下的相关网页：
<http://www.beyotime.com/support/western.htm>。

相关产品：

产品编号	产品名称	包装
P0023A-100ml	Western一抗稀释液	100ml
P0023A-500ml	Western一抗稀释液	500ml
P0023B-100ml	Western封闭液	100ml
P0023B-500ml	Western封闭液	500ml
P0023D-100ml	Western二抗稀释液	100ml
P0023D-500ml	Western二抗稀释液	500ml
P0220	QuickBlock™封闭液(PBS)	100ml
P0222	QuickBlock™封闭液(PBSTw)	100ml
P0226	QuickBlock™封闭液(PBSTx)	100ml
P0228	QuickBlock™封闭液(TBS)	100ml
P0231	QuickBlock™封闭液(TBSTw)	100ml
P0233	QuickBlock™封闭液(TBSTx)	100ml
P0235	QuickBlock™封闭液(10X)	100ml
P0252-100ml	QuickBlock™ Western封闭液	100ml
P0252-500ml	QuickBlock™ Western封闭液	500ml
P0256-100ml	QuickBlock™ Western一抗稀释液	100ml
P0256-500ml	QuickBlock™ Western一抗稀释液	500ml
P0258-100ml	QuickBlock™ Western二抗稀释液	100ml
P0258-500ml	QuickBlock™ Western二抗稀释液	500ml

使用本产品的文献：

1. Wang J, Feng S, Li M, Liu Y, Yan J, Tang Y, Du D, Chen F Increased Expression of Kv10.2 in the Hippocampus Attenuates Valproic Acid-Induced Autism-Like Behaviors in Rats. *Neurochem Res* 2019 Dec
2. Wang X, Meng K, He Y, Wang H, Zhang Y, Quan F Melatonin Stimulates STAR Expression and Progesterone Production via Activation of the PI3K/AKT Pathway in Bovine Theca Cells. *Int J Biol Sci* 2019 Jan 1
3. Liao H, Liu S, Wang H, Su H, Liu Z Enhanced antifungal activity of bovine lactoferrin-producing probiotic *Lactobacillus casei* in the murine model of vulvovaginal candidiasis. *BMC Microbiol* 2019 Jan 8
4. Li M, Wang J, Liu Y, Wu J, Li F, Feng S, Du D, Chen F Central changes in the Kv10.2 potassium channel in stress-induced hypertension rats. *Neuroreport* 2019 Jun 12
5. Pan Z, Fan Z, Ma J, Liu H, Shen L, He B, Zhang M Profiling and functional characterization of circulation lncRNAs that are associated with coronary atherosclerotic plaque stability. *Am J Transl Res* 2019 Jun 15
6. Liu Z, Tian X, Liu W, Xian Y, Chen W, Chen H, Zhou R Development of two antigen-binding fragments to a conserved linear epitope of human adenovirus and their application in immunofluorescence. *PLoS One* 2019 Jun 26
7. Xiao Q, Zeng JH, Zhou H, Qiu QH, Ke B, Deng L, Hu ZM, Roh J, Dai M. Expression and effects of leukemia inhibitory factor on nucleus pulposus degeneration. *Mol Med Rep* 2019 Mar
8. Ran K, Yang Z, Zhao Y, Wang X. Transmural pressure drives proliferation of human arterial smooth muscle cells via mechanism associated with NADPH oxidase and Survivin. *Microvasc Res* 2019 Nov
9. Tan Y, Jin Y, Wang Q, Huang J, Wu X, Ren Z Perilipin 5 Protects against Cellular

- Oxidative Stress by Enhancing Mitochondrial Function in HepG2 Cells. *Cells* 2019 Oct 11
10. Deng Y, Han Q, Mei S, Li H, Yang F, Wang J, Ge S, Jing X, Xu H, Zhang T Cyclin-dependent kinase subunit 2 overexpression promotes tumor progression and predicts poor prognosis in uterine leiomyosarcoma. *Oncol Lett* 2019 Sep
 11. Yang X, Zhang C, Tie H, Luo J, Wang Y, Wu Q miR-760 exerts an antioncogenic effect in esophageal squamous cell carcinoma by negatively driving fat metabolism via targeting c-Myc. *J Cell Biochem* 2020 Apr
 12. Chufan Yan, Jianguye Zhao, Yu Qin, Fangkun Zhao, Liyang Ji, Jinsong Zhang Overexpression of ATG4a promotes autophagy and proliferation, and inhibits apoptosis in lens epithelial cells via the AMPK and Akt pathways *Mol Med Rep* 2020 Aug;22(2):1295-1302.
 13. Yongjie Xiong, Qirun Yin, Jing Li, Shaojun He Oxidative Stress and Endoplasmic Reticulum Stress Are Involved in the Protective Effect of Alpha Lipoic Acid Against Heat Damage in Chicken Testes *Animals (Basel)* 2020 Feb 27;10(3):384.
 14. Feng Min, Xian Jie Jia, Qin Gao, Fang Niu, Zhi Yuan Hu, Ya Ling Han, Hong Jie Shi, Ying Yu Remote ischemic post-conditioning protects against myocardial ischemia/reperfusion injury by inhibiting the Rho-kinase signaling pathway *Exp Ther Med* 2020 Jan;19(1):99-106.
 15. Xiangru Zheng, Jinlai Wei, Wenjun Li, Xiaoli Li, Wuyi Wang, Jinbao Guo, Zhongxue Fu PRDX2 removal inhibits the cell cycle and autophagy in colorectal cancer cells *Aging (Albany NY)* 2020 Jul 20;12(16):16390-16409.
 16. Hui Hou, Ruirui Yang, Xiaohong Liu, Xiaolong Wu, Sulin Zhang, Kaixian Chen, Mingyue Zheng Discovery of triazoloxinoxaline as novel STING agonists via structure-based virtual screening *Bioorg Chem* 2020 Jul;100:103958.
 17. Xin Dang, Nai-Kei Wong, Yongli Xie, Vengatesen Thiyagarajan, Fan Mao, Xiangyu Zhang, Yue Lin, Zhiming Xiang, Jun Li, Shu Xiao, Zohaib Noor, Yuanqiu He, Yang Zhang, Ziniu Yu Autophagy Dually Induced by AMP Surplus and Oxidative Stress Enhances Hemocyte Survival and Bactericidal Capacity via AMPK Pathway in *Crassostrea hongkongensis* *Front Cell Dev Biol* 2020 Jun 3;8:411.
 18. Jiangbei Yuan, Zihan Zheng, Liting Wang, Haiying Ran, Xiangyu Tang, Xiaodong Xie, Fei Li, Fang Liu, Xiaoyang Wang, Jiale Zhang, Junying Zhang, Yi Huang, Xuefeng Xia, Ying Wan The Dynll1-Cox4i1 Complex Regulates Intracellular Pathogen Clearance via Release of Mitochondrial Reactive Oxygen Species *Infect Immun* 2020 Mar 23;88(4):e00738-19.
 19. Feng-Zhi Xin, Ze-Hua Zhao, Rui-Nan Zhang, Qin Pan, Zi-Zhen Gong, Chao Sun, Jian-Gao Fan Folic acid attenuates high-fat diet-induced steatohepatitis via deacetylase SIRT1-dependent restoration of PPAR α *World J Gastroenterol* 2020 May 14;26(18):2203-2220.
 20. Niya Long, Shuo Peng, Liangzhao Chu, Jun Jia, Minghao Dong, Jian Liu Paclitaxel inhibits the migration of CD133+ U251 malignant glioma cells by reducing the expression of glycolytic enzymes *Exp Ther Med* 2020 Nov;20(5):72.
 21. Jiechun Huang, Fangrui Wang, Xiaotian Sun, Xianglin Chu, Rongrong Jiang, Yiqing Wang, Liewen Pang Myocardial infarction cardiomyocytes-derived exosomal miR-328-3p promote apoptosis via Caspase signaling *Am J Transl Res* 2021 Apr 15;13(4):2365-2378.
 22. Chao Zhou, Mengyu Liu, Xiang Mei, Qian Li, Wenjuan Zhang, Ping Deng, Zhixin He, Yu Xi, Tong Tong, Hui Feng Pi, Yonghui Lu, Chunhai Chen, Lei Zhang, Zhengping Yu, Zhou Zhou, Mindi He Histone hypoacetylation contributes to neurotoxicity induced by chronic nickel exposure in vivo and in vitro *Sci Total Environ* 2021 Aug 20;783:147014.
 23. Chuan-Suo Zhang, Qian Han, Zhao-Wei Song, Hong-Yan Jia, Tian-Peng Shao, Yan-Peng Chen Hydrogen gas post-conditioning attenuates early neuronal pyroptosis in a rat model of subarachnoid hemorrhage through the mitoKATP signaling pathway *Exp Ther Med* 2021 Aug;22(2):836.
 24. Qian Zhou, Hong Cao, Xiaoyi Hang, Huamin Liang, Miaomiao Zhu, Yixian Fan, Jiawei Shi, Nianguo Dong, Ximiao He Midkine Prevents Calcification of Aortic Valve Interstitial Cells via Intercellular Crosstalk *Front Cell Dev Biol* 2021 Dec 15;9:794058.
 25. Xin Wang, Wenting Zhang, Jingtao Na, Yanping Huo, Yacheng Wang, Ketong Liu Spironolactone Inhibits Cardiomyocyte Hypertrophy by Regulating the Ca²⁺/Calcineurin/p-NFATc3 Pathway *J Healthc Eng* 2021 Dec 16;2021:3843830.
 26. Yuhan Shu, Xin Sun, Guiqin Ye, Mengting Xu, Zhipan Wu, Caixia Wu, Shouxin Li, Jingkui Tian, Haote Han, Jianbin Zhang DHOK Exerts Anti-Cancer Effect Through Autophagy Inhibition in Colorectal Cancer *Front Cell Dev Biol* 2021 Dec 17;9:760022.
 27. Weibo Dai, Xinyi Zhan, Weijie Peng, Xin Liu, Weiwen Peng, Quanxi Mei, Xianjing Hu Ficus pandurata Hance Inhibits Ulcerative Colitis and Colitis-Associated Secondary Liver Damage of Mice by Enhancing Antioxidation Activity *Oxid Med Cell Longev* 2021 Dec 18;2021:2617881.
 28. Bin Wang, Ping Yu, Wei Lin, Zhaohui Zhai MicroRNA-21-5p Reduces Hypoxia/Reoxygenation-Induced Neuronal Cell Damage through Negative Regulation of CPEB3 *Anal Cell Pathol (Amst)* 2021 Dec 2;2021:5543212.
 29. Wenting Zhang, Xin Wang, Jing Li, Mingyuan Xu, Xiaolu Ren, Huiying Liu, Lu Xu, Jun Shao Astragaloside IV Reduces OxLDL-Induced BNP Overexpression by Regulating HDAC *J Healthc Eng* 2021 Dec 3;2021:3433615.
 30. Mengmeng Li, Hongjian Lu, Xueyan Wang, Chengwei Duan, Xiangyang Zhu, Yi Zhang, Xin Ge, Feng Ji, Xueqin Wang, Jianbin Su, Dongmei Zhang Pyruvate kinase M2 (PKM2) interacts with activating transcription factor 2 (ATF2) to bridge glycolysis and pyroptosis in microglia *Mol Immunol* 2021 Dec;140:250-266.
 31. Li-Tian Ma, Yang Bai, Jie Li, Yu Qiao, Yang Liu, Jin Zheng Elemene Emulsion Injection Administration Reduces Neuropathic Pain by Inhibiting Astrocytic NDRG2 Expression within Spinal Dorsal Horn *Chin J Integr Med* 2021 Dec;27(12):912-918.
 32. Yi Wu, Heming Cui, Yuying Zhang, Ping Yu, Yuangeng Li, Dan Wu, Yan Xue, Wenwen Fu Inonotus obliquus extract alleviates myocardial ischemia/reperfusion injury by suppressing endoplasmic reticulum stress *Mol Med Rep* 2021 Jan;23(1):77.
 33. Bin Zeng, Zhiwei Sun, Qiting Zhao, Doudou Liu, Hao Chen, Xiaoshuang Li, H Rosie Xing, Jianyu Wang SEC23A Inhibit Melanoma Metastatic through Secretory PF4 Cooperation with SPARC to Inhibit MAPK Signaling Pathway *Int J Biol Sci* 2021 Jul 13;17(12):3000-3012.
 34. Li-Min Zhang, Dong-Xue Zhang, Wei-Chao Zheng, Jin-Shu Hu, Lan Fu, Yan Li, Yue Xin, Xu-Peng Wang CORM-3 exerts a neuroprotective effect in a rodent model of traumatic brain injury via the bidirectional gut-brain interactions *Exp Neurol* 2021 Jul;341:113683.
 35. Keke Wu, Biao Li, Qiuzhen Lin, Wanghan Xu, Wanyun Zuo, Jiayi Li, Na Liu, Tao Tu, Baojian Zhang, Yichao Xiao, Qiming Liu Nicotinamide mononucleotide attenuates isoproterenol-induced cardiac fibrosis by regulating oxidative stress and Smad3 acetylation *Life Sci* 2021 Jun 1;274:119299.
 36. Wenjin Chen, Weiqiang Jia, Cuiying Wu, Lihua Chen, Kai Sun, Ji Wang, Boyun Ding, Ning Liu, Ruxiang Xu The Neurogenic Compound P7C3 Regulates the Aerobic Glycolysis by Targeting Phosphoglycerate Kinase 1 in Glioma *Front Oncol* 2021 Jun 18;11:644492.
 37. Jingfei Zheng, Xuehe Li, Weili Yang, Fang Zhang Dihydroartemisinin regulates apoptosis, migration, and invasion of ovarian cancer cells via mediating RECK *J Pharmacol Sci* 2021 Jun;146(2):71-81.
 38. Yanjie Tan, Yi Jin, Sheng Wang, Jianhua Cao, Zhuqing Ren The RNA surveillance factor UPF1 regulates the migration and adhesion of porcine skeletal muscle satellite cells *J Muscle Res Cell Motil* 2021 Jun;42(2):203-217.
 39. Dongmei Ye, Yuxuan Li, Heliang Zhang, Zhiwei Zhou, Yujie Tang, Peng Wu, Qiang Zhao, Zhiwei Zhang Silencing PRSS1 suppresses the growth and proliferation

- of gastric carcinoma cells via the ERK pathway *Int J Biol Sci* 2021 Mar 1;17(4):957-971.
40. Donghua Du, Wenting Lv, Rina Su, Chunwei Yu, Xiaoxia Jing, Nuwenqimuge Bai, Surong Hasi Hydrolyzed camel whey protein alleviated heat stress-induced hepatocyte damage by activated Nrf2/HO-1 signaling pathway and inhibited NF- κ B/NLRP3 axis *Cell Stress Chaperones* 2021 Mar;26(2):387-401.
 41. Cai-Chou Zhao, Hao Guo, Ying Wang, Jiu-Hong Li Comprehensive upstream and downstream regulatory analyses identify miR-675-3p as a potential prognostic biomarker in melanoma *Hum Cell* 2021 Mar;34(2):654-666.
 42. Shan-Shan Li, Ling-Ling Xie, Zhuang-Zhuang Li, Yong-Jian Fan, Man-Man Qi, Yan-Guo Xi Androgen is responsible for enhanced susceptibility of melatonin against traumatic brain injury in females *Neurosci Lett* 2021 May 1;752:135842.
 43. Qudong Lu, Yang Yang, Hengshuai Zhang, Cheng Chen, Jiang Zhao, Zhenxing Yang, Yi Fan, Longkun Li, Huan Feng, Jingzhen Zhu, Shanhong Yi Activation of GPR18 by Resolvin D2 Relieves Pain and Improves Bladder Function in Cyclophosphamide-Induced Cystitis Through Inhibiting TRPV1 *Drug Des Devel Ther* 2021 Nov 15;15:4687-4699.
 44. Jiaqi Chen, Qiaoya Pan, Yang Bai, Xuepeng Chen, Yi Zhou Hydroxychloroquine Induces Apoptosis in Cholangiocarcinoma via Reactive Oxygen Species Accumulation Induced by Autophagy Inhibition *Front Mol Biosci* 2021 Sep 10;8:720370.
 45. Guohai Xie, Xinyi Zheng, Zhong Zheng, Ruoyu Wu, Zhixian Yao, Wenjie Huang, Feng Sun, Xingyu Mu, Ke Wu, Junhua Zheng The ceRNA PVT1 inhibits proliferation of ccRCC cells by sponging miR-328-3p to elevate FAM193B expression *Aging (Albany NY)* 2021 Sep 13;13(17):21712-21728.
 46. Song Yu, Bingxuan Jia, Na Liu, Dianzhen Yu, Shuo Zhang, Aibo Wu Fumonisin B1 triggers carcinogenesis via HDAC/PI3K/Akt signalling pathway in human esophageal epithelial cells *Sci Total Environ* 2021 Sep 15;787:147405.
 47. Hang-Xing Yu, Wei Lin, Kang Yang, Li-Juan Wei, Jun-Li Chen, Xin-Yue Liu, Ke Zhong, Xin Chen, Ming Pei, Hong-Tao Yang Transcriptome-Based Network Analysis Reveals Hirudin Potentiates Anti-Renal Fibrosis Efficacy in UOU Rats *Front Pharmacol* 2021 Sep 21;12:741801.
 48. Lishan Zhang, Shanhui Ge, Wanmei He, Qingui Chen, Caixia Xu, Mian Zeng Ghrelin protects against lipopolysaccharide-induced acute respiratory distress syndrome through the PI3K/AKT pathway *J Biol Chem* 2021 Sep;297(3):101111.
 49. Si-Min Zhang, Chuan-Yuan Wei, Qiang Wang, Lu Wang, Lu Lu, Fa-Zhi Qi M2-polarized macrophages mediate wound healing by regulating connective tissue growth factor via AKT, ERK1/2, and STAT3 signaling pathways *Mol Biol Rep* 2021 Sep;48(9):6443-6456.
 50. Yinpei Luo, Hong Yang, Xiaojing Yan, Yaran Wu, Guoliang Wei, Xiaoying Wu, Xuelong Tian, Ying Xiong, Guangyan Wu, Huizhong Wen Transcranial Direct Current Stimulation Alleviates Neurovascular Unit Dysfunction in Mice With Preclinical Alzheimer's Disease *Front Aging Neurosci* 2022 Apr 14:857415.
 51. Xiaohong Sun, Keda Zhu, Chengcheng Feng, Jie Zhu, Shuangshuang Chen, Wenkai Tang, Zhifang Wang, Long Xiao, Hong Li, Dechun Geng, Zhirong Wang Paeoniflorin Ameliorates Hyperprolactinemia-Induced Inhibition of Osteoblastogenesis by Suppressing the NF- κ B Signaling Pathway *Int J Endocrinol* 2022 Apr 15:2022:4572033.
 52. Wenwen Su, Leilei Wu, Qichao Liang, Xiaoyue Lin, Xiaoyi Xu, Shikai Yu, Yitong Lin, Jiadong Zhou, Yang Fu, Xiaoyan Gao, Bo Zhang, Li Li, Dan Li, Yongkui Yin, Gaochen Song Extraction Optimization, Structural Characterization, and Anti-Hepatoma Activity of Acidic Polysaccharides From *Scutellaria barbata* D. Don *Front Pharmacol* 2022 Apr 4:13:827782.
 53. Wenqiang Zhu, Chen Ding, Piaopiao Huang, Juanli Ran, Pingan Lian, Yaxin Tang, Wen Dai, Xiansheng Huang Metformin Ameliorates Hepatic Steatosis induced by olanzapine through inhibiting LXR α /PCSK9 pathway *Sci Rep* 2022 Apr 4;12(1):5639.
 54. Guanghai Guo, Jin Dong Diosmetin attenuates oxidative stress-induced damage to lens epithelial cells via the mitogen-activated protein kinase (MAPK) pathway *Bioengineered* 2022 Apr;13(4):11072-11081.
 55. Ling-Ling Xie, Shan-Shan Li, Yong-Jian Fan, Man-Man Qi, Zhuang-Zhuang Li Melatonin alleviates traumatic brain injury-induced anxiety-like behaviors in rats: Roles of the protein kinase A/cAMP-response element binding signaling pathway *Exp Ther Med* 2022 Apr;23(4):248.
 56. Xinyi Zhan, Weijie Peng, Zhuqiang Wang, Xin Liu, Weibo Dai, Quanxi Mei, Xianjing Hu Polysaccharides from Garlic Protect against Liver Injury in DSS-Induced Inflammatory Bowel Disease of Mice via Suppressing Pyroptosis and Oxidative Damage *Oxid Med Cell Longev* 2022 Aug 16:2022:2042163.
 57. Cong Pang, Sen Gao, Xun-Zhi Liu, Xiao-Jian Li, Zheng Peng, Hua-Sheng Zhang, Yan Zhou, Xiang-Xin Chen, Tao Tao, Yue Lu, Wei Li, Chun-Hua Hang Astrocytic CD24 Protects Neuron from Recombinant High-Mobility Group Box 1 Protein(rHMGBl)-Elicited Neuronal Injury *Brain Sci* 2022 Aug 23;12(9):1119.
 58. Tianyang Lan, Kang Zhang, Feifei Lin, Qifu He, Shenghui Wu, Zhiming Xu, Yong Zhang, Fusheng Quan Effects of MICU1-Mediated Mitochondrial Calcium Uptake on Energy Metabolism and Quality of Vitri-fied-Thawed Mouse Metaphase II Oocytes *Int J Mol Sci* 2022 Aug 3;23(15):8629.
 59. Fan Wang, Yishan Lu, Junming Cao Dynamics impacts of oxytetracycline on growth performance, intestinal health and antibiotic residue of grouper in exposure and withdrawal treatment *Ecotoxicol Environ Saf* 2022 Dec 1:247:114203.
 60. Tao Ma, Yan Chen, Zhi-Gang Yi, Jia Liu, Yan-Hong Li, Jun Bai, Wen-Ting Tie, Mei Huang, Xiao-Feng Zhu, Ji Wang, Juan Du, Xiu-Qin Zuo, Qin Li, Fan-Li Lin, Liu Tang, Jing Guo, Hong-Wen Xiao, Qian Lei, Xiao-Li Ma, Li-Juan Li, Lian-Sheng Zhang NORAD promotes multiple myeloma cell progression via BMP6/P-ERK1/2 axis *Cell Signal* 2022 Dec:100:110474.
 61. Zhuotong Zeng, Yaoyao Wang, Yangfan Xiao, Jie Zheng, Ruizhen Liu, Xinglan He, Jiangfan Yu, Bingsi Tang, Xiangning Qiu, Rui Tang, Yaqian Shi, Rong Xiao Overexpression of OASL upregulates TET1 to induce aberrant activation of CD4+ T cells in systemic sclerosis via IRF1 signaling *Arthritis Res Ther* 2022 Feb 19;24(1):50.
 62. Kaixi Wang, Jiafan Wu, Shuna Chen, Hangye Zhao, Puming He, Youying Tu, Bo Li Transcriptome analysis provides insight into the anti-diabetic mechanism of theaflavins in high-fat diet and streptozotocin-induced mice *Food Funct* 2022 Feb 21;13(4):2033-2043.
 63. Zenglin Cui, Yuwei Li, Gaorui Liu, Yanmeng Jiang miR-103a-3p Silencing Ameliorates Calcium Oxalate Deposition in Rat Kidney by Activating the UMOD/TRPV5 Axis *Dis Markers* 2022 Feb 23:2022:2602717.
 64. Huawei Li, Xiaoling Chen, Zhiqing Huang, Daiwen Chen, Bing Yu, Yuheng Luo, Jun He, Ping Zheng, Jie Yu, Hong Chen Ellagic acid enhances muscle endurance by affecting the muscle fiber type, mitochondrial biogenesis and function *Food Funct* 2022 Feb 7;13(3):1506-1518.
 65. Xiaoqin Tan, Chunpu Li, Ruirui Yang, Sen Zhao, Fei Li, Xutong Li, Lifan Chen, Xiaozhe Wan, Xiaohong Liu, Tianbiao Yang, Xiaochu Tong, Tingyang Xu, Rongrong Cui, Hualiang Jiang, Sulin Zhang, Hong Liu, Mingyue Zheng Discovery of Pyrazolo[3,4-d]pyridazinone Derivatives as Selective DDR1 Inhibitors via Deep Learning Based Design, Synthesis, and Biological Evaluation *J Med Chem* 2022 Jan 13;65(1):103-119.
 66. Yihan Wang, Xiang Yuan, Malik Ahsan Ali, Ziyue Qin, Yan Zhang, Changjun Zeng piR-121380 Is Involved in Cryo-Capacitation and Regulates Post-Thawed Boar Sperm Quality Through Phosphorylation of ERK2 via Targeting PTPN7 *Front Cell Dev Biol* 2022 Jan 26;9:792994.
 67. Penghui Xu, Xing Zhang, Jiacheng Cao, Jing Yang, Zetian Chen, Weizhi Wang, Sen Wang, Lu Zhang, Li Xie, Lang Fang, Yiwen Xia, Zhe Xuan, Jialun Lv, Hao Xu,

- Zekuan Xu The novel role of circular RNA ST3GAL6 on blocking gastric cancer malignant behaviours through autophagy regulated by the FOXp2/MET/mTOR axis *Clin Transl Med* 2022 Jan;12(1):e707.
68. Min Zhang, Shaoyang Sun, Lei Wang, Xu Wang, Tianhui Chen, Zexu Chen, Yongxiang Jiang Zonular defects in lox1-deficient zebrafish *Clin Exp Ophthalmol* 2022 Jan;50(1):62-73.
69. Yuanyuan Wu, Wensi Zhu, Ainiwaer Rouzi, Lin Tong, Linxiao Han, Juan Song, Jianwen Ding, Yu Yan, Miao Li, Ting Pan, Jie Liu, Qin Wang, Yuanlin Song, Jie Shen, Jian Zhou The traditional Chinese patented medicine Qingke Pingchuan granules alleviate acute lung injury by regenerating club cells *Pulm Circ* 2022 Jul 1;12(3):e12138.
70. Jun Yan Xu, Yu Yan Xiong, Rui Jie Tang, Wen Yang Jiang, Yu Ning, Zhao Ting Gong, Pei Sen Huang, Gui Hao Chen, Jun Xu, Chun Xiao Wu, Meng Jin Hu, Jing Xu, Yi Xu, Cun Rong Huang, Chen Jin, Xiao Tong Lu, Hai Yan Qian, Xiang Dong Li, Yue Jin Yang Interleukin-5-induced eosinophil population improves cardiac function after myocardial infarction *Cardiovasc Res* 2022 Jul 20;118(9):2165-2178.
71. Sheng He, Lu Li, Haifan Chen, Xiaoli Hu, Wendi Wang, Hui Zhang, Ruiping Wei, Xiaoxiao Zhang, Yaosheng Chen, Xiaohong Liu PRRSV Infection Induces Gasdermin D-Driven Pyroptosis of Porcine Alveolar Macrophages through NLRP3 Inflammasome Activation *J Virol* 2022 Jul 27;96(14):e0212721.
72. Huijie Zhang, Zhixin He, Ping Deng, Muxue Lu, Chao Zhou, Lingling Yang, Zhengping Yu PIN1-mediated ROS production is involved in antagonism of N-acetyl-L-cysteine against arsenic-induced hepatotoxicity *Toxicol Res (Camb)* 2022 Jul 8;11(4):628-643.
73. Yanqi Zhong, Yang Zhang, Weifang Liu, Yin Zhao, Li Zou, Xiaoxia Liu TLR4 Modulates Senescence and Paracrine Action in Placental Mesenchymal Stem Cells via Inhibiting Hedgehog Signaling Pathway in Preeclampsia *Oxid Med Cell Longev* 2022 Jun 14;2022:7202837.
74. Ruru Gao, Qiong Luo, Yang Li, Liming Song, Junnan Stephen Cai, Ying Xiong, Fei Yan, Jianhua Liu Biosynthetic Nanobubble-Mediated CRISPR/Cas9 Gene Editing of Cdh2 Inhibits Breast Cancer Metastasis *Pharmaceutics* 2022 Jun 30;14(7):1382.
75. Lishan Huang, Zhou Chen, Ruiyu Chen, Lu Lin, Lingjia Ren, Meilian Zhang, Libin Liu Increased fatty acid metabolism attenuates cardiac resistance to β -adrenoceptor activation via mitochondrial reactive oxygen species: A potential mechanism of hypoglycemia-induced myocardial injury in *Redox Biol* 2022 Jun;52:102320.
76. Shuai Zhang, Mengyu Xing, Gaojun Chen, Lei Tong, Haili Zhang, Dongshu Du Up-regulation of miR-335 and miR-674-3p in the rostral ventrolateral medulla contributes to stress-induced hypertension *J Neurochem* 2022 Jun;161(5):387-404.
77. Xing Wei, Andrew Chia Hao Chang, Haishuang Chang, Shan Xu, Yilin Xue, Yuanxin Zhang, Ming Lei, Alex Chia Yu Chang, Qingyong Zhang Hypoglycemia-Exacerbated Mitochondrial Connexin 43 Accumulation Aggravates Cardiac Dysfunction in Diabetic Cardiomyopathy *Front Cardiovasc Med* 2022 Mar 16;9:800185.
78. Yang Bai, Jiaqi Chen, Weijian Hu, Lei Wang, Yulian Wu, Shi'an Yu Silibinin Therapy Improves Cholangiocarcinoma Outcomes by Regulating ERK/Mitochondrial Pathway *Front Pharmacol* 2022 Mar 23;13:847905.
79. Lu Jiang, Chunhua Chi, Fang Yuan, Meiqi Lu, Dongqing Hu, Lin Wang, Xiaoming Liu Anti-inflammatory effects of anemonin on acute ulcerative colitis via targeted regulation of protein kinase C- θ *Chin Med* 2022 Mar 28;17(1):39.
80. Siyuan Chen, Yajie Xue, Yutian Shen, Hao Ju, Xiaodong Zhang, Jinsong Liu, Yongxia Wang Effects of different selenium sources on duodenum and jejunum tight junction network and growth performance of broilers in a model of fluorine-induced chronic oxidative stress *Poult Sci* 2022 Mar;101(3):101664.
81. Zheng-Rong Gao, Qiong Liu, Jie Zhao, Ya-Qiong Zhao, Li Tan, Shao-Hui Zhang, Ying-Hui Zhou, Yun Chen, Yue Guo, Yun-Zhi Feng A comprehensive analysis of the circRNA-miRNA-mRNA network in osteocyte-like cell associated with Mycobacterium leprae infection *PLoS Negl Trop Dis* 2022 May 2;16(5):e0010379.
82. Zeyu Liu, Moli Huang, Yue Hong, Shaoyang Wang, Yongle Xu, Cheng Zhong, Jingyuan Zhang, Zhengping Zhuang, Shan Shan, Tao Ren Isovalerylsipramycin I suppresses non-small cell lung carcinoma growth through ROS-mediated inhibition of PI3K/AKT signaling pathway *Int J Biol Sci* 2022 May 21;18(9):3714-3730.
83. Fang Hong, Guiyan He, Manqi Zhang, Boyang Yu, Chengzhi Chai The Establishment of a Mouse Model of Recurrent Primary Dysmenorrhea *Int J Mol Sci* 2022 May 30;23(11):6128.
84. Cong Zhang, Lijin Zeng, Guoyi Cai, Yuanling Zhu, Yan Xiong, Hong Zhan, Zhen Yang miR-340-5p Alleviates Oxidative Stress Injury by Targeting MyD88 in Sepsis-Induced Cardiomyopathy *Oxid Med Cell Longev* 2022 May 4;2022:2939279.
85. Tong Yi, Wenxin Ding, Yuanzhen Hao, Lifeng Cen, Jiyang Li, Xunlong Shi, Ting Wang, Daofeng Chen, Haiyan Zhu Neutrophil extracellular traps mediate severe lung injury induced by influenza A virus H1N1 in mice coinfecting with Staphylococcus aureus *Microb Pathog* 2022 May;166:105558.
86. Changgui Wu, Shaohua Chen, Yang Liu, Bo Kong, Wei Yan, Tao Jiang, Hao Tian, Zhaoyi Liu, Qi Shi, Yongjun Wang, Qianqian Liang, Xiaobing Xi, Hao Xu Cynarin suppresses gouty arthritis induced by monosodium urate crystals *Bioengineered* 2022 May;13(5):11782-11793.
87. Jia-Hui Sun, Ming Huang, Zhou Fang, Tian-Xiao Li, Ting-Ting Wu, Yi Chen, Da-Ping Quan, Ying-Ying Xu, Yu-Ming Wang, Yi Yang, Jian-Long Zou Nerve bundle formation during the promotion of peripheral nerve regeneration: collagen VI-neural cell adhesion molecule 1 interaction *Neural Regen Res* 2022 May;17(5):1023-1033.
88. Shan Shan, Yeping Yang, Jilan Jiang, Bingxin Yang, Yisai Yang, Feng Sun, Junyu Zhang, Yu Lin, Hong Xu Extracellular vesicle-derived long non-coding RNA as circulating biomarkers for endometriosis *Reprod Biomed Online* 2022 May;44(5):923-933.
89. Yangjun Yang, Xi Li, Zonghan Liu, Xinyu Ruan, Huihui Wang, Qiang Zhang, Lu Cao, Luchen Song, Yinghong Chen, Yi Sun Moderate Treadmill Exercise Alleviates NAFLD by Regulating the Biogenesis and Autophagy of Lipid Droplet Nutrients *2022 Nov 20;14(22):4910.*
90. Zhilan Peng, Jialong Gao, Weimin Su, Wenhong Cao, Guoping Zhu, Xiaoming Qin, Chaohua Zhang, Yi Qi Purification and Identification of Peptides from Oyster (*Crassostrea hongkongensis*) Protein Enzymatic Hydrolysates and Their Anti-Skin Photoaging Effects on UVB-Irradiated HaCaT Cells *Mar Drugs* 2022 Nov 28;20(12):749.
91. Yu Liu, Yan Zhang, Huanhuan Zhu, Wenzhi Shen, Zheng Chen, Jian Bai, Tian Shuang, Qi Chen Aucubin administration suppresses STING signaling and mitigated high-fat diet-induced atherosclerosis and steatohepatitis in LDL receptor deficient mice *Food Chem Toxicol* 2022 Nov;169:113422.
92. Leyu Li, Guoxin Huang, Tingbo Chen, Hui Lin, Ruiyan Xu, Jinyan Cheng, Ying Hu, Weibo Dai, Gengting Dong Fufang Fanshiliu Decoction Revealed the Antidiabetic Effect through Modulating Inflammatory Response and Gut Microbiota Composition *Evid Based Complement Alternat Med* 2022 Oct 10;2022:3255401.
93. Xiaoyong Chen, Ziwei Li, Shuaiwei Wang, Guangzhi Tong, Keyuan Chen, Yan Zhao Proteomic analysis reveals zinc-finger CCHC-type containing protein 3 as a factor inhibiting virus infection by promoting innate signaling *Virus Res* 2022 Oct 2;319:198876.
94. Mu-Zi Li, Xiao-Yang Wen, Xiao-Qiang Liu, Yu-Qing Wang, Lei Yan LPS-Induced Activation of the cGAS-STING Pathway is Regulated by Mitochondrial Dysfunction and Mitochondrial DNA Leakage in Endometritis *J Inflamm Res* 2022

- Oct 5:15:5707-5720.
95. Yueshuai Liu, Hongxiang Ding, Yuze Yang, Yan Liu, Xin Cao, Tao Feng Progesterone Induces Apoptosis and Steroidogenesis in Porcine Placental Trophoblasts *Animals (Basel)* 2022 Oct 8;12(19):2704.
 96. Hai-Qun Dong, Shi-Jing Liang, Yu-Ling Xu, Yi Dai, Na Sun, Dong-Hong Deng, Peng Cheng Liproxstatin-1 induces cell cycle arrest, apoptosis, and caspase-3/GSDME-dependent secondary pyroptosis in K562 cells *Int J Oncol* 2022 Oct;61(4):119.
 97. Junwen Wang, Xintong Zhu, Limeng Dai, Ziyi Wang, Xingying Guan, Xiaoyin Tan, Jia Li, Mao Zhang, Yun Bai, Hong Guo Supt16 haploinsufficiency causes neurodevelopment disorder by disrupting MAPK pathway in neural stem cells *Hum Mol Genet* 2023 Feb 19;32(5):860-872.
 98. Wen Yang, Chunwang Jia, Long Liu, Yu Fu, Yawei Wu, Zhicheng Liu, Ruixuan Yu, Xiaojie Ma, Ao Gong, Fangming Liu, Yanni Xia, Yong Hou, Yuhua Li, Lei Zhang Hypoxia-Inducible Factor-1 α Protects Against Intervertebral Disc Degeneration Through Antagonizing Mitochondrial Oxidative Stress Inflammation *Inflammation* 2023 Feb;46(1):270-284.
 99. Fanghua Chen, Yin Li, Ling Aye, Yingcheng Wu, Liangqing Dong, Zijian Yang, Qiang Gao, Shu Zhang FUT8 is regulated by miR-122-5p and promotes malignancies in intrahepatic cholangiocarcinoma via PI3K/AKT signaling *Cell Oncol (Dordr)* 2023 Feb;46(1):79-91.
 100. Fan Gu, Hailin Wu, Zhuo Huang, Fei Wang, Ruihuan Yang, Zhuan Bian, Miao He The effects of dimethyl fumarate on cytoplasmic LPS-induced noncanonical pyroptosis in periodontal ligament fibroblasts and dental pulp cells *Int Endod J* 2023 Jul;56(7):869-880.
 101. Weijie Zhong, Juan Cheng, Xiaosheng Yang, Wenwu Liu, Yi Li Heliox Preconditioning Exerts Neuroprotective Effects on Neonatal Ischemia/Hypoxia Injury by Inhibiting Necroptosis Induced by Ca²⁺ Elevation *Transl Stroke Res* 2023 Jun;14(3):409-424.
 102. Ye Kuang, Yun Cheng, Jia Wang, Hongyan Li, Xianghong Cao, Yang Wang KIAA1429 mediates epithelial mesenchymal transition in sorafenib-resistant hepatocellular carcinoma through m6A methylation modification *Cancer Med* 2023 Mar;12(6):7222-7233.
 103. Qiong Fang, Zhiying Li, Ye Xue, Xin Zong, Wenshuang Ma, Guangmin Xi, Xiao Feng Zhang, Zuwei Li Embelin Enhances the Sensitivity of Renal Cancer Cells to Axitinib by Inhibiting HIF Signaling Pathway *Anticancer Agents Med Chem* 2023;23(7):807-816.
 104. Ding M1, 2, 3, Ning J2, Feng N2, Li Z2, Liu Z2, Wang Y2, Wang Y2, Li X3, Huo C3, Jia X3, Xu R3, Fu F2, Wang X3, Pei J2. Dynamin - related protein 1 - mediated mitochondrial fission contributes to post - traumatic cardiac dysfunction in rats and the protective effect of melatonin. *J Pineal Res* 2018 Jan;64(1). doi: 10.1111/jpi.12447.
 105. Ding M1, 2, 3, Ning J2, Feng N2, Li Z2, Liu Z2, Wang Y2, Wang Y2, Li X3, Huo C3, Jia X3, Xu R3, Fu F2, Wang X3, Pei J2. Dynamin - related protein 1 - mediated mitochondrial fission contributes to post - traumatic cardiac dysfunction in rats and the protective effect of melatonin. *J Pineal Res* 2018 Jan;64(1). doi: 10.1111/jpi.12447.
 106. Ding M1, 2, 3, Ning J2, Feng N2, Li Z2, Liu Z2, Wang Y2, Wang Y2, Li X3, Huo C3, Jia X3, Xu R3, Fu F2, Wang X3, Pei J2. Dynamin - related protein 1 - mediated mitochondrial fission contributes to post - traumatic cardiac dysfunction in rats and the protective effect of melatonin. *J Pineal Res* 2018 Jan;64(1). doi: 10.1111/jpi.12447.
 107. Zhao Z1, Zhao Y1, Zhuang XY2, Lo WC3, Baker MAB4, Lo CJ5, Bai F Frequent pauses in Escherichia coli flagella elongation revealed by single cell real-time fluorescence imaging. *Nature Communications* 2018 May 14;9(1):1885.
 108. Zhao Z1, Zhao Y1, Zhuang XY2, Lo WC3, Baker MAB4, Lo CJ5, Bai F Frequent pauses in Escherichia coli flagella elongation revealed by single cell real-time fluorescence imaging. *Nature Communications* 2018 May 14;9(1):1885.
 109. Zhao Z1, Zhao Y1, Zhuang XY2, Lo WC3, Baker MAB4, Lo CJ5, Bai F Frequent pauses in Escherichia coli flagella elongation revealed by single cell real-time fluorescence imaging. *Nature Communications* 2018 May 14;9(1):1885.
 110. Pei-PeiZhao, Hao-RanHu, Jia-YuLiu, Qin-FeiKe, Xiao-YuanPeng, HaoDing, Ya-PingGuo. Gadolinium phosphate/chitosan scaffolds promote new bone regeneration via Smad/Runx2 pathway. *CHEM ENG J* 2019 March 1;359:1120-1129.
 111. Pei-PeiZhao, Hao-RanHu, Jia-YuLiu, Qin-FeiKe, Xiao-YuanPeng, HaoDing, Ya-PingGuo. Gadolinium phosphate/chitosan scaffolds promote new bone regeneration via Smad/Runx2 pathway. *CHEM ENG J* 2019 March 1;359:1120-1129.
 112. Pei-PeiZhao, Hao-RanHu, Jia-YuLiu, Qin-FeiKe, Xiao-YuanPeng, HaoDing, Ya-PingGuo. Gadolinium phosphate/chitosan scaffolds promote new bone regeneration via Smad/Runx2 pathway. *CHEM ENG J* 2019 March 1;359:1120-1129.
 113. Pei-PeiZhao, Hao-RanHu, Jia-YuLiu, Qin-FeiKe, Xiao-YuanPeng, HaoDing, Ya-PingGuo. Gadolinium phosphate/chitosan scaffolds promote new bone regeneration via Smad/Runx2 pathway. *CHEM ENG J* 2019 March 1;359:1120-1129.
 114. Long J1, 2, Xiao Y3, Liu L4, Cao Y5, 6. The adverse vascular effects of multi-walled carbon nanotubes (MWCNTs) to human vein endothelial cells (HUVECs) in vitro: role of length of MWCNTs. *J NANOBIOTECHNOL*2017 Nov 10;15(1):80.
 115. Long J1, 2, Xiao Y3, Liu L4, Cao Y5, 6. The adverse vascular effects of multi-walled carbon nanotubes (MWCNTs) to human vein endothelial cells (HUVECs) in vitro: role of length of MWCNTs. *J NANOBIOTECHNOL*2017 Nov 10;15(1):80.
 116. Long J, Xiao Y, Liu L, Cao Y. The adverse vascular effects of multi-walled carbon nanotubes (MWCNTs) to human vein endothelial cells (HUVECs) in vitro: role of length of MWCNTs. *J NANOBIOTECHNOL*2017 Nov 10;15(1):80.
 117. Long J, Xiao Y, Liu L, Cao Y. The adverse vascular effects of multi-walled carbon nanotubes (MWCNTs) to human vein endothelial cells (HUVECs) in vitro: role of length of MWCNTs. *J NANOBIOTECHNOL*2017 Nov 10;15(1):80.
 118. Long J, Xiao Y, Liu L, Cao Y. The adverse vascular effects of multi-walled carbon nanotubes (MWCNTs) to human vein endothelial cells (HUVECs) in vitro: role of length of MWCNTs. *J NANOBIOTECHNOL*2017 Nov 10;15(1):80.
 119. Long J, Xiao Y, Liu L, Cao Y. The adverse vascular effects of multi-walled carbon nanotubes (MWCNTs) to human vein endothelial cells (HUVECs) in vitro: role of length of MWCNTs. *J NANOBIOTECHNOL*2017 Nov 10;15(1):80.
 120. Long J, Xiao Y, Liu L, Cao Y. The adverse vascular effects of multi-walled carbon nanotubes (MWCNTs) to human vein endothelial cells (HUVECs) in vitro: role of length of MWCNTs. *J NANOBIOTECHNOL*2017 Nov 10;15(1):80.
 121. Hu H1, Zhao P2, Liu J2, Ke Q2, Zhang C1, Guo Y3, Ding H4. Lanthanum phosphate/chitosan scaffolds enhance cytocompatibility and osteogenic efficiency via the Wnt/ β -catenin pathway. *J NANOBIOTECHNOL*2018 Nov 29;16(1):98.
 122. Hu H1, Zhao P2, Liu J2, Ke Q2, Zhang C1, Guo Y3, Ding H4. Lanthanum phosphate/chitosan scaffolds enhance cytocompatibility and osteogenic efficiency via the Wnt/ β -catenin pathway. *J NANOBIOTECHNOL*2018 Nov 29;16(1):98.
 123. Hu H1, Zhao P2, Liu J2, Ke Q2, Zhang C1, Guo Y3, Ding H4. Lanthanum phosphate/chitosan scaffolds enhance cytocompatibility and osteogenic efficiency via the Wnt/ β -catenin pathway. *J NANOBIOTECHNOL*2018 Nov 29;16(1):98.
 124. Chen GL1, 2, Ye T1, Chen HL3, Zhao ZY4, Tang WQ4, Wang LS3, Xia J1, 4. Xanthine dehydrogenase downregulation promotes TGF β signaling and cancer stem cell-related gene expression in hepatocellular carcinoma. *Oncogenesis* 2017 Sep 25;6(9):e382.
 125. Chen GL1, 2, Ye T1, Chen HL3, Zhao ZY4, Tang WQ4, Wang LS3, Xia J1, 4.

- Xanthine dehydrogenase downregulation promotes TGF β signaling and cancer stem cell-related gene expression in hepatocellular carcinoma. *Oncogenesis* 2017 Sep 25;6(9):e382.
126. Chen GL, Ye T1, Chen HL3, Zhao ZY4, Tang WQ4, Wang LS3, Xia JL1, 4. Xanthine dehydrogenase downregulation promotes TGF β signaling and cancer stem cell-related gene expression in hepatocellular carcinoma. *Oncogenesis* 2017 Sep 25;6(9):e382.
127. Fu C, Yin D, Nie H, Sun D. Notoginsenoside R1 Protects HUVEC Against Oxidized Low Density Lipoprotein (Ox-LDL)-Induced Atherogenic Response via Down-Regulating miR-132. *CELL PHYSIOL BIOCHEM* 2018;51(4):1739-1750.
128. Fu C, Yin D, Nie H, Sun D. Notoginsenoside R1 Protects HUVEC Against Oxidized Low Density Lipoprotein (Ox-LDL)-Induced Atherogenic Response via Down-Regulating miR-132. *CELL PHYSIOL BIOCHEM* 2018;51(4):1739-1750.
129. Fu C, Yin D, Nie H, Sun D. Notoginsenoside R1 Protects HUVEC Against Oxidized Low Density Lipoprotein (Ox-LDL)-Induced Atherogenic Response via Down-Regulating miR-132. *CELL PHYSIOL BIOCHEM* 2018;51(4):1739-1750.
130. Tengfei Liu, Xiaoyan Liu, Hui Xiong, Cheng Xu, Jianxu Yao, Xiumei Zhu, Jianping Zhou and Jing Yao. Mechanisms of TPGS and its derivatives inhibiting P-glycoprotein efflux pump and application for reversing multidrug resistance in hepatocellular carcinoma. *POLYM CHEM-UK* 2018;9:1827-1839.
131. Tengfei Liu, Xiaoyan Liu, Hui Xiong, Cheng Xu, Jianxu Yao, Xiumei Zhu, Jianping Zhou and Jing Yao. Mechanisms of TPGS and its derivatives inhibiting P-glycoprotein efflux pump and application for reversing multidrug resistance in hepatocellular carcinoma. *POLYM CHEM-UK* 2018;9:1827-1839.
132. Tengfei Liu, Xiaoyan Liu, Hui Xiong, Cheng Xu, Jianxu Yao, Xiumei Zhu, Jianping Zhou and Jing Yao. Mechanisms of TPGS and its derivatives inhibiting P-glycoprotein efflux pump and application for reversing multidrug resistance in hepatocellular carcinoma. *POLYM CHEM-UK* 2018;9:1827-1839.
133. Tengfei Liu, Xiaoyan Liu, Hui Xiong, Cheng Xu, Jianxu Yao, Xiumei Zhu, Jianping Zhou and Jing Yao. Mechanisms of TPGS and its derivatives inhibiting P-glycoprotein efflux pump and application for reversing multidrug resistance in hepatocellular carcinoma. *POLYM CHEM-UK* 2018;9:1827-1839.
134. Yang HW1, Liu XY2, Shen ZF1, Yao W1, Gong XB2, Huang HX3, Ding GH1. An investigation of the distribution and location of mast cells affected by the stiffness of substrates as a mechanical niche. *Int J Biol Sci* 2018 Jun 22;14(9):1142-1152.
135. Yang HW1, Liu XY2, Shen ZF1, Yao W1, Gong XB2, Huang HX3, Ding GH1. An investigation of the distribution and location of mast cells affected by the stiffness of substrates as a mechanical niche. *Int J Biol Sci* 2018 Jun 22;14(9):1142-1152.
136. Yang HW1, Liu XY2, Shen ZF1, Yao W1, Gong XB2, Huang HX3, Ding GH1. An investigation of the distribution and location of mast cells affected by the stiffness of substrates as a mechanical niche. *Int J Biol Sci* 2018 Jun 22;14(9):1142-1152.
137. Jin X1, Yu Y1, Zou Q2, Wang M1, Cui Y1, Xie J1, Wang Z1. MicroRNA - 105 promotes epithelial - mesenchymal transition of nonsmall lung cancer cells through upregulating Mcl - 1. *J Cell Biochem* 2018 Oct 14. doi: 10.1002/jcb.27873.
138. Jin X1, Yu Y1, Zou Q2, Wang M1, Cui Y1, Xie J1, Wang Z1. MicroRNA - 105 promotes epithelial - mesenchymal transition of nonsmall lung cancer cells through upregulating Mcl - 1. *J Cell Biochem* 2018 Oct 14. doi: 10.1002/jcb.27873.
139. Jin X1, Yu Y1, Zou Q2, Wang M1, Cui Y1, Xie J1, Wang Z1. MicroRNA - 105 promotes epithelial - mesenchymal transition of nonsmall lung cancer cells through upregulating Mcl - 1. *J Cell Biochem* 2018 Oct 14. doi: 10.1002/jcb.27873.
140. Ding Y1, Zheng Y1, Liu T1, Chen T1, Wang C1, Sun Q1, Hua M2, Hua T3. Changes in GABAergic markers accompany degradation of neuronal function in the primary visual cortex of senescent rats. *SCI REP-UK* 2017 Nov 2;7(1):14897.
141. Ding Y1, Zheng Y1, Liu T1, Chen T1, Wang C1, Sun Q1, Hua M2, Hua T3. Changes in GABAergic markers accompany degradation of neuronal function in the primary visual cortex of senescent rats. *SCI REP-UK* 2017 Nov 2;7(1):14897.
142. Ding Y, Zheng Y, Liu T, Chen T, Wang C, Sun Q, Hua M, Hua T. Changes in GABAergic markers accompany degradation of neuronal function in the primary visual cortex of senescent rats. *SCI REP-UK* 2017 Nov 2;7(1):14897.
143. Guo D1, Hu X2, Zhang H3, Lu C4, Cui G5, Luo X6. Orientin and neuropathic pain in rats with spinal nerve ligation. *Int Immunopharmacol* 2018 May;58:72-79.
144. Guo D1, Hu X2, Zhang H3, Lu C4, Cui G5, Luo X6. Orientin and neuropathic pain in rats with spinal nerve ligation. *Int Immunopharmacol* 2018 May;58:72-79.
145. Guo D1, Hu X2, Zhang H3, Lu C4, Cui G5, Luo X6. Orientin and neuropathic pain in rats with spinal nerve ligation. *Int Immunopharmacol* 2018 May;58:72-79.
146. Guo D1, Hu X2, Zhang H3, Lu C4, Cui G5, Luo X6. Orientin and neuropathic pain in rats with spinal nerve ligation. *Int Immunopharmacol* 2018 May;58:72-79.
147. Li XH1, Wang HP1, Tan J1, Wu YD1, Yang M1, Mao CZ2, Gao SF2, Li H2, Chen H3, Cai WB4. Loss of pigment epithelium-derived factor leads to ovarian oxidative damage accompanied by diminished ovarian reserve in mice. *Life Sci* 2019 Jan 1;216:129-139.
148. Li XH1, Wang HP1, Tan J1, Wu YD1, Yang M1, Mao CZ2, Gao SF2, Li H2, Chen H3, Cai WB4. Loss of pigment epithelium-derived factor leads to ovarian oxidative damage accompanied by diminished ovarian reserve in mice. *Life Sci* 2019 Jan 1;216:129-139.
149. Li XH1, Wang HP1, Tan J1, Wu YD1, Yang M1, Mao CZ2, Gao SF2, Li H2, Chen H3, Cai WB4. Loss of pigment epithelium-derived factor leads to ovarian oxidative damage accompanied by diminished ovarian reserve in mice. *Life Sci* 2019 Jan 1;216:129-139.
150. Li XH1, Wang HP1, Tan J1, Wu YD1, Yang M1, Mao CZ2, Gao SF2, Li H2, Chen H3, Cai WB4. Loss of pigment epithelium-derived factor leads to ovarian oxidative damage accompanied by diminished ovarian reserve in mice. *Life Sci* 2019 Jan 1;216:129-139.
151. Liu L1, Zhang Y1, Chang X1, Li R2, Wu C1, Tang L3, Zhou Z4. Fluorochloridone perturbs blood-testis barrier/Sertoli cell barrier function through Arp3-mediated F-actin disruption. *Toxicol Lett* 2018 Oct 1;295:277-287.
152. Liu L1, Zhang Y1, Chang X1, Li R2, Wu C1, Tang L3, Zhou Z4. Fluorochloridone perturbs blood-testis barrier/Sertoli cell barrier function through Arp3-mediated F-actin disruption. *Toxicol Lett* 2018 Oct 1;295:277-287.
153. Liu L, Zhang Y, Chang X, Li R, Wu C, Tang L, Zhou Z. Fluorochloridone perturbs blood-testis barrier/Sertoli cell barrier function through Arp3-mediated F-actin disruption. *Toxicol Lett* 2018 Oct 1;295:277-287.
154. Liu L, Zhang Y, Chang X, Li R, Wu C, Tang L, Zhou Z. Fluorochloridone perturbs blood-testis barrier/Sertoli cell barrier function through Arp3-mediated F-actin disruption. *Toxicol Lett* 2018 Oct 1;295:277-287.
155. Liu L, Zhang Y, Chang X, Li R, Wu C, Tang L, Zhou Z. Fluorochloridone perturbs blood-testis barrier/Sertoli cell barrier function through Arp3-mediated F-actin disruption. *Toxicol Lett* 2018 Oct 1;295:277-287.
156. Liu L, Zhang Y, Chang X, Li R, Wu C, Tang L, Zhou Z. Fluorochloridone perturbs blood-testis barrier/Sertoli cell barrier function through Arp3-mediated F-actin disruption. *Toxicol Lett* 2018 Oct 1;295:277-287.
157. Liu L1, Zhang Y1, Chang X1, Li R2, Wu C1, Tang L3, Zhou Z4. Fluorochloridone perturbs blood-testis barrier/Sertoli cell barrier function through Arp3-mediated F-actin disruption. *Toxicol Lett* 2018 Oct 1;295:277-287.
158. Lin L, Wang L, Liu Y, Xu C, Tu Y, Zhou J. Non-thermal plasma inhibits tumor

- growth and proliferation and enhances the sensitivity to radiation in vitro and in vivo. *Oncol Rep* 2018 Dec;40(6):3405-3415.
159. Wu K, Mu XY, Jiang JT, Tan MY, Wang RJ, Zhou WJ, Wang X, He YY, Li MQ, Liu ZH miRNA-26a-5p and miR-26b-5p inhibit the proliferation of bladder cancer cells by regulating PDCD10. *Oncol Rep* 2018 Dec;40(6):3523-3532.
160. Wu K1, Mu XY2, Jiang JT2, Tan MY2, Wang RJ2, Zhou WJ3, Wang X2, He YY4, Li MQ5, Liu ZH2. miRNA-26a-5p and miR-26b-5p inhibit the proliferation of bladder cancer cells by regulating PDCD10. *Oncol Rep* 2018 Dec;40(6):3523-3532.
161. Wu K1, Mu XY2, Jiang JT2, Tan MY2, Wang RJ2, Zhou WJ3, Wang X2, He YY4, Li MQ5, Liu ZH2. miRNA-26a-5p and miR-26b-5p inhibit the proliferation of bladder cancer cells by regulating PDCD10. *Oncol Rep* 2018 Dec;40(6):3523-3532.
162. Wu K1, Mu XY2, Jiang JT2, Tan MY2, Wang RJ2, Zhou WJ3, Wang X2, He YY4, Li MQ5, Liu ZH2. miRNA-26a-5p and miR-26b-5p inhibit the proliferation of bladder cancer cells by regulating PDCD10. *Oncol Rep* 2018 Dec;40(6):3523-3532.
163. Wu K1, Mu XY2, Jiang JT2, Tan MY2, Wang RJ2, Zhou WJ3, Wang X2, He YY4, Li MQ5, Liu ZH2. miRNA-26a-5p and miR-26b-5p inhibit the proliferation of bladder cancer cells by regulating PDCD10. *Oncol Rep* 2018 Dec;40(6):3523-3532.
164. Lin L, Wang L, Liu Y, Xu C, Tu Y, Zhou J Non-thermal plasma inhibits tumor growth and proliferation and enhances the sensitivity to radiation in vitro and in vivo. *Oncol Rep* 2018 Dec;40(6):3405-3415.
165. Lin L, Wang L, Liu Y, Xu C, Tu Y, Zhou J Non-thermal plasma inhibits tumor growth and proliferation and enhances the sensitivity to radiation in vitro and in vivo. *Oncol Rep* 2018 Dec;40(6):3405-3415.
166. Wu K, Mu XY, Jiang JT, Tan MY, Wang RJ, Zhou WJ, Wang X, He YY, Li MQ, Liu ZH miRNA-26a-5p and miR-26b-5p inhibit the proliferation of bladder cancer cells by regulating PDCD10. *Oncol Rep* 2018 Dec;40(6):3523-3532.
167. Wu K, Mu XY, Jiang JT, Tan MY, Wang RJ, Zhou WJ, Wang X, He YY, Li MQ, Liu ZH miRNA-26a-5p and miR-26b-5p inhibit the proliferation of bladder cancer cells by regulating PDCD10. *Oncol Rep* 2018 Dec;40(6):3523-3532.
168. Lin L, Wang L, Liu Y, Xu C, Tu Y, Zhou J Non-thermal plasma inhibits tumor growth and proliferation and enhances the sensitivity to radiation in vitro and in vivo. *Oncol Rep* 2018 Dec;40(6):3405-3415.
169. Wu K, Mu XY, Jiang JT, Tan MY, Wang RJ, Zhou WJ, Wang X, He YY, Li MQ, Liu ZH miRNA-26a-5p and miR-26b-5p inhibit the proliferation of bladder cancer cells by regulating PDCD10. *Oncol Rep* 2018 Dec;40(6):3523-3532.
170. Zhang W1, Chen W1, Li Z1, Ma L1, Yu J1, Wang H1, Liu Z1, Xu B1. Identification and Characterization of Three New Cytochrome P450 Genes and the Use of RNA Interference to Evaluate Their Roles in Antioxidant Defense in *Apis cerana cerana* Fabricius. *Front Physiol* 2018 Nov 15;9:1608.
- 注: 更多使用本产品的文献请参考产品网页。

Version 2024.03.12