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- Yunpeng Zhao, Cheng Qiu, Wenhan Wang, Jiangfan Peng, Xiang Cheng, Yangtao Shanguan, Mingyang Xu, Jiayi Li, Ruize Qu, Xiaomin Chen, Suyi Jia, Dan Luo, Long Liu, Peng Li, Fengjin Guo, Krasimir Vasilev, Liang Liu, John Hayball, Shuli Dong, Xin Pan, Yuhua Li, Linlin Guo, Lei Cheng, Weiwei Li . Cortistatin protects against intervertebral disc degeneration through targeting mitochondrial ROS-dependent NLRP3 inflammasome activation Theranostics. 2020 May 27;10(15):7015-7033.
- Dong Liang, Chen Yu, Zhao Ma, Xingye Yang, Zhenzhen Li, Xuhui Dong, Xiaojun Qin, Lupei Du, Minyong Li . Identification of anthelmintic parbendazole as a therapeutic molecule for HNSCC through connectivity map-based drug repositioning Acta Pharm Sin B. 2022 May;12(5):2429-2442.
- Haofeng Zheng, Yunnan Zhang, Lei Li, Rui Zhang, Zihuan Luo, Zhe Yang, Yongrong Ye, Jiannan He, Qiquan Sun . Depletion of Toll-Like Receptor-9 Attenuates Renal Tubulointerstitial Fibrosis After Ischemia-Reperfusion Injury Front Cell Dev Biol. 2021 Feb 12;9:641527.
- Bing Liu, Qiang-Ming Li, Zhen-Zi Shang, Xue-Qiang Zha, Li-Hua Pan, Jian-Ping Luo . Anti-gastric cancer activity of cultivated Dendrobium huoshanense stem polysaccharide in tumor-bearing mice: Effects of molecular weight and O-acetyl group Int J Biol Macromol. 2021 Dec 1;192:590-599.
- Samantha E C Wolff, Xiao-Lan Wang, Han Jiao, Jia Sun, Andries Kalsbeek, Chun-Xia Yi, Yuanqing Gao . The Effect of Rev-erba Agonist SR9011 on the Immune Response and Cell Metabolism of Microglia Front Immunol. 2020 Sep 25;11:550145.
- Haofeng Zheng, Yunnan Zhang, Jiannan He, Zhe Yang, Rui Zhang, Lei Li, Zihuan Luo, Yongrong Ye, Qiquan Sun . Hydroxychloroquine Inhibits Macrophage Activation and Attenuates Renal Fibrosis After Ischemia-Reperfusion Injury Front Immunol. 2021 Apr 14;12:645100.
- Yufang Hou, Rixin Zhang, Jinbao Zong, Weiqi Wang, Mingxuan Zhou, Zheng Yan, Tiegang Li, Wenqiang Gan, Silin Lv, Zifan Zeng, Min Yang . Comprehensive Analysis of a Cancer-Immunity Cycle-Based Signature for Predicting Prognosis and Immunotherapy Response in Patients With Colorectal Cancer Front Immunol. 2022 May 31;13:892512.
- Mingxuan Zhou, Silin Lv, Yufang Hou, Rixin Zhang, Weiqi Wang, Zheng Yan, Tiegang Li, Wenqiang Gan, Zifan Zeng, Fang Zhang, Min Yang . Characterization of sialylation-related long noncoding RNAs to develop a novel signature for predicting prognosis, immune landscape, and chemotherapy response in colorectal cancer Front Immunol. 2022 Oct 18;13:994874.
- Chao Song, Biao Yu, Junjun Wang, Xinmiao Ji, Lei Zhang, Xiaofei Tian, Xin Yu, Chuanlin Feng, Xinyu Wang, Xin Zhang . Moderate Static Magnet Fields Suppress Ovarian Cancer Metastasis via ROS-Mediated Oxidative Stress Oxid Med Cell Longev. 2021 Dec 7;2021:7103345.
- Si-Lin Lv, Zi-Fan Zeng, Wen-Qiang Gan, Wei-Qi Wang, Tie-Gang Li, Yu-Fang Hou, Zheng Yan, Ri-Xin Zhang, Min Yang . Lp-PLA2 inhibition prevents Ang II-induced cardiac inflammation and fibrosis by blocking macrophage NLRP3 inflammasome activation Acta Pharmacol Sin. 2021 Dec;42(12):2016-2032.
- Li Z, Wang Y . Long non-coding RNA FTH1P3 promotes the metastasis and aggressiveness of non-small cell lung carcinoma by inducing epithelial-mesenchymal transition. INT J CLIN EXP PATHO. 2019 Oct 1 12(10):3782-3790. eCollection 2019.
- Kaiyan Lou, Ping Huang, Huijuan Ma, Xiaolei Wang, Huan Xu, Wei Wang . Orlistat increases arsenite tolerance in THP-1 derived macrophages through the up-regulation of ABCA1 Drug Chem Toxicol. 2019 Oct 31;1:9.
- Yu W, Yu Z, Wu D, Zhang J, Zhu Y, Zhang Y, Ning H, Wang M, Zhang J, Zhao J . Lens-specific conditional knockout of Msx2 in mice leads to ocular anterior segment dysgenesis via activation of a calcium signaling pathway. Lab Invest. 2019 Nov 99(11):1714-1727.
- Zhang J, Zheng Y, Luo Y, Du Y, Zhang X, Fu J . Curcumin inhibits LPS-induced neuroinflammation by promoting microglial M2 polarization via TREM2/TLR4/ NF-κB pathways in BV2 cells. Mol Immunol. 2019 Dec 116:29-37.
- Jia J, Nie L, Liu Y . Butyrate alleviates inflammatory response and NF-κB activation in human degenerated intervertebral disc tissues. Int Immunopharmacol. 2020 Jan 78:106004.
- Jian Pu, Jianchu Wang, Wenchuan Li, Yuan Lu, Xianjian Wu, Xidai Long, Chunying Luo, Huamei Wei . hsa_circ_0000092 promotes hepatocellular carcinoma progression through up-regulating HN1 expression by binding to microRNA-338-3p J Cell Mol Med. 2020 Feb 20. doi: 10.1111/jcmm.15010.
- Xiaolei Wang, Kaiyan Lou, Xiaodong Song, Huijuan Ma, Xinyi Zhou, Huan Xu, Wei Wang . Mebendazole is a potent inhibitor to chemoresistant T cell acute lymphoblastic leukemia cells TOXICOL APPL PHARM . 2020 Jun 1;396:115001.
- Lin Song, Jiaqi Cui, Ning Wang, Rui Wang, Jianqun Yan, Bo Sun . Maternal exercise during gestation and lactation decreases high-fat diet preference by altering central reward system gene expression in adult female offspring from high-fat fed dams Behav Brain Res. 2020 Jul 15;390:112660.
- Wenhan Wang, Ruitong Yang, Minfa Zhang, Jiayi Li, Jiangfan Peng, Mingyang Xu, Yunpeng Zhao, Hao Li, Xin Pan . Glycitin Suppresses Cartilage Destruction of Osteoarthritis in Mice Inflammation. 2020 Aug;43(4):1312-1322.
- Yufang Hou, Weiqi Wang, Zifan Zeng, Wenqiang Gan, Silin Lv, Tiegang Li, Zheng Yan, Rixin Zhang, Min Yang . High SEMA4C expression promotes the epithelial-mesenchymal transition and predicts poor prognosis in colorectal carcinoma AGING-US. 2020 Nov 7;12(21):21992-22018.
- Chao Tian, Shuai Li, Lang He, Xiaobo Han, Feng Tang, Rongqi Huang, Zuoxian Lin, Sihao Deng, Junjie Xu, Hualin Huang, Huifang Zhao, Zhiyuan Li .

- Transient receptor potential ankyrin 1 contributes to the lysophosphatidylcholine-induced oxidative stress and cytotoxicity in OLN-93 oligodendrocyte CELL STRESS CHAPERON. 2020 Nov;25(6):955-968.
22. Rixin Zhang, Tiegang Li, Weiqi Wang, Wenqiang Gan, Silin Lv, Zifan Zeng, Yufang Hou, Zheng Yan, Min Yang . Indoleamine 2, 3-Dioxygenase 1 and CD8 Expression Profiling Revealed an Immunological Subtype of Colon Cancer With a Poor Prognosis Front Oncol. 2020 Dec 7;10:594098.
 23. Xiaoqun Lv, Jinguo Zhang, Jun Zhang, Wencai Guan, Weifang Ren, Yujuan Liu, Guoxiong Xu . A Negative Feedback Loop Between NAMPT and TGF- β Signaling Pathway in Colorectal Cancer Cells Onco Targets Ther. 2021 Jan 8;14:187-198.
 24. Wei Wei, Xiaoguang Wei, Minyu Zhang, Cheng Peng . Ultrasound microbubble-mediated miR-150-5p inhibits gastric cancer cell growth by targeting the expression of NR2F2 J RECEPT SIG TRANSD. 2021 Feb 12;1-10.
 25. Yulin Mo, Wenxiu Lai, Ying Zhong, Zhuoqing Hu, Meigui You, Minqun Du, Pan Wang, Xinyou Wu, Cailing Chen, Huanmin He, Zhimin Gao, Yaping Xu, Dongtao Wang, Liao Cui, Yajun Yang . TXNIP contributes to bone loss via promoting the mitochondrial oxidative phosphorylation during glucocorticoid-induced osteoporosis Life Sci. 2021 Feb 1;266:118938.
 26. Xin Liu, Mingyang Xiao, Liang Zhang, Liuli Li, Guolian Zhu, Erdong Shen, Mingyue Lv, Xiaobo Lu, Zhe Sun . The m6A methyltransferase METTL14 inhibits the proliferation, migration, and invasion of gastric cancer by regulating the PI3K/AKT/mTOR signaling pathway J Clin Lab Anal. 2021 Mar;35(3):e23655.
 27. Shuang Yao, Fan Zheng, Yang Yu, Yuxia Zhan, Ning Xu, Guanghua Luo, Lu Zheng . Apolipoprotein M promotes cholesterol uptake and efflux from mouse macrophages FEBS Open Bio. 2021 Jun;11(6):1607-1620.
 28. Jianwei Cui, Chengping Li, Xiao Cui, Xueyan Liu, Chaoqun Meng, Guoli Zhou . Shortening of HO1 3'UTRs by Alternative Polyadenylation Suppresses Adipogenesis in 3T3-L1 J Agric Food Chem. 2021 Jul 21;69(28):8038-8049.
 29. Zhenwan Li, Jin Ma, Zhongsheng Kuang, Yong Jiang . β -Asarone Attenuates A β -Induced Neuronal Damage in PC12 Cells Overexpressing APPswe by Restoring Autophagic Flux Front Pharmacol. 2021 Jul 28;12:701635.
 30. Xuelu Chen, Shengfu Li, Dan Long, Juan Shan, Youping Li . Rapamycin facilitates differentiation of regulatory T cells via enhancement of oxidative phosphorylation Cell Immunol. 2021 Jul;365:104378.
 31. Rong Li, Sihao Hou, Ming Zou, Kunqi Ye, Li Xiang . miR-543 impairs cell proliferation, migration, and invasion in breast cancer by suppressing VCAN Biochem Biophys Res Commun. 2021 Sep 17;570:191-198.
 32. Jun-Hua Nie, Tao Yang, Hong Li, Hai-Shan Ye, Guo-Qing Zhong, Ting-Ting Li, Chi Zhang, Wen-Han Huang, Jin Xiao, Zhi Li, Jian-Li He, Bo-Le Du, Yu Zhang, Jia Liu . Identification of GPC3 mutation and upregulation in a multidrug resistant osteosarcoma and its spheroids as therapeutic target J Bone Oncol. 2021 Sep 20;30:100391.
 33. Jinqi Wang, Xuanhao Zeng, Yating Liu, Weiling Lian, Haozhen Lv, Kelu Wei, Qi Zhang, Jinhua Xu . Human embryonic stem cell-derived melanocytes exhibit limited immunogenicity Biochem Biophys Res Commun. 2021 Oct 8;573:151-157.
 34. Zhi-Yun Liu, Chun-Guang Zhao . miR-532-3p inhibits the progression of tongue squamous cell carcinoma by targeting podoplanin Chin Med J (Engl). 2021 Dec 8;134(24):2999-3008.
 35. Wenyun Li, Zhuo Sun, Min Wu, Zequn Deng, Mengman Zheng, Zhichao Kuang, Yuwei Liu, Gengsheng He . Deficiency of coiled-coil domain containing 80 increases plasma cholesterol by decreasing fecal sterols excretion in hypercholesterolemic mice J Nutr Biochem. 2021 Dec;98:108868.
 36. Chen Wang, Dongpeng Zheng, Fanglin Weng, Yongzeng Jin, Ling He . Sodium butyrate ameliorates the cognitive impairment of Alzheimer's disease by regulating the metabolism of astrocytes Psychopharmacology (Berl). 2022 Jan;239(1):215-227.
 37. Zheng Pang, Yan Xu, Qingjun Zhu . Early Growth Response 1 Suppresses Macrophage Phagocytosis by Inhibiting NRF2 Activation Through Upregulation of Autophagy During Pseudomonas aeruginosa Infection Front Cell Infect Microbiol. 2022 Jan 12;11:773665.
 38. Chen Zhao, Wenyuan Pu, Mengyuan Niu, Junaid Wazir, Shiyu Song, Lulu Wei, Li Li, Zhonglan Su, Hongwei Wang . Respiratory exposure to PM2.5 soluble extract induced chronic lung injury by disturbing the phagocytosis function of macrophage Environ Sci Pollut Res Int. 2022 Feb;29(10):13983-13997.
 39. Tiegang Li, Weiqi Wang, Wenqiang Gan, Silin Lv, Zifan Zeng, Yufang Hou, Zheng Yan, Rixin Zhang, Min Yang . Comprehensive bioinformatics analysis identifies LAPTM5 as a potential blood biomarker for hypertensive patients with left ventricular hypertrophy Aging (Albany NY). 2022 Feb 14;14(3):1508-1528.
 40. Ling Li, Tanggang Deng, Lin Zhang, Yanpeng Wang, Yangyang Zhou, Yan Liu, Hui Li, Jing Dai, Yuan Yang, Neng Ling, Huimin Liu, Jing Liu, Lanqin Cao, Min Xu, Mao Ye . ERK-mediated Cytoplasmic Retention of USP11 Contributes to Breast Cancer Cell Proliferation by Stabilizing Cytoplasmic p21 Int J Biol Sci. 2022 Mar 21;18(6):2568-2582.
 41. Hao Qi, Lulu Wang, Minghui Zhang, Zhe Wang, Xuejun Gao, Meng Li . Methionine and leucine induce ARID1A degradation to promote mTOR expression and milk synthesis in mammary epithelial cells J Nutr Biochem. 2022 Mar;101:108924.
 42. Qi Wu, Zhiyuan Zhang, Meiling Ji, Tao Yan, Yudong Jiang, Yijiao Chen, Jiang Chang, Jicheng Zhang, Dong Tang, Dexiang Zhu, Ye Wei . The Establishment and Experimental Verification of an lncRNA-Derived CD8+ T Cell Infiltration ceRNA Network in Colorectal Cancer Clin Med Insights Oncol. 2022 Apr 22;16:11795549221092218.
 43. Hai Zhu, Wen Jiang, Haixing Zhu, Jinwei Hu, Bingge Tang, Zhiqiang Zhou, Xinyang He . Elevation of ADAM12 facilitates tumor progression by enhancing metastasis and immune infiltration in gastric cancer Int J Oncol. 2022 May;60(5):51.
 44. Xi-Chen Wu, Zhi-Heng Zhu, Jun-Peng Zhang, Fu-Ming Shao, Jing-Qiu Peng, Yan Chen, Xue-Zong Wang, Wen-Yao Li, Yue-Long Cao, Wei Feng, Jian-Guang Xu, Dao-Fang Ding . Identification of thrombin as a key regulator of chondrocyte catabolic activity through RNA-Seq and experimental verification Gene. 2022 May 20;823:146327.
 45. Xue Wang, Chaofeng Yu, Xiaomei Liu, Jiasong Yang, Yuliang Feng, Yajun Wu, Yali Xu, Yihua Zhu, Wensheng Li . Fenofibrate Ameliorated Systemic and Retinal Inflammation and Modulated Gut Microbiota in High-Fat Diet-Induced Mice Front Cell Infect Microbiol. 2022 Jun 2;12:839592.
 46. Xunshan Ren, Bin Li, Changgeng Xu, Huangming Zhuang, Tianrun Lei, Fuze Jiang, Panghu Zhou . High expression of Piezo1 induces senescence in

chondrocytes through calcium ions accumulation Biochem Biophys Res Commun. 2022 Jun 4;607:138-145.

47. Fangting Liu, Zhehao Shi, Wenming Bao, Jiuyi Zheng, Kaiyu Chen, Zhihui Lin, Hao-Nan Song, Xin Luo, Qiantong Dong, Lei Jiang, Yi Wang, Gang Chen, Xiaolei Chen . ZIC2 promotes colorectal cancer growth and metastasis through the TGF- β signaling pathway Exp Cell Res. 2022 Jun 15;415(2):113118.
48. Wenxin Wang, Jie Geng, Xiaohan Wu, Jianguang Zhang, Chenna Zheng, Huachun Rao, Tianyu Li, Yong Diao, Huiyong Yang . Preparation of a miR-155-activating nucleic acid nanoflower to study the molecular mechanism of miR-155 in inflammation Mol Med. 2022 Jun 17;28(1):66.
49. Jiajun Chen, Jixin Qin, Jin Liu . Elucidation of the mechanism of miR-122-5p in mediating FOXO3 injury and apoptosis of mouse cochlear hair cells induced by hydrogen peroxide Exp Ther Med. 2022 Jun;23(6):435.
50. Lu-Yao Tong, Yong-Bing Deng, Wei-Hong Du, Wen-Zhu Zhou, Xin-Yu Liao, Xue Jiang . Clemastine Promotes Differentiation of Oligodendrocyte Progenitor Cells Through the Activation of ERK1/2 via Muscarinic Receptors After Spinal Cord Injury Front Pharmacol. 2022 Jul 5;13:914153.
51. Ya-Hui Xu, Yi-Fan Feng, Rong Zou, Fei Yuan, Yuan-Zhi Yuan . Silencing of YAP attenuates pericyte-myofibroblast transition and subretinal fibrosis in experimental model of choroidal neovascularization Cell Biol Int. 2022 Aug;46(8):1249-1263.
52. Lin Song, Jiaqi Cui, Shuyuan Hu, Rui Wang, Hongbao Li, Bo Sun . Maternal Treatment with Metformin Persistently Ameliorates High-Fat Diet-Induced Metabolic Symptoms and Modulates Gut Microbiota in Rat Offspring Nutrients. 2022 Sep 1;14(17):3612.
53. Haocun Kong, Luxi Yu, Caiming Li, Xiaofeng Ban, Zhengbiao Gu, Zhaofeng Li . Short-Clustered Maltodextrin Activates Ileal Glucose-Sensing and Induces Glucagon-like Peptide 1 Secretion to Ameliorate Glucose Homeostasis in Type 2 Diabetic Mice J Agric Food Chem. 2022 Oct 5;70(39):12604-12619.
54. Yuting Li, Ziqiong Yang, Xiangfen Li, Yang Yu, Xiaofeng Li, Peng Chen, Bing Li, Xiaoxiao Wang, Shou-Dong Ye . Prdm14 promotes mouse ESC self-renewal and PGCLC specification through enhancement of Stat3 activity iScience. 2022 Oct 6;25(11):105293.
55. Qi Hao, Lulu Wang, Minghui Zhang, Zhe Wang, Meng Li, Xuejun Gao . Taurine stimulates protein synthesis and proliferation of C2C12 myoblast cells through the PI3K-ARID4B-mTOR pathway Br J Nutr. 2022 Nov 28;128(10):1875-1886.
56. Meng-Di Lu, Hong Li, Jun-Hua Nie, Sheng Li, Hai-Shan Ye, Ting-Ting Li, Mo-Li Wu, Jia Liu . Dual Inhibition of BRAF-MAPK and STAT3 Signaling Pathways in Resveratrol-Suppressed Anaplastic Thyroid Cancer Cells with BRAF Mutations Int J Mol Sci. 2022 Nov 19;23(22):14385.
57. Yin Zhao, Shengyan Huang, Xirong Tan, Liufen Long, Qingmei He, Xiaoyu Liang, Jiewen Bai, Qingjie Li, Jiayi Lin, Yingqin Li, Na Liu, Jun Ma, Yupei Chen . N6 -Methyladenosine-Modified CBX1 Regulates Nasopharyngeal Carcinoma Progression Through Heterochromatin Formation and STAT1 Activation Adv Sci (Weinh). 2022 Dec;9(36):e2205091.

Version 2024.03.12