



碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301 或 800-8283301
订货 e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

β-半乳糖苷酶报告基因检测试剂盒

产品编号	产品名称	包装
RG0036	β-半乳糖苷酶报告基因检测试剂盒	200次

产品简介：

- β-半乳糖苷酶报告基因检测试剂盒(β-galactosidase Assay Kit)是一种用于β-半乳糖苷酶报告基因检测的试剂盒。β-半乳糖苷酶是一种常用的报告基因分子，经常与荧光素酶报告基因一起转染细胞，被用作荧光素酶报告基因检测的内参照，从而消除因为质粒的转染效率不同而带来的误差。
- 碧云天生产的β-半乳糖苷酶报告基因检测试剂盒，是以无色的ONPG为底物，在β-半乳糖苷酶的催化下生成黄色的onitrophenol，然后通过酶标仪或分光光度计在420nm波长附近测定吸光度，从而实现对β-半乳糖苷酶活性的测定。
- 本试剂盒所提供的裂解液适合于动物细胞的裂解，通常不适合于细菌、真菌和植物细胞的裂解。
- 如果使用96孔板检测，足够测定200个样品。

包装清单：

产品编号	产品名称	包装
RG0036-1	报告基因细胞裂解液	100ml
RG0036-2	β-半乳糖苷酶检测试剂	10ml
RG0036-3	β-半乳糖苷酶反应终止液	35ml
—	说明书	1份

保存条件：

报告基因细胞裂解液4°C保存3个月有效，-20°C保存至少一年有效。β-半乳糖苷酶检测试剂-20°C避光保存至少6个月有效，-80°C避光保存至少一年有效。β-半乳糖苷酶反应终止液可以室温保存，也可以4°C或-20°C保存。

注意事项：

- 最佳测定波长为420nm。在410-430nm范围内测定，具有相近的灵敏度。
- 本检测试剂盒和Promega荧光素酶报告基因的裂解液兼容。用Promega荧光素酶报告基因裂解液裂解的样品可以用本试剂盒测定。
- 本产品仅限于专业人员的科学的研究用，不得用于临床诊断或治疗，不得用于食品或药品，不得存放于普通住宅内。
- 为了您的安全和健康，请穿实验服并戴一次性手套操作。

使用说明：

1. 裂解细胞：按如下方式加入报告基因细胞裂解液，充分裂解细胞。细胞裂解后可以立即测定β-半乳糖苷酶活性，也可以先冻存，待以后再测定。冻存的样品需融解，并达到室温后才可以开始测定。建议使用12、24或48孔板进行报告基因检测。

器皿类型	96孔板	48孔板	24孔板	12孔板	6孔板
细胞裂解液 (微升/孔)	100	150	200	300	500

2. 融解β-半乳糖苷酶检测试剂，并达到室温。充分融解β-半乳糖苷酶反应终止液，如有难溶结晶，请在37°C水浴溶解。
3. 在96孔板中，分别加入5-50微升样品，然后加入报告基因裂解液至最后体积为50微升。由于不同的细胞体系，不同的转染方法的转染效率相差很大。在第一次测定时推荐取两个样品作预实验。在96孔板中，分别加入0微升，5微升，10微升，25微升，50微升样品，然后加入报告基因裂解液至最后体积为50微升。
4. 在每个孔中加入β-半乳糖苷酶检测试剂50微升，混合后，盖上96孔板的盖子。可以用parafilm或保鲜膜封住96孔板防止蒸发。在37°C放置30分钟或直至样品孔内出现浅黄色。通常大约3小时后，吸光度达到平台，孵育更长时间吸光度不会明显升高。
5. 加入150微升β-半乳糖苷酶反应终止液终止反应，混匀，尽量避免气泡。
6. 按仪器操作说明书开启酶标仪或分光光度计，将测定波长设定为420nm，测定吸光度。

使用本产品的文献：

1. Qi L, Yue L, Feng D, Qi F, Li J, Dong X. . Genome-wide mRNA processing in methanogenic archaea reveals post-transcriptional regulation of ribosomal protein synthesis. Nucleic Acids Res. 2017 Jul;45(12):7285-7298.
2. Zu-Guo Zheng, Si-Tong Zhu, Hui-Min Cheng, Xin Zhang, Gang Cheng,

- Pyone Myat Thu, Supeng Perry Wang, Hui-Jun Li, Ming Ding, Lei Qiang, Xiao-Wei Chen, Qing Zhong, Ping Li, Xiaojun Xu . Discovery of a potent SCAP degrader that ameliorates HFD-induced obesity, hyperlipidemia and insulin resistance via an autophagy-independent lysosomal pathway Autophagy. 2021 Jul;17(7):1592-1613.

3. Zhou Y, Xiang Z, Li D, Han X. . Regulation of microcystin-LR-induced toxicity in mouse spermatogonia by miR-96. *Environ Sci Technol*. 2014 Jun 3;48(11):6383-90.
4. Jin Wang, Jiayi Xie, Daosong Wang, Xue Han, Minqi Chen, Guojun Shi, Linjia Jiang, Meng Zhao . CXCR4high megakaryocytes regulate host-defense immunity against bacterial pathogens Elife. 2022 Jul 29:11:e78662.
5. Wang W, Vellaisamy K, Li G, Wu C, Ko CN, Leung CH, Ma DL. . Development of a Long-Lived Luminescence Probe for Visualizing β -Galactosidase in OvarianCarcinoma Cells. *Anal Chem*. 2017 Nov 7;89(21):11679-11684.
6. Lijun Zheng, Chenyu Luo, Nan Yang, Hailong Pei, Mintao Ji, Yinyin Shu, Zhisen Zhang, Shuai Dong, Xiuxiu Wang, Xuemei Li, Wensheng Zhang, Yan Wang, Bingyan Li, Tom K Hei, Guangming Zhou, Lei Chang . Ionizing radiation-induced long noncoding RNA CRYBG3 regulates YAP/TAZ through mechanotransduction *Cell Death Dis*. 2022 Mar 4;13(3):209.
7. Ren X, Chen Z, Liu Y, Zhang H, Zhang M, Liu Q, Hong X, Zhu JK, Gong Z . ABO3, a WRKY transcription factor, mediates plant responses to abscisic acid and droughttolerance in *Arabidopsis*. *Plant J*. 2010 Aug;63(3):417-29.
8. Liu Y, Liu W, Hu C, Xue Z, Wang G, Ding B, Luo H, Tang L, Kong X, Chen X, Liu N, Ding Y, Jin Y . MiR-17 modulates osteogenic differentiation through a coherent feed-forward loop inmesenchymal stem cells isolated from periodontal ligaments of patients with periodontitis. *Stem Cells*. 2011 Nov;29(11):1804-16.
9. Wang D, Liu C, Wang Y, Wang W, Wang K, Wu X, Li Z, Zhao C, Li L, Peng L . Impact of miR-26b on cardiomyocyte differentiation in P19 cells through regulating canonical/non-canonical Wnt signalling. *CELL PROLIFERAT*. 2017 Dec;50(6).
10. Chen H, Sun JG, Cao XW, Ma XG, Xu JP, Luo FK, Chen ZT . Preliminary validation of ERBB2 expression regulated by miR-548d-3p and miR-559. *BIOCHEM BIOPH RES CO*. 2009 Aug 7;385(4):596-600.
11. Yang Z, Zhang Y, Yang Y, Sun L, Han D, Li H, Wang C . Pharmacological and toxicological target organelles and safe use of single-walled carbonnanotubes as drug carriers in treating Alzheimer disease. *Nanomedicine*. 2010 Jun;6(3):427-41.
12. Yin K, Deng X, Mo ZC, Zhao GJ, Jiang J, Cui LB, Tan CZ, Wen GB, Fu Y, Tang CK . Tristetraprolin-dependent post-transcriptional regulation of inflammatory cytokine mRNAexpression by apolipoprotein A-I: role of ATP-binding membrane cassette transporter A1 and signal transducer and activator of transcription 3. *J Biol Chem*. 2011 Apr 22;286(16):13834-45.
13. Cao X, Yang M, Wei RC, Zeng Y, Gu JF, Huang WD, Yang DQ, Li HL, Ding M, Wei N, Zhang KJ, Xu B, Liu XR, Qian QJ, Liu XY . Cancer targeting Gene-Viro-Therapy of liver carcinoma by dual-regulated oncolytic adenovirusarmed with TRAIL gene. *Gene Ther*. 2011 Aug;18(8):765-77.
14. Luo M, Weng Y, Tang J, Hu M, Liu Q, Jiang F, Yang D, Liu C, Zhan X, Song P, Bai H, Li B, Shi Q . MicroRNA-450a-3p represses cell proliferation and regulates embryo development by regulating Bub1 expressionin mouse. *PLoS One*. 2012;7(10):e47914.
15. Li L, Shi JY, Zhu GQ, Shi B . MiR-17-92 cluster regulates cell proliferation and collagen synthesis by targeting TGFB pathway inmouse palatal mesenchymal cells. *J Cell Biochem*. 2012 Apr;113(4):1235-44.
16. Yin J, Sun L, Dong Y, Chi X, Zhu W, Qi SH, Gao H . Expression of blaA underlies unexpected ampicillin-induced cell lysis of *Shewanella oneidensis*. *PLoS One*. 2013;8(3):e60460.
17. Zhang H, Fu H, Wang J, Sun L, Jiang Y, Zhang L, Gao H . Impacts of nitrate and nitrite on physiology of *Shewanella oneidensis*. *PLoS One*. 2013 Apr 23;8(4):e62629.
18. Shi M, Wu L, Xia Y, Chen H, Luo Q, Sun L, Gao H . Exoprotein production correlates with morphotype changes of nonmotile *Shewanella oneidensis* mutants. *J Bacteriol*. 2013 Apr;195(7):1463-74.
19. Zhou L, Zheng H, Tang Y, Yu W, Gong Q . Eugenol inhibits quorum sensing at sub-inhibitory concentrations. *Biotechnol Lett*. 2013 Apr;35(4):631-7.
20. Bai R, Yi S, Zhang X, Liu H, Fang X . Role of ICAM-1 polymorphisms (G241R, K469E) in mediating its single-molecule binding ability: atomic force microscopy measurements on living cells. *BIOCHEM BIOPH RES CO*. 2014 Jun 13;448(4):372-8.
21. Shi M, He X, Wei W, Wang J, Zhang T, Shen X . Tenascin-C induces resistance to apoptosis in pancreatic cancer cell through activation of ERK/NF- κ B pathway. *Apoptosis*. 2015 Jun;20(6):843-57.
22. Cheng S, Liang X, Wang Y, Jiang Z, Liu Y, Hou W, Li S, Zhang J, Wang Z . The circadian Clock gene regulates acrosin activity of sperm through serine protease inhibitor A3K. *EXP BIOL MED*. 2015 Jun;20(6):843-57.
23. Li L, Zheng M, Ma H, Gong S, Ai G, Liu X, Li J, Wang K, Dong X . Significant performance enhancement of a UASB reactor by using acyl homoserine lactones to facilitate the long filaments of *Methanosaeta harundinacea* 6Ac. *APPL MICROBIOL BIOT*. 2015 Aug;99(15):6471-80.
24. Zhang R, Weng Y, Li B, Jiang Y, Yan S, He F, Chen X, Deng F, Wang J, Shi Q . BMP9-induced osteogenic differentiation is partially inhibited by miR-30a in the mesenchymal stem cell line C3H10T1/2. *J Mol Histol*. 2015 Oct;46(4-5):399-407.
25. Zhang R, Yan S, Wang J, Deng F, Guo Y, Li Y, Fan M, Song Q, Liu H, Weng Y, Shi Q . MiR-30a regulates the proliferation, migration, and invasion of human osteosarcoma by targeting Runx2. *TUMOR BIOL* . 2016 Mar;37(3):3479-88.
26. Wang XJ, Wang J, Wang YY, Guo YJ, Chu BB, Yang GY . Sus scrofa matrix attachment region enhances expression of the PiggyBac system transfected intoHEK293T cells. *Biotechnol Lett*. 2016 Jun;38(6):949-58.
27. Wang J, Xin B, Wang H, He X, Wei W, Zhang T, Shen X . Gastrin regulates ABCG2 to promote the migration, invasion and side populations in pancreatic cancercells via activation of NF- κ B signaling. *Exp Cell Res*. 2016 Aug 1;346(1):74-84.
28. Xin B, He X, Wang J, Cai J, Wei W, Zhang T, Shen X . Nerve growth factor regulates CD133 function to promote tumor cell migration and invasion via activatingERK1/2 signaling in pancreatic cancer. *Pancreatology*. 2016 Nov - Dec;16(6):1005-14.
29. Wang H, Cai J, Du S, Guo Z, Xin B, Wang J, Wei W, Shen X . Fractalkine/CX3CR1 induces apoptosis resistance and proliferation through the activation of the AKT/NF- κ B cascade in pancreatic cancer cells. *Cell Biochem Funct*. 2017 Aug;35(6):315-326.
30. Yang N, Li Y, Wang G, Ding Y, Jin Y, Xu Y . Tumor necrosis factor-a suppresses adipogenic and osteogenic differentiation of humanperiodontal ligament stem cell by inhibiting miR-21/Spry1 functional axis. *Differentiation*. 2017 Sep - Oct;97:33-43.
31. Liu H, Li G, Zhang B, Sun D, Wu J, Chen F, Kong F, Luan Y, Jiang W, Wang R, Xue X . Suppression of the NF- κ B signaling pathway in colon cancer cells by the natural compound Riccardin D from Dumortierahirsute. *Mol Med Rep*. 2018 Apr;17(4):5837-5843.
32. Li S, Sun W, Zheng H, Tian F . Microrna-145 accelerates the inflammatory reaction through activation of NF- κ B signaling in atherosclerosis cells and mice. *Biomed Pharmacother*. 2018 Jul;103:851-857.
33. Liu Z, Hu K, Tang Y, Li H, Tang H, Hu X, Ma X, Liu Z . SmpB

- down-regulates proton-motive force for the persister tolerance to aminoglycosides in *Aeromonas veronii*. *BIOCHEM BIOPH RES CO.* 2018 Dec 9;507(1-4):407-413.
34. Miaomia Chi, Ying Liu, Hongling Wen, Li Zhao, Yanyan Song, Na Liu, Lianli Chi, Zhiyu Wang . Role of the head-stalk linker region (aa115-122) of Newcastle disease virus haemagglutinin-neuraminidase in regulating fusion triggering J Gen Virol. 2020 Jan;101(1):5-20.
 35. Ying Liu, Miaomia Chi, Hongling Wen, Li Zhao, Yanyan Song, Na Liu, Lianli Chi, Zhiyu Wang . The DI-DII linker of human parainfluenza virus type 3 fusion protein is critical for the virus Virus Genes. 2020 Feb;56(1):37-48.
 36. Xiaoxue Shen, Zhibin Zhao, Bin Yang . MicroRNA-155 promotes apoptosis of colonic smooth muscle cells and aggravates colonic dysmotility by targeting IGF-1 Exp Ther Med. 2020 Apr;19(4):2725-2732.
 37. Feng Wang, Ke Fan, Ying Zhao, Mei-Lin Xie . Apigenin attenuates TGF-β1-stimulated cardiac fibroblast differentiation and extracellular matrix production by targeting miR-155-5p/c-Ski/Smad pathway J Ethnopharmacol. 2021 Jan 30;265:113195.
 38. Yanan Huang, Yaqing Liu, Yanguo Li, Ying Liu, Chi Zhang, Hongling Wen, Li Zhao, Yanyan Song, Liyang Wang, Zhiyu Wang . Role of key amino acids in the transmembrane domain of the Newcastle disease virus fusion protein Biosci Trends. 2021 Mar 15;15(1):16-23.
 39. Wang Feng, Zhao Ying, Fan Ke, Xie Mei-Lin . Apigenin suppresses TGF-β1-induced cardiac fibroblast differentiation and collagen synthesis through the downregulation of HIF-1α expression by miR-122-5p Phytomedicine. 2021 Mar;83:153481.
 40. Feng Wang, Jun Zhang, Guanghao Niu, Jiayi Weng, Qian Zhang, Meilin Xie, Chunjian Li, Kangyun Sun . Apigenin inhibits isoproterenol-induced myocardial fibrosis and Smad pathway in mice by regulating oxidative stress and miR-122-5p/155-5p expressions Drug Dev Res. 2022 Jun;83(4):1003-1015.
 41. Yihui Zhang, Kang Zhan, Zixuan Hu, Guoqi Zhao . SV40T/E6E7-Induced Proliferation Is Involved in the Activity of E2F3 in Bovine Mammary Epithelial Cells Animals (Basel). 2022 Jul 12;12(14):1790.
 42. Mengdi Yang, Dan Li, Zhiyuan Jiang, Changcan Li, Suyuan Ji, Jing Sun, Yujie Chang, Shunyi Ruan, Zhiyu Wang, Rui Liang, Xueyu Dai, Bin Li, Hui Zhao . TGF-β-Induced FLRT3 Attenuation Is Essential for Cancer-Associated Fibroblast-Mediated Epithelial-Mesenchymal Transition in Colorectal Cancer Mol Cancer Res. 2022 Aug 5;20(8):1247-1259.
 43. Guangzhi Li, Dongmei Zhang, Shizhen Wang, Ni Jie, Yongting Qin . Protective Role of Taraxasterol against Cardiovascular Aging and Aging-Induced Desensitization of Insulin Signaling Front Biosci (Landmark Ed). 2022 Nov 22;27(11):311.

Version 2024.03.12