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1. Tao Wan, Jiafeng Zhong, Qi Pan, Tianhua Zhou, Yuan Ping, Xiangrui Liu . Exosome-mediated delivery of Cas9 ribonucleoprotein complexes for tissue-specific gene therapy of liver diseases Sci Adv. 2022 Sep 16;8(37):eabp9435.
2. Pu L, Geng Y, Liu S, Chen J, Luo K, Wang G, Gu Z. . Electroneutralized amphiphilic triblock copolymer with a peptide dendron for efficient muscular gene delivery. ACS APPL MATER INTER. 2014 Sep 10;6(17):15344-51.
3. Nianyin Lv, Yufang Zhao, Xiaoyi Liu, Lusha Ye, Zihao Liang, Yanhua Kang, Yeping Dong, Wei Wang, Narasaiah Kolliputi, Liyun Shi . Dysfunctional telomeres through mitostress-induced cGAS/STING activation to aggravate immune senescence and viral pneumonia Aging Cell. 2022 Apr;21(4):e13594.
4. Lei R, Zhang K, Liu K, Shao X, Ding Z, Wang F, Hong Y, Zhu M, Li H, Li H. . Transferrin receptor facilitates TGF-β and BMP signaling activation to control craniofacial morphogenesis. Cell Death Dis. 2016 Jun 30;7(6):e2282.
5. Liping Chu, Yuxiu Qu, Yang An, Linjun Hou, Juewan Li, Weijia Li, Gaofeng Fan, Bao-Liang Song, En Li, Liye Zhang, Wei Qi . Induction of senescence-associated secretory phenotype underlies the therapeutic efficacy of PRC2 inhibition in cancer Cell Death Dis. 2022 Feb 15;13(2):155.
6. Hongwei Liao, Xiang Li, Lianzheng Zhao, Yalong Wang, Xiaodan Wang, Ye Wu, Xin Zhou, Wei Fu, Lei Liu, Hong-Gang Hu, Ye-Guang Chen . A PROTAC peptide induces durable β-catenin degradation and suppresses Wnt-dependent intestinal cancer Cell Discov. 2020 Jun 9;6:35.
7. Shen T, Li S, Cai LD, Liu JL, Wang CY, Gan WJ, Li XM, Wang JR, Sun LN, Deng M, Liu YH, Li JM . Erbin exerts a protective effect against inflammatory bowel disease by suppressing autophagic cell death. ONCOTARGET. 2018 Jan 4;9(15):12035-12049.
8. Zhang L, Gao X, Men K, Wang B, Zhang S, Qiu J, Huang M, Gou M, Huang N, Qian Z, Zhao X, Wei Y. . Gene therapy for C-26 colon cancer using heparin-polyethyleneimine nanoparticle-mediatedsurvivin T34A. INT J NANOMED. 2011;6:2419-27.
9. Zhen-Hua Chen, Si Li, Mingrui Xu, Candace C Liu, Hongying Ye, Ben Wang, Qing-Feng Wu . Single-cell transcriptomic profiling of the hypothalamic median eminence during aging J Genet Genomics. 2022 Jun;49(6):523-536.

10. Zhang Q, Yu K, Yuan W, Xue F, Sun H, Zhu H. . Purification of large-scale plasmid DNAs by selective precipitation with cetyltrimethylammonium bromide. *生物工程学报*. 2008 Dec;24(12):2117-21.
11. Li J, Zeng W, Huang Y, Zhang Q, Hu P, Rabkin SD, Liu R. . Treatment of breast cancer stem cells with oncolytic herpes simplex virus. *Cancer Gene Ther*. 2012 Oct;19(10):707-1.
12. Wang XJ, Li Y, Huang H, Zhang XJ, Xie PW, Hu W, Li DD, Wang SQ. . A simple and robust vector-based shRNA expression system used for RNA interference. *PLoS One*. 2013;8(2):e56110.
13. Wang J, Dan G, Shangguan T, Hao H, Tang R, Peng K, Zhao J, Sun H, Zou Z. . miR-198 Represses the Proliferation of HaCaT Cells by Targeting Cyclin D2. *Int J Mol Sci*. 2015 Jul 27;16(8):17018-28.
14. Tan W, Gu Z, Shen B, Jiang J, Meng Y, Da Z, Liu H, Tao T, Cheng C. . PTEN/Akt - p27kip1 signaling promote the BM - MSCs senescence and apoptosis in SLE patients. *J Cell Biochem*. 2015 Aug;116(8):1583-94.
15. Zhang Y, Huang X, Zhao XY, Hu YJ, Sun HY, Kong WJ. . Role of the Ubiquitin C-Terminal Hydrolase L1-Modulated Ubiquitin Proteasome System in Auditory Cortex Senescence. *ORL J OTO-RHINO-LARY*. 2017;79(3):153-163.
16. Gu C, Li Y, Liu J, Ying X, Liu Y, Yan J, Chen C, Zhou H, Cao L, Ma Y. . LncRNA mediated SIRT1/FoxO3a and SIRT1/p53 signaling pathways regulate type II alveolar epithelial cell senescence in patients with chronic obstructive pulmonary disease. *Mol Med Rep*. 2017 May;15(5):3129-3134.
17. Cheng L, Jiang H, Fan J, Wang J, Hu P, Ruan Y, Liu R. . A novel oncolytic herpes simplex virus armed with the carboxyl-terminus of murine MyD116 has enhanced anti-tumour efficacy against human breast cancer cells. *Oncol Lett*. 2018 May;15(5):7046-7052.
18. Li Z, Liu T, Ma J, Guo Q, Ma L, Lv Q, Jiang Y, Wei C, Zhang J. . TGF- β induces corneal endothelial senescence via increase of mitochondrial reactive oxygen species in chronic corneal allograft failure. *AGING-US*. 2018 Nov 28;10(11):3474-3485.
19. Hai-Long Zhang, Bing Zhao, Pin Yang, Yin-Quan Du, Wei Han, Jianming Xu, Dong-Min Yin. . Steroid Receptor Coactivator 3 Regulates Synaptic Plasticity and Hippocampus-dependent Memory *Neurosci Bull*. 2021 Dec;37(12):1645-1657.
20. Jia Ye, Zhihui Jin, Sen Chen, Weichun Guo. . Uridine relieves MSCs and chondrocyte senescence in vitro and exhibits the potential to treat osteoarthritis in vivo *Cell Cycle*. 2022 Jan;21(1):33-48.
21. Yanshuang Zhang, Tingyu Wu, Fan Li, Yirui Cheng, Qing Han, Xin Lu, Shun Lu, Weiliang Xia. . FGF19 Is Coamplified With CCND1 to Promote Proliferation in Lung Squamous Cell Carcinoma and Their Combined Inhibition Shows Improved Efficacy *Front Oncol*. 2022 Apr 7;12:846744.
22. Niran Fang, Bin Yang, Ting Xu, Yanpeng Li, Huimin Li, Hanghui Zheng, Aiguo Zhang, Ruiai Chen. . Expression and Immunogenicity of Recombinant African Swine Fever Virus Proteins Using the Semliki Forest Virus *Front Vet Sci*. 2022 May 9;9:870009.
23. Han Tang, Han Du, Xielan Kuang, Hao Huang, Jingshu Zeng, Chongde Long, Binbin Zhu, Licheng Fu, Hua Wang, Qingjiong Zhang, Shuibin Lin, Jianhua Yan, Huangxuan Shen. . Arbutin Protects Retinal Pigment Epithelium Against Oxidative Stress by Modulating SIRT1/FOXO3a/PGC-1 α/β Pathway *Front Genet*. 2022 Aug 16;13:922807.
24. Pingcui Xu, Bingfeng Lin, Xuehui Deng, Kai Huang, Yan Zhang, Nani Wang. . VDR activation attenuates osteoblastic ferroptosis and senescence by stimulating the Nrf2/GPX4 pathway in age-related osteoporosis *Free Radic Biol Med*. 2022 Nov 20;193(Pt 2):720-735.

Version 2024.03.12