



碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301 或 800-8283301
 订货 e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

增强型ATP检测试剂盒

产品编号	产品名称	包装
S0027	增强型ATP检测试剂盒	200次

产品简介:

- 增强型ATP检测试剂盒(Enhanced ATP Assay Kit)可以用于检测普通溶液、细胞或组织内的ATP(adenosine 5'-triphosphate)水平。细胞和组织样品一步裂解即可完成样品制备,检测灵敏度高达0.1nM,化学发光可以持续稳定达30分钟。本试剂盒的多项指标优于许多国外的同类产品。并且获得的样品还可以同时进行Western检测。
- ATP,作为最重要的能量分子,在细胞的各种生理、病理过程中起着重要作用。ATP水平的改变,会影响细胞的功能。通常细胞在凋亡、坏死或处于一些毒性状态下,ATP水平会下降;而高葡萄糖刺激等可上调某些细胞的细胞内ATP水平。通常ATP水平的下降表明线粒体的功能受损或下降,在细胞凋亡时ATP水平的下降通常和线粒体的膜电位下降同时发生。
- 本试剂盒根据萤火虫萤光素酶(firefly luciferase,也称荧光素酶)催化萤光素产生萤光时需要ATP提供能量研制而成。当萤火虫萤光素酶和萤光素都过量时,在一定的浓度范围内萤光的产生和ATP的浓度成正比。这样就可以高灵敏地检测溶液中的ATP浓度。
- **样品制备简单。**本试剂盒提供了可以用于细胞和组织裂解的ATP检测裂解液,简单裂解后即可用于ATP检测。无需进行高氯酸或三氯乙酸(TCA)抽提,或样品裂解后的煮沸等繁琐操作。
- **灵敏度高,线性范围宽,在0.1nM至10 μ M范围内有良好检测效果。**本试剂盒在样品体积为100微升时可以检测浓度低达0.1nM的ATP。而常规的细胞或组织裂解液中ATP的浓度为0.1-1 μ M,一些常见细胞的细胞内ATP水平约为10nmol/mg蛋白。并且本试剂盒的检测浓度范围非常宽,检测上限可以高达10 μ M,并在0.1nM-10 μ M范围内可以形成良好的标准曲线。
- **读数稳定。**本试剂盒进行了特殊的优化设计,使检测ATP时的化学发光非常稳定。对于ATP标准曲线的检测结果显示,在开始反应后10分钟内化学发光无明显下降,开始反应后30分钟内化学发光的下降不超过10%。
- **制备的样品兼容性好。**使用本试剂盒中的ATP检测裂解液裂解获得的细胞或组织样品,不仅可以用于ATP检测,还可用于检测蛋白浓度、进行SDS-PAGE或一些常规的较易溶解蛋白的Western检测。
- **使用方便快捷。**通常10-20个样品可以在30-60分钟内测定完毕。
- 碧云天的三款ATP检测试剂盒的主要特点和差异如下:

产品编号	S0026	S0026B	S0027
产品名称	ATP检测试剂盒	ATP检测试剂盒	增强型ATP检测试剂盒
检测灵敏度	++++	+++	+++++
特殊样品兼容性	++++	+++++	++++
信号值	++++	+++	+++++
信号稳定性	+++++	++++	+++++
检测下限(ATP)	1nM	5nM	0.1nM
检测上限(ATP)	10 μ M	10 μ M	10 μ M
线性范围(ATP)	1nM-10 μ M	10nM-10 μ M	0.1nM-10 μ M
推荐指数	++++	+++	+++++

- 一个包装的本试剂盒至少可以检测200个样品。

包装清单:

产品编号	产品名称	包装
S0027-1	ATP检测试剂	4ml(1ml/管,共4管)
S0027-2	ATP检测试剂稀释液	20ml
S0027-3	ATP标准溶液(0.5mM)	0.1ml
S0027-4	ATP检测裂解液	100ml
—	说明书	1份

保存条件:

-20 $^{\circ}$ C保存,半年有效。-70 $^{\circ}$ C保存,一年有效。ATP检测试剂需避光保存。

注意事项:

- ATP检测试剂中含有萤光素酶，反复冻融会导致其逐渐失活。尽管经测试ATP检测试剂反复冻融5次对于其检测效果无明显影响，为取得良好的使用效果，建议用户使用时的冻融次数不宜超过3次。ATP检测试剂稀释成ATP检测工作液后，最好一次用完，不宜冻存后再使用。
- ATP，特别是裂解后样品中的ATP在室温不太稳定，需在4°C或冰上操作。ATP在冰上可以稳定长达6小时。
- 本试剂盒需使用luminometer，即化学发光仪(检测萤光素酶报告基因时所用的仪器)。如果没有luminometer，也可以使用液闪仪。液闪仪的测定效果取决于液闪仪的检测灵敏度和检测精度。
- 使用可检测化学发光的多功能酶标仪时，推荐使用孔和孔之间不透光的96孔白板或黑板。如使用普通的透明96孔板，须特别注意在检测孔之间设置间隔孔，以减少邻近孔之间的相互干扰。对于透明96孔板，一个发光孔可以使上下或左右邻近孔的RLU值升高该孔的10-20%左右，使上下或左右间隔一个孔的邻近孔的RLU值升高该孔的1%-4%左右；对于相同的样品，底部不透光的96孔白板的化学发光读数可以达到透明96孔板的5-10倍左右，达到底部透光96孔白板读数的3倍左右(实测数据会因96孔板、检测仪器和样品等的不同而存在差异)。
- 本试剂盒提供的ATP检测裂解液可以有效裂解并释放常见的培养细胞和组织中的ATP。对于一些特殊的组织或样品，如果发现检测出来的ATP水平显著低于预期水平，可以在裂解样品后并且在离心前，取部分样品煮沸2分钟以充分释放ATP。煮沸后样品中的蛋白会变性，从而会在后续的离心步骤中被沉淀，因此煮沸的样品不能用于蛋白浓度测定、SDS-PAGE和Western检测。可以使用剩余的部分样品进行蛋白浓度测定、SDS-PAGE和Western检测。
- 本产品仅限于专业人员的科学研究用，不得用于临床诊断或治疗，不得用于食品或药品，不得存放于普通住宅内。
- 为了您的安全和健康，请穿实验服并戴一次性手套操作。

使用说明：

1. 样品测定的准备：(注意：样品裂解需在4°C或冰上操作)

对于贴壁细胞：

吸除培养液，按照6孔板每孔加入200微升裂解液的比例(即相当于细胞培养液量2毫升的1/10)加入裂解液，裂解细胞。裂解细胞时为了裂解充分，可以使用移液器进行反复吹打或晃动培养板使裂解液充分接触并裂解细胞。通常细胞在接触裂解液后会立即裂解。裂解后4°C 12000g离心5分钟，取上清，用于后续的测定。

对于悬浮细胞：

用离心管离心沉淀细胞，弃上清，轻轻弹散细胞，按照6孔板每孔的细胞量加入200微升裂解液的比例加入裂解液，裂解细胞。裂解细胞时为了裂解充分可以弹击离心管管底或适当Vortex使裂解液充分接触并裂解细胞。通常细胞在接触裂解液后会立即裂解。裂解后4°C 12000g离心5分钟，取上清，用于后续的测定。

对于组织样品：

按照每20毫克组织加入约100-200微升裂解液的比例加入裂解液，然后用玻璃匀浆器或其它匀浆设备进行匀浆。充分匀浆可以确保组织被完全裂解。裂解后4°C 12000g离心5分钟，取上清，用于后续的测定。

2. 标准曲线测定的准备：

冰浴上融解待用试剂，把ATP标准溶液用ATP检测裂解液稀释成适当的浓度梯度。具体的浓度需根据样品中ATP的浓度而定。初次检测可以检测0.01、0.03、0.1、0.3、1、3和10 μ M这几个浓度，在后续的实验中，可以根据样品中ATP的浓度对标准品的浓度范围进行适当调整。

3. ATP检测工作液的配制：

按照每个样品或标准品需100微升ATP检测工作液的比例配制适当量的ATP检测工作液。把待用试剂在冰浴上融解。取适量的ATP检测试剂，按照1:4的比例用ATP检测试剂稀释液稀释ATP检测试剂。例如1毫升ATP检测试剂加入4毫升ATP检测试剂稀释液配制成5毫升ATP检测工作液。稀释后的ATP检测试剂即为用于后续实验的ATP检测工作液。ATP检测工作液可在冰浴上暂时保存。

4. ATP浓度的测定：

- a. 加100微升ATP检测工作液到检测孔或检测管内。室温放置3-5分钟，以使本底性的ATP全部被消耗掉，从而降低本底。可以一次性把10-20个检测孔或检测管分别加上100微升ATP检测工作液，从而节省时间。
- b. 在检测孔或检测管内加上20微升样品或标准品，迅速用枪(微量移液器)混匀，至少间隔2秒后，用化学发光仪(luminometer)或液闪仪测定RLU值或CPM。(注：样品的体积可以自行在10-100微升范围内调节。如果样品中的ATP浓度比较低则可以加入100微升样品，如果样品中ATP浓度比较高则可以加入较小体积的样品，同时标准品也需要使用相同的体积。如果样品中ATP的浓度特别高，可以用ATP检测裂解液稀释样品后再测定。本试剂盒在加入10-100微升标准品时，在0.1nM-10 μ M的浓度范围内有很好的线性关系(参考图1)。




图1. 本产品对ATP标准品的检测效果。图中数据为20微升标准品并减去空白对照后的数据。实测数据会因检测仪器等的不同而存在差异，图中数据仅供参考。

- c. 根据标准曲线计算出样品中ATP的浓度。
- d. 为了消除样品制备时由于蛋白量的差异而造成的误差，可以用碧云天生产的BCA蛋白浓度测定试剂盒(P0009/P0010/P0010S/P0011/P0012/P0012S)测定样品中的蛋白浓度。然后把ATP的浓度换算成nmol/mg蛋白的形式。

常见问题:

1. Luminometer和荧光分光光度计有何不同?

荧光分光光度计检测的样品本身不能发光，样品需要由特定波长的激发光激发，然后才能产生荧光并被荧光分光光度计检测。Luminometer，即化学发光检测仪，检测的样品本身可以发光，不需要激发光进行激发。也就是说Luminometer是检测化学发光(萤光)的仪器。有些型号的荧光分光光度计也具有luminometer的功能，即也可以检测化学发光。您所使用的荧光分光光度计能否用于化学发光的测定请仔细阅读该仪器的说明书。

2. 可以进行萤光素酶报告基因检测的仪器是否就可以用于本试剂盒的ATP检测?

是。萤光素酶报告基因的检测原理和本增强型ATP检测试剂盒的原理相同，可以用相同的仪器测定，并且可以选择相同的测定参数，例如检测前等待时间为2秒，检测时间为10秒等。

3. 多功能酶标仪是否可以用于本试剂盒的ATP检测?

不一定。如果该多功能酶标仪具有检测化学发光的功能，即具有luminometer的功能，就可以用于本试剂盒的检测，否则就不能了。

相关产品:

产品编号	产品名称	包装
S0026	ATP检测试剂盒	200次
S0026B	ATP检测试剂盒	200次
S0027	增强型ATP检测试剂盒	200次

使用本产品的文献:

1. Zhou J, Zhang F, Meng H, Zhang Y, Li Y. Introducing extra NADPH consumption ability significantly increases the photosynthetic efficiency and biomass production of cyanobacteria. *Metab Eng* 2016 Nov;38:217-227.
2. Xue RQ, Sun L, Yu XJ, Li DL, Zang WJ. Vagal nerve stimulation improves mitochondrial dynamics via an M3 receptor/CaMKKβ/AMPK pathway in isoproterenol-induced myocardial ischaemia. *J Cell Mol Med* 2017 Jan;21(1):58-71.
3. Zhang J, Cai Q, Jiang M, Liu Y, Gu H, Guo J, Sun H, Fang J, Jin L. *Exp Gerontol* 2017 Mar;89:45-56.
4. Li Y, Ye Z, Lai W, Rao J, Huang W, Zhang X, Yao Z, Lou T. *Front Pharmacol* 2017 Apr 5;8:178.
5. Zhou W, Tian D, He J, Zhang L, Tang X, Zhang L, Wang Y, Li L, Zhao J, Yuan X, Peng S. *Environ Pollut* 2017 Jul;226:412-425.
6. Zhou W, Tian D, He J, Zhang L, Tang X, Zhang L, Wang Y, Li L, Zhao J, Yuan X, Peng S. Exposure scenario: Another important factor determining the toxic effects of PM2.5 and possible mechanisms involved. *Environ Pollut* 2017 Jul;226:412-425.
7. Zhang J, Cai Q, Jiang M, Liu Y, Gu H, Guo J, Sun H, Fang J, Jin L. Mesencephalic astrocyte-derived neurotrophic factor alleviated 6-OHDA-induced cell damage via ROS-AMPK/mTOR mediated autophagic inhibition. *Exp Gerontol* 2017 Mar;89:45-56.
8. Xue RQ, Sun L, Yu XJ, Li DL, Zang WJ. Vagal nerve stimulation improves mitochondrial dynamics via an M3 receptor/CaMKKβ/AMPK pathway in isoproterenol-induced myocardial ischaemia. *J Cell Mol Med* 2017 Jan;21(1):58-71.
9. Li HJ, Sun XM, Li ZK, Yin QW, Pang H, Pan JJ, Li X, Chen W. LncRNA UCA1 Promotes Mitochondrial Function of Bladder Cancer via the MiR-195/ARL2 Signaling Pathway. *CELL PHYSIOL BIOCHEM* 2017;43(6):2548-2561.
10. Li Y, Ye Z, Lai W, Rao J, Huang W, Zhang X, Yao Z, Lou T. Activation of Sirtuin 3 by Silybin Attenuates Mitochondrial Dysfunction in Cisplatin-induced Acute Kidney Injury. *Front Pharmacol* 2017 Apr 5;8:178.
11. Zhao R, Lv M, Li Y, Sun M, Kong W, Wang L, Song S, Fan C, Jia L, Qiu S, Sun Y, Song H, Hao R. Stable Nanocomposite Based on PEGylated and Silver Nanoparticles Loaded Graphene Oxide for Long-Term Antibacterial Activity. *ACS APPL MATER INTER* 2017 May 10;9(18):15328-15341.
12. Chen Y, Liang P, Huang Y, Li M, Zhang X, Ding C, Feng J, Zhang Z, Zhang X, Gao Y, Zhang Q, Cao S, Zheng H, Liu D, Songyang Z, Huang J. Glycerol kinase-like proteins cooperate with Pld6 in regulating sperm mitochondrial sheath formation and male fertility. *Cell Discov* 2017 Aug 22;3:17030.

13. Long X, Pan Q, Wang C, Wang H, Li H, Li X. Microbial fuel cell-photoelectrocatalytic cell combined system for the removal of azo dyewastewater. *BIORESOURCE TECHNOL* 2017 Nov;244(Pt 1):182-191.
14. Zhou DC, Su YH, Jiang FQ, Xia JB, Wu HY, Chang ZS, Peng WT, Song GH, Park KS, Kim SK, Cai DQ, Zheng L, Qi XF. CpG oligodeoxynucleotide preconditioning improves cardiac function after myocardial infarction via modulation of energy metabolism and angiogenesis. *J Cell Physiol* 2018 May;233(5):4245-4257.
15. Cheng C, Wang T, Song Z, Peng L, Gao M, Hermine O, Rousseaux S, Khochbin S, Mi JQ, Wang J. Induction of autophagy and autophagy-dependent apoptosis in diffuse large B-cell lymphoma by a new antimalarial artemisinin derivative, SM1044. *CANCER MED-US* 2018 Feb;7(2):380-396.
16. Li H, Yang XL, Song HL, Zhang S, Long XZ. Effects of direct current on *Klebsiella* spp. viability and corresponding resistance gene expression in simulative bio-electrochemical reactors. *Chemosphere* 2018 Apr;196:251-259.
17. Zhang Y, Wei J, Xu J, Leong WS, Liu G, Ji T, Cheng Z, Wang J, Lang J, Zhao Y, You L, Zhao X, Wei T, Anderson GJ, Qi S, Kong J, Nie G, Li S. Suppression of Tumor Energy Supply by Liposomal Nanoparticle-Mediated Inhibition of Aerobic Glycolysis. *ACS APPL MATER INTER* 2018 Jan 24;10(3):2347-2353.
18. Wang C, Li K, Li T, Chen Z, Wen Y, Liu X, Jia X, Zhang Y, Xu Y, Han M, Komatsu N, Zhao L, Chen X. Monocyte-mediated chemotherapy drug delivery in glioblastoma. *NANOMEDICINE-UK* 2018 Jan;13(2):157-178.
19. Suo Y, Fu H, Ren M, Liao Z, Ma Y, Wang J. Enhanced butyric acid production in *Clostridium tyrobutyricum* by overexpression of rate-limiting enzymes in the Embden-Meyerhof-Parnas pathway. *J Biotechnol* 2018 Apr 20;272-273:14-21.
20. Zhao H, Wang Q, Liu C, Shang Y, Wen F, Wang F, Liu W, Xiao W, Li W. A Role for the Respiratory Chain in Regulating Meiosis Initiation in *Saccharomyces cerevisiae*. *Genetics* 2018 Mar;208(3):1181-1194.
21. Hu G, Zhou J, Chen X, Qian Y, Gao C, Guo L, Xu P, Chen W, Chen J, Li Y, Liu L. Engineering synergetic CO₂-fixing pathways for malate production. *Metab Eng* 2018 May;47:496-504.
22. Zhou DC, Su YH, Jiang FQ, Xia JB, Wu HY, Chang ZS, Peng WT, Song GH, Park KS, Kim SK, Cai DQ, Zheng L, Qi XF. CpG oligodeoxynucleotide preconditioning improves cardiac function after myocardial infarction via modulation of energy metabolism and angiogenesis. *J Cell Physiol* 2018 May;233(5):4245-4257.
23. Li TF, Li K, Zhang Q, Wang C, Yue Y, Chen Z, Yuan SJ, Liu X, Wen Y, Han M, Komatsu N, Xu YH, Zhao L, Chen X. Dendritic cell-mediated delivery of doxorubicin-polyglycerol-nanodiamond composites elicits enhanced anti-cancer immune response in glioblastoma. *Biomaterials* 2018 Oct;181:35-52.
24. Ji F, Zhao C, Wang B, Tang Y, Miao Z, Wang Y. The role of 5-hydroxymethylcytosine in mitochondria after ischemic stroke. *J Neurosci Res* 2018 Oct;96(10):1717-1726.
25. Deyu H, Luqing C, Xianglian L, Pu G, Qirong L, Xu W, Zonghui Y. Protective mechanisms involving enhanced mitochondrial functions and mitophagy against T-2 toxin-induced toxicities in GH3 cells. *Toxicol Lett* 2018 Oct 1;295:41-53.
26. Zhou Y, Lian S, Zhang J, Lin D, Huang C, Liu L, Chen Z. Mitochondrial Perturbation Contributing to Cognitive Decline in Streptozotocin-Induced Type 1 Diabetic Rats. *CELL PHYSIOL BIOCHEM* 46(4):1668-1682. 2018
27. Zhao Y, Du R, Zhou T, Yang D, Huang Y, Wang Y, Huang J, Ma X, He F, Qiu J, Wang G. Arsenic Trioxide-Coated Stent Is an Endothelium-Friendly Drug Eluting Stent. *Adv Healthc Mater* 7(15):e1800207. 2018 Aug
28. Huang LL, Long RT, Jiang GP, Jiang X, Sun H, Guo H, Liao XH. Augmenter of liver regeneration promotes mitochondrial biogenesis in renal ischemia-reperfusion injury. *Apoptosis* 23(11-12):695-706. 2018 Dec
29. Cao K, Lai F, Zhao XL, Wei QX, Miao XY, Ge R, He QY, Sun X. The mechanism of iron-compensation for manganese deficiency of *Streptococcus pneumoniae*. *J Proteomics* 184:62-70. 2018 Jul 30
30. Dai J, Huang Q, Niu K, Wang B, Li Y, Dai C, Chen Z, Tao K, Dai J. Sestrin 2 confers primary resistance to sorafenib by simultaneously activating AKT and AMPK in hepatocellular carcinoma. *CANCER MED-US* 7(11):5691-5703. 2018 Nov
31. Middepogu A, Hou J, Gao X, Lin D. Effect and mechanism of TiO₂ nanoparticles on the photosynthesis of *Chlorella pyrenoidosa*. *ECOTOX ENVIRON SAFE* 161:497-506. 2018 Oct
32. Zhang Y, Liu Y, Bi X, Hu C, Ding F, Ding W. Therapeutic Approaches in Mitochondrial Dysfunction, Inflammation, and Autophagy in Uremic Cachexia: Role of Aerobic Exercise. *MEDIAT INFLAMM* 2019;2789014. 2019 Aug 18
33. Wang HW, Zhang Y, Tan PP, Jia LS, Chen Y, Zhou BH. Mitochondrial respiratory chain dysfunction mediated by ROS is a primary point of fluoride-induced damage in Hepa1-6 cells. *Environ Pollut* 255(Pt 3):113359. 2019 Dec
34. Meng Q, Peng X, Zhao S, Xu T, Wang S, Liu Q, Cai R, Fan Y. Hypoxic storage of erythrocytes slows down storage lesions and prolongs shelf-life. *J Cell Physiol* 234(12):22833-22844. 2019 Dec
35. Yu HF, Duan CC, Yang ZQ, Wang YS, Yue ZP, Guo B. HB-EGF Ameliorates Oxidative Stress-Mediated Uterine Decidualization Damage. *Oxid Med Cell Longev* 2019;6170936. 2019 Dec 2
36. Han T, Yan J, Chen H, Ji Y, Chen J, Cui J, Shen W, Zou J. HIF-1 α contributes to tube malformation of human lymphatic endothelial cells by upregulating VEGFR-3. *Int J Oncol* 54(1):139-151. 2019 Jan
37. Cheng Z, Zhang M, Hu J, Lin J, Feng X, Wang S, Wang T, Gao E, Wang H, Sun D. Cardiac-specific Mst1 deficiency inhibits ROS-mediated JNK signalling to alleviate Ang II-induced cardiomyocyte apoptosis. *J Cell Mol Med* 23(1):543-555. 2019 Jan
38. Liu H, Huang H, Li R, Bi W, Feng L, E L, Hu M, Wen W. Mitophagy protects SH-SY5Y neuroblastoma cells against the TNF α -induced inflammatory injury: Involvement of microRNA-145 and Bnip3. *Biomed Pharmacother* 109:957-968. 2019 Jan
39. Xiao Y, Zhang J, Shu X, Bai L, Xu W, Wang A, Chen A, Tu WY, Wang J, Zhang K, Luo B, Shen C. Loss of mitochondrial protein CHCHD10 in skeletal muscle causes neuromuscular junction impairment. *Hum Mol Genet* pii: ddz154. 2019 Jul 1
40. Li H, Wang C, He T, Zhao T, Chen YY, Shen YL, Zhang X, Wang LL. Mitochondrial Transfer from Bone Marrow Mesenchymal Stem Cells to Motor Neurons in Spinal Cord Injury Rats via Gap Junction. *Theranostics* 9(7):2017-2035. 2019 Mar 17
41. Zhao Y, Liu S, Zhou L, Li X, Meng Y, Li Y, Li L, Jiao B, Bai L, Yu Y, Zhang S, Li W, Hoffman AR, Hu JF, Cui J. Aberrant shuttling of long noncoding RNAs during the mitochondria-nuclear crosstalk in hepatocellular carcinoma cells. *Am J Cancer Res* 9(5):999-1008. eCollection 2019. 2019 May 1
42. Jiang Q, Zhang C, Wang H, Peng T, Zhang L, Wang Y, Han W, Shi C. Mitochondria-Targeting Immunogenic Cell Death Inducer Improves the Adoptive T-Cell Therapy Against Solid Tumor. *Front Oncol* 9:1196. 2019 Nov 12
43. Zhou TT, Zhao T, Ma F, Zhang YN, Jiang J, Ruan Y, Yan QY, Wang GH, Ren J, Guan XW, Guo J, Zhao YH, Ye JM, Hu LH, Chen J, Shen X. Small molecule IVQ, as a prodrug of gluconeogenesis inhibitor QVO, efficiently ameliorates glucose homeostasis in type 2 diabetic mice. *Acta Pharmacol Sin* 40(9):1193-1204. 2019 Sep
44. Zhang J, Si H, Li B, Zhou X, Zhang J. Myriscignan Exhibits Activities Against *Toxoplasma gondii* RH Strain by Triggering Mitochondrial Dysfunction. *Front Microbiol* 10:2152. 2019 Sep 18
45. Feng C, Li D, Chen M, Jiang L, Liu X, Li Q, Geng C, Sun X, Yang G, Zhang L, Yao X. Citreoviridin induces myocardial apoptosis through PPAR- γ -mTORC2-mediated autophagic pathway and the protective effect of thiamine and selenium. *CHEM-BIOL INTERACT* 311:108795. 2019 Sep 25
46. Chunyan Cheng, Tao Wang, Zhiqun Song, Lijun Peng, Mengqing Gao, Olivier Hermine, Sophie Rousseaux, Saadi Khochbin, Jian-Qing Mi, Jin Wang. Induction of Autophagy and Autophagy-Dependent Apoptosis in Diffuse Large B-cell Lymphoma by a New Antimalarial Artemisinin Derivative, SM1044. *CANCER MED-US* 2018 Feb;7(2):380-396.
47. Kui Zhu, Shang Chen, Tatyana A Syssoeva, Lingchong You. Universal Antibiotic Tolerance Arising From Antibiotic-Triggered Accumulation of Pyocyanin in

- Pseudomonas Aeruginosa* PLoS Biol 2019 Dec 16;17(12):e3000573.
48. Zhou T, Chang L, Luo Y, Zhou Y, Zhang J. Mst1 inhibition attenuates non-alcoholic fatty liver disease via reversing Parkin-related mitophagy. Redox Biol 21:101120. 2019 Feb
 49. Yuan Liu, Yuqian Jia, Kangni Yang, Ziwen Tong, Jingru Shi, Ruichao Li, Xia Xiao, Wenkai Ren, Rüdiger Hardeland, Russel J Reiter, Zhiqiang Wang Melatonin overcomes MCR-mediated colistin resistance in Gram-negative pathogens Theranostics 2020 Aug 29;10(23):10697-10711.
 50. Min Wang, Rui-Ying Wang, Jia-Hui Zhou, Xue-Heng Xie, Gui-Bo Sun, Xiao-Bo Sun Calenduloside E Ameliorates Myocardial Ischemia-Reperfusion Injury through Regulation of AMPK and Mitochondrial OPA1 Oxid Med Cell Longev 2020 Aug 31;2020:2415269.
 51. Binwu Hu, Shuo Zhang, Weijian Liu, Peng Wang, Songfeng Chen, Xiao Lv, Deyao Shi, Kaige Ma, Baichuan Wang, Yongchao Wu, Zengwu Shao Inhibiting Heat Shock Protein 90 Protects Nucleus Pulposus-Derived Stem/Progenitor Cells From Compression-Induced Necroptosis and Apoptosis Front Cell Dev Biol 2020 Aug 7;8:685.
 52. Deqin Kong, Rui Liu, Jiangzheng Liu, Qingbiao Zhou, Jiaxin Zhang, Wenli Li, Hua Bai, Chunxu Hai Cubic Membranes Formation in Synchronized Human Hepatocellular Carcinoma Cells Reveals a Possible Role as a Structural Antioxidant Defense System in Cell Cycle Progression Front Cell Dev Biol 2020 Dec 14;8:617406.
 53. Luyao Zhang, Zichuan Wang, Tengfei Lu, Lin Meng, Yan Luo, Xiangwei Fu, Yunpeng Hou Mitochondrial Ca²⁺ Overload Leads to Mitochondrial Oxidative Stress and Delayed Meiotic Resumption in Mouse Oocytes Front Cell Dev Biol 2020 Dec 15;8:580876.
 54. Meng Peng, Yuan Liu, Xiang-Qin Zhang, Ya-Wei Xu, Yin-Tao Zhao, Hai-Bo Yang CTRP5-Overexpression Attenuated Ischemia-Reperfusion Associated Heart Injuries and Improved Infarction Induced Heart Failure Front Pharmacol 2020 Dec 22;11:603322.
 55. Huiting Zhang, Ke Yan, Lumin Sui, Junyu Nie, Kexin Cui, Jiahao Liu, Hengye Zhang, Xiaogan Yang, Kehuan Lu, Xingwei Liang Constant light exposure causes oocyte meiotic defects and quality deterioration in mice Environ Pollut 2020 Dec;267:115467.
 56. Xiang Lin, Yongdong Dai, Xiaomei Tong, Wenzhi Xu, Qianmeng Huang, Xiaoying Jin, Chao Li, Feng Zhou, Hanjin Zhou, Xiaona Lin, Dong Huang, Songying Zhang Excessive oxidative stress in cumulus granulosa cells induced cell senescence contributes to endometriosis-associated infertility Redox Biol 2020 Feb;30:101431.
 57. Yanping Liu, Chunli Yang, Xueqin Feng, Linglu Qi, Jun Guo, Dan Zhu, Phung N Thai, Yingying Zhang, Pengjie Zhang, Miao Sun, Juanxiu Lv, Lubo Zhang, Zhice Xu, Xiyuan Lu Prenatal High-Salt Diet-Induced Metabolic Disorders via Decreasing Peroxisome Proliferator-Activated Receptor Gamma Coactivator 1 α in Adult Male Rat Offspring Mol Nutr Food Res 2020 Jul;64(14):e2000196.
 58. Shuang Nie, Xuettian Qian, Mengyue Shi, Hongzhen Li, Chunyan Peng, Xiwei Ding, Shu Zhang, Bin Zhang, Guifang Xu, Ying Lv, Lei Wang, Helmut Friess, Bo Kong, Xiaoping Zou, Shanshan Shen ALDH1A3 Accelerates Pancreatic Cancer Metastasis by Promoting Glucose Metabolism Front Oncol 2020 Jun 16;10:915.
 59. Chunyi Tong, Xianghua Zhong, Yuejun Yang, Xu Liu, Guowei Zhong, Chang Xiao, Bin Liu, Wei Wang, Xiaoping Yang PB@PDA@Ag nanosystem for synergistically eradicating MRSA and accelerating diabetic wound healing assisted with laser irradiation Biomaterials 2020 Jun;243:119936.
 60. Jing Wang, Huanying Deng, Jixiang Zhang, Dandan Wu, Jiao Li, Jingjing Ma, Weiguo Dong α -Hederin induces the apoptosis of gastric cancer cells accompanied by glutathione decrement and reactive oxygen species generation via activating mitochondrial dependent pathway Phytother Res 2020 Mar;34(3):601-611.
 61. Yonglian Zeng, Zhenya Guo, Zhigao Hu, Mingjiang Liu, Yubing Chen, Shilian Chen, Bo Peng, Peng Zhang, Zhan Wu, Hongliang Luo, Fudi Zhong, Keqing Jiang, Yi Lu, Guandou Yuan, Songqing He FGD1 exhibits oncogenic properties in hepatocellular carcinoma through regulating cell morphology, autophagy and mitochondrial function Biomed Pharmacother 2020 May;125:110029.
 62. Wei Yang, Shuai Liu, Yunlei Li, Yujie Wang, Yao Deng, Weimin Sun, Hualan Huang, Junmou Xie, Andong He, Honglv Chen, Ailin Tao, Jie Yan Pyridoxine induces monocyte-macrophages death as specific treatment of acute myeloid leukemia Cancer Lett 2020 Nov 1;492:96-105.
 63. Can Hu, Ge Zhou, Kun Liu, Wenjun Yin, Lingyun Zhou, Jianglin Wang, Linhua Chen, Shanru Zuo, Yueliang Xie, Xiacong Zuo CaMKII as a key regulator of contrast-induced nephropathy through mPTP opening in HK-2 cells Cell Signal 2020 Nov;75:109734.
 64. Xiuge Gao, Xiangchun Ruan, Hui Ji, Lin Peng, Yawei Qiu, Dan Yang, Xinhao Song, Chunlei Ji, Dawei Guo, Shanxiang Jiang Maduramicin triggers methuosis-like cell death in primary chicken myocardial cells Toxicol Lett 2020 Oct 15;333:105-114.
 65. Zhi-Teng Chen, Hai-Feng Zhang, Meng Wang, Shao-Hua Wang, Zhu-Zhi Wen, Qing-Yuan Gao, Mao-Xiong Wu, Wen-Hao Liu, Yong Xie, Jing-Ting Mai, Ying Yang, Jing-Feng Wang, Yang-Xin Chen Long non-coding RNA Linc00092 inhibits cardiac fibroblast activation by altering glycolysis in an ERK-dependent manner Cell Signal 2020 Oct;74:109708.
 66. Yi-Zhi Ren, Ben-Zheng Zhang, Xiao-Jing Zhao, Zhi-Yuan Zhang Resolvin D1 ameliorates cognitive impairment following traumatic brain injury via protecting astrocytic mitochondria J Neurochem 2020 Sep;154(5):530-546.
 67. Nanshan Song, Hong Zhu, Rong Xu, Jiaqi Liu, Yinquan Fang, Jing Zhang, Jianhua Ding, Gang Hu, Ming Lu Induced Expression of kir6.2 in A1 Astrocytes Propagates Inflammatory Neurodegeneration via Drp1-dependent Mitochondrial Fission Front Pharmacol 2021 Jan 28;11:618992.
 68. Zihao Dai, Zongren Wang, Kai Lei, Junbin Liao, Zhenwei Peng, Manxia Lin, Ping Liang, Jie Yu, Sui Peng, Shuling Chen, Ming Kuang Irreversible electroporation induces CD8⁺ T cell immune response against post-ablation hepatocellular carcinoma growth Cancer Lett 2021 Apr 10;503:1-10.
 69. Wenhua Kuang, Jie Zhang, Zhou Lan, R N V Krishna Deepak, Chao Liu, Zhilong Ma, Lili Cheng, Xinbin Zhao, Xianbin Meng, Weihua Wang, Xueying Wang, Lina Xu, Yupei Jiao, Qi Luo, Ziyi Meng, Kehkooi Kee, Xiaohui Liu, Haiteng Deng, Wei Li, Hao Fan, Ligong Chen SLC22A14 is a mitochondrial riboflavin transporter required for sperm oxidative phosphorylation and male fertility Cell Rep 2021 Apr 20;35(3):109025.
 70. Tuersunayi Abudurehman, Jing Xia, Ming-Hao Li, Hang Zhou, Wei-Wei Zheng, Neng Zhou, Rong-Yi Shi, Jian-Min Zhu, Li-Ting Yang, Li Chen, Liang Zheng, Kai Xue, Kai Qing, Cai-Wen Duan CDK7 Inhibitor THZ1 Induces the Cell Apoptosis of B-Cell Acute Lymphocytic Leukemia by Perturbing Cellular Metabolism Front Oncol 2021 Apr 6;11:663360.
 71. Xiaoyong Huang, Xiaole Zhao, Kui Zhu, Shuangyang Ding, Bing Shao Sodium dehydroacetate exposure decreases locomotor persistence and hypoxia tolerance in zebrafish Environ Res 2021 Apr;195:110276.
 72. Cuijuan Qian, Shihang Chen, Sen Li, Yichao Wang, Jun Yao Circ-0000003 regulates glutamine metabolism and tumor progression of tongue squamous cell carcinoma via the miR-330-3p/GLS axis Oncol Rep 2021 Apr;45(4):45.
 73. Jingwu Zhu, Jiang Tian, Chao Yang, Junpeng Chen, Lihuang Wu, Mengni Fan, Xiaojun Cai L-Arg-Rich Amphiphilic Dendritic Peptide as a Versatile NO Donor for NO/Photodynamic Synergistic Treatment of Bacterial Infections and Promoting Wound Healing Small 2021 Aug;17(32):e2101495.
 74. Youkun Bi, Xuejun Guo, Mengqi Zhang, Keqi Zhu, Chentao Shi, Baoqi Fan, Yanyun Wu, Zhiguang Yang, Guangju Ji Bone marrow derived-mesenchymal stem cell improves diabetes-associated fatty liver via mitochondria transformation in mice Stem Cell Res Ther 2021 Dec 11;12(1):602.
 75. Zhihong Zheng, Jude Juventus Aweya, Shiyuan Bao, Defu Yao, Shengkang Li, Ngoc Tuan Tran, Hongyu Ma, Yueling Zhang The Microbial Composition of Penaeid Shrimps' Hepatopancreas Is Modulated by Hemocyanin J Immunol 2021 Dec 1;207(11):2733-2743.
 76. Yun Cao, Zhaowei Chen, Jijia Hu, Jun Feng, Zijiang Zhu, Yanqin Fan, Qiaoxuan Lin,

- Guohua Ding Mfn2 Regulates High Glucose-Induced MAMs Dysfunction and Apoptosis in Podocytes via PERK Pathway *Front Cell Dev Biol* 2021 Dec 20;9:769213.
77. Fuli Zheng, Zhousong Luo, Xinpei Lin, Wei Wang, Michael Aschner, Ping Cai, Yuan-Liang Wang, Wenya Shao, Guangxia Yu, Zhenkun Guo, Siying Wu, Huangyuan Li Intercellular transfer of mitochondria via tunneling nanotubes protects against cobalt nanoparticle-induced neurotoxicity and mitochondrial damage *Nanotoxicology* 2021 Dec;15(10):1358-1379.
78. Meimei Zheng, Xu Chen, Yiqiang Cui, Wen Li, Haiqian Dai, Qiuling Yue, Hao Zhang, Ying Zheng, Xuejiang Guo, Hui Zhu TULP2, a New RNA-Binding Protein, Is Required for Mouse Spermatid Differentiation and Male Fertility *Front Cell Dev Biol* 2021 Feb 18;9:623738.
79. Wei Yang, Jianhua Liu, Lin Hou, Qingmin Chen, Yahui Liu Shikonin differentially regulates glucose metabolism via PKM2 and HIF1 α to overcome apoptosis in a refractory HCC cell line *Life Sci* 2021 Jan 15;265:118796.
80. Liang Sun, Xiaoqin Yang, Xiaoheng Huang, Yizhou Yao, Xiangyu Wei, Shugao Yang, Diyuan Zhou, Wei Zhang, Zhimin Long, Xiaoyan Xu, Xinguo Zhu, Songbing He, Xiong Su 2-Hydroxylation of Fatty Acids Represses Colorectal Tumorigenesis and Metastasis via the YAP Transcriptional Axis *Cancer Res* 2021 Jan 15;81(2):289-302.
81. Shiyuan Fu, Xinyu Zhao, Zhou Zhou, Mengyan Li, Liang Zhu Effective removal of odor substances using intimately coupled photocatalysis and biodegradation system prepared with the silane coupling agent (SCA)-enhanced TiO₂ coating method *Water Res* 2021 Jan 1;188:116569.
82. Xiaoqiu Zhu, Yiyi Yao, Mingyan Guo, Jin Li, Pengfeng Yang, Hui Xu, Daowei Lin Sevoflurane increases intracellular calcium to induce mitochondrial injury and neuroapoptosis *Toxicol Lett* 2021 Jan 1;336:11-20.
83. Yuan Liu, Kangni Yang, Yuqian Jia, Jingru Shi, Ziwen Tong, Zhiqiang Wang Thymine Sensitizes Gram-Negative Pathogens to Antibiotic Killing *Front Microbiol* 2021 Jan 28;12:622798.
84. Yawei Wang, Binlin Tang, Lei Long, Peng Luo, Wei Xiang, Xueru Li, Huilan Wang, Qingzhi Jiang, Xu Tan, Shenglin Luo, Huijuan Li, Ziwen Wang, Zelin Chen, Yu Leng, Zhongyong Jiang, Yang Wang, Le Ma, Rui Wang, Chunyu Zeng, Zujian Liu, Yu Wang, Hongming Miao, Chunmeng Shi Improvement of obesity-associated disorders by a small-molecule drug targeting mitochondria of adipose tissue macrophages *Nat Commun* 2021 Jan 4;12(1):102.
85. Yi Li, Yu-Fan Feng, Xiao-Tian Liu, Yu-Chen Li, Hui-Min Zhu, Meng-Ru Sun, Ping Li, Baolin Liu, Hua Yang Songorine promotes cardiac mitochondrial biogenesis via Nrf2 induction during sepsis *Redox Biol* 2021 Jan;38:101771.
86. Xing Liu, Lingyan Huang, Christopher Rensing, Jie Ye, Kenneth H Neelson, Shungui Zhou Syntrophic interspecies electron transfer drives carbon fixation and growth by *Rhodospseudomonas palustris* under dark, anoxic conditions *Sci Adv* 2021 Jul 2;7(27):eabh1852.
87. Yuan Liu, Kangni Yang, Yuqian Jia, Jingru Shi, Ziwen Tong, Dan Fang, Bingqing Yang, Chengrui Su, Ruichao Li, Xia Xiao, Zhiqiang Wang Gut microbiome alterations in high-fat-diet-fed mice are associated with antibiotic tolerance *Nat Microbiol* 2021 Jul;6(7):874-884.
88. Yu Zhang, Zhengzhe Zhang, Yinguang Chen Biochar Mitigates N₂O Emission of Microbial Denitrification through Modulating Carbon Metabolism and Allocation of Reducing Power *Environ Sci Technol* 2021 Jun 15;55(12):8068-8078.
89. Yue Xiao, Bao Yuan, Weiyi Hu, Jiajia Qi, Hao Jiang, Boxing Sun, Jiabao Zhang, Shuang Liang Tributyltin Oxide Exposure During in vitro Maturation Disrupts Oocyte Maturation and Subsequent Embryonic Developmental Competence in Pigs *Front Cell Dev Biol* 2021 Jun 28;9:683448.
90. Qi Xu, Demin Cheng, Guanru Li, Yi Liu, Ping Li, Wenqing Sun, Dongyu Ma, Chunhui Ni CircHIPK3 regulates pulmonary fibrosis by facilitating glycolysis in miR-30a-3p/FOXK2-dependent manner *Int J Biol Sci* 2021 Jun 4;17(9):2294-2307.
91. Shang Wu, Yang Yang, Tingwei Wang, Jiadi Sun, Yinzhi Zhang, Jian Ji, Xiulan Sun Effects of acid, alkaline, cold, and heat environmental stresses on the antibiotic resistance of the *Salmonella enterica* serovar Typhimurium *Food Res Int* 2021 Jun;144:110359.
92. Yunzhen Xie, Jiaping Zheng, Shiqi Li, Huiying Li, Yu Zhou, Wenrong Zheng, Meilian Zhang, Libin Liu, Zhou Chen GLP-1 improves the neuronal supportive ability of astrocytes in Alzheimer's disease by regulating mitochondrial dysfunction via the cAMP/PKA pathway *Biochem Pharmacol* 2021 Jun;188:114578.
93. Ye Wu, Dong Lu, Yixin Jiang, Jinmei Jin, Sanhong Liu, Lili Chen, Hong Zhang, Yudong Zhou, Hongzhan Chen, Dale G Nagle, Xin Luan, Weidong Zhang Stapled Wasp Venom-Derived Oncolytic Peptides with Side Chains Induce Rapid Membrane Lysis and Prolonged Immune Responses in Melanoma *J Med Chem* 2021 May 13;64(9):5802-5815.
94. Weiyi Hu, Yan Zhang, Dali Wang, Tingting Yang, Jiajia Qi, Yonghong Zhang, Hao Jiang, Jiabao Zhang, Boxing Sun, Shuang Liang Iron Overload-Induced Ferroptosis Impairs Porcine Oocyte Maturation and Subsequent Embryonic Developmental Competence in vitro *Front Cell Dev Biol* 2021 May 28;9:673291.
95. Yan Huang, Gang Wang, Zhan Zhou, Zhengshan Tang, Ningning Zhang, Xiaoyan Zhu, Xin Ni Endogenous Hydrogen Sulfide Is an Important Factor in Maintaining Arterial Oxygen Saturation *Front Pharmacol* 2021 May 31;12:677110.
96. Yuefei Fang, Yang He, Canhao Wu, Meng Zhang, Zeyun Gu, Jiaxin Zhang, Ergang Liu, Qin Xu, Akmal M Asrorov, Yongzhuo Huang Magnetism-mediated targeting hyperthermia-immunotherapy in "cold" tumor with CSF1R inhibitor *Theranostics* 2021 May 3;11(14):6860-6872.
97. Jiaping Zheng, Yunzhen Xie, Lingjia Ren, Liqin Qi, Li Wu, Xiaodong Pan, Jianxing Zhou, Zhou Chen, Libin Liu GLP-1 improves the supportive ability of astrocytes to neurons by promoting aerobic glycolysis in Alzheimer's disease *Mol Metab* 2021 May;47:101180.
98. Tingpeng Yang, Yanzhi Wang, Weijie Liao, Shikuan Zhang, Songmao Wang, Naihan Xu, Weidong Xie, Cheng Luo, Yangyang Wang, Ziqiang Wang, Yaou Zhang Down-regulation of EPB41L4A-AS1 mediated the brain aging and neurodegenerative diseases via damaging synthesis of NAD⁺ and ATP *Cell Biosci* 2021 Nov 10;11(1):192.
99. Yixian Zhou, Boyi Niu, Yiting Zhao, Jintao Fu, Ting Wen, Kaixin Liao, Guilan Quan, Xin Pan, Chuanbin Wu Multifunctional nanoreactors-integrated microneedles for cascade reaction-enhanced cancer therapy *J Control Release* 2021 Nov 10;339:335-349.
100. Beibei Xie, Tingting Liu, Shuang Chen, Yan Zhang, Dongxian He, Qian Shao, Zhen Zhang, Chenhui Wang Combination of DNA demethylation and chemotherapy to trigger cell pyroptosis for inhalation treatment of lung cancer *Nanoscale* 2021 Nov 18;13(44):18608-18615.
101. Simeng Zhang, Zhongyan Hua, Gen Ba, Ning Xu, Jianing Miao, Guifeng Zhao, Wei Gong, Zhihui Liu, Carol J Thiele, Zhijie Li Antitumor effects of the small molecule DMAMCL in neuroblastoma via suppressing aerobic glycolysis and targeting PFKFB3 *Cancer Cell Int* 2021 Nov 24;21(1):619.
102. Hua Yu, Ying Bai, Jing Qiu, Xiaomei He, Junzhi Xiong, Qian Dai, Xingmin Wang, Yuanyuan Li, Halei Sheng, Rong Xin, Lu Jiang, Qiaoqiao Li, Defeng Li, Hong Zhang, Le Zhang, Qian Chen, Jin Peng, Xiaomei Hu, Kebin Zhang *Pseudomonas aeruginosa* PerV Enhances the Nitric Oxide-Mediated Tumoricidal Activity of Tumor-Associated Macrophages via a TLR4/PI3K/AKT/mTOR-Glycolysis-Nitric Oxide Circuit *Front Oncol* 2021 Nov 25;11:736882.
103. Muhan Jing, Shanshan Wang, Ding Li, Zeyu Wang, Ziwen Li, Yichen Lu, Tong Sun, Chen Qiu, Fang Chen, Haijuan Yu, Wei Zhang Lorcaserin Inhibit Glucose-Stimulated Insulin Secretion and Calcium Influx in Murine Pancreatic Islets *Front Pharmacol* 2021 Nov 5;12:761966.
104. Ye He, Jin Leng, Ke Li, Kun Xu, Chuanchuan Lin, Zhang Yuan, Rui Zhang, Danyang Wang, Bailong Tao, Tony Jun Huang, Kaiyong Cai A multifunctional hydrogel coating to direct fibroblast activation and infected wound healing via simultaneously controllable photobiomodulation and photodynamic therapies

- Biomaterials 2021 Nov;278:121164.
105. An-Te Ou, Jia-Xin Zhang, Yue-Fei Fang, Rong Wang, Xue-Ping Tang, Peng-Fei Zhao, Yu-Ge Zhao, Meng Zhang, Yong-Zhuo Huang Disulfiram-loaded lactoferrin nanoparticles for treating inflammatory diseases *Acta Pharmacol Sin* 2021 Nov;42(11):1913-1920.
 106. Nanshan Song, Yinqian Fang, Hong Zhu, Jiaqi Liu, Siyuan Jiang, Sifan Sun, Rong Xu, Jianhua Ding, Gang Hu, Ming Lu Kir6.2 is essential to maintain neurite features by modulating PM20D1-reduced mitochondrial ATP generation *Redox Biol* 2021 Nov;47:102168.
 107. Liu Wang, Guangzhe Li, Lei Cao, Yi Dong, Yang Wang, Shisheng Wang, Yueqing Li, Xiuhan Guo, Yi Zhang, Fangfang Sun, Xuemei Du, Jiangan Su, Qing Li, Xiaojun Peng, Kun Shao, Weijie Zhao An ultrasound-driven immune-boosting molecular machine for systemic tumor suppression *Sci Adv* 2021 Oct 22;7(43):eabj4796.
 108. Wanbing He, Jieping Huang, Yang Liu, Changming Xie, Kun Zhang, Xinhong Zhu, Jie Chen, Hui Huang Deletion of soluble epoxide hydrolase suppressed chronic kidney disease-related vascular calcification by restoring Sirtuin 3 expression *Cell Death Dis* 2021 Oct 23;12(11):992.
 109. Lin Meng, Hongmei Hu, Zhiqiang Liu, Luyao Zhang, Qingrui Zhuan, Xue Li, Xiangwei Fu, Shien Zhu, Yunpeng Hou The Role of Ca²⁺ in Maturation and Reprogramming of Bovine Oocytes: A System Study of Low-Calcium Model *Front Cell Dev Biol* 2021 Oct 26;9:746237.
 110. Lei Shi, Li-Li Gao, Shi-Zhong Cai, Qian-Wei Xiong, Zhou-Rui Ma A novel selective mitochondrial-targeted curcumin analog with remarkable cytotoxicity in glioma cells *Eur J Med Chem* 2021 Oct 5;221:113528.
 111. Zeyu Zhang, Yarui Ma, Xiaolei Guo, Yingxi Du, Qing Zhu, Xiaobing Wang, Changzhu Duan FDX1 can Impact the Prognosis and Mediate the Metabolism of Lung Adenocarcinoma *Front Pharmacol* 2021 Oct 8;12:749134.
 112. Lina Liang, Yan Peng, Liyan Qiu Mitochondria-targeted vitamin E succinate delivery for reversal of multidrug resistance *J Control Release* 2021 Sep 10;337:117-131.
 113. Mingqiang Ye, Yujing Kong, Cuiling Zhang, Yifei Lv, Shasha Cheng, Dongyan Hou, Yuezhong Xian Near-Infrared Light Controllable DNA Walker Driven by Endogenous Adenosine Triphosphate for in Situ Spatiotemporal Imaging of Intracellular MicroRNA *ACS Nano* 2021 Sep 28;15(9):14253-14262.
 114. Qing Wu, Ming Zhao, Xi He, Runqing Xue, Dongling Li, Xiaojiang Yu, Shengpeng Wang, Weijin Zang Acetylcholine reduces palmitate-induced cardiomyocyte apoptosis by promoting lipid droplet lipolysis and perilipin 5-mediated lipid droplet-mitochondria interaction *Cell Cycle* 2021 Sep;20(18):1890-1906.
 115. Shan Zhou, Xian Ling, Ping Meng, Ye Liang, Kunyu Shen, Qinyu Wu, Yunfang Zhang, Qiyang Chen, Shuangqin Chen, Youhua Liu, Lili Zhou Cannabinoid receptor 2 plays a central role in renal tubular mitochondrial dysfunction and kidney ageing *J Cell Mol Med* 2021 Sep;25(18):8957-8972.
 116. Yan Guo, Suxing Jin, Hao Yuan, Tao Yang, Kun Wang, Zijian Guo, Xiaoyong Wang DNA-Unresponsive Platinum(II) Complex Induces ERS-Mediated Mitophagy in Cancer Cells *J Med Chem* 2022 Jan 13;65(1):520-530.
 117. Shuang Qiu, Zhengshuai Wu, Zhipeng Chen, Abdul-Wahab Abbew, Jinxiang Li, Shijian Ge Microalgal Activity and Nutrient Uptake from Wastewater Enhanced by Nanoscale Zerovalent Iron: Performance and Molecular Mechanism *Environ Sci Technol* 2022 Jan 4;56(1):585-594.
 118. Ke Gong, Qingqi Hong, Huiwen Wu, Fen Wang, Linhai Zhong, Lingjia Shen, Piao Xu, Wenqing Zhang, Hanwei Cao, Yan-Yan Zhan, Tianhui Hu, Xiaoting Hong Gap junctions mediate glucose transfer to promote colon cancer growth in three-dimensional spheroid culture *Cancer Lett* 2022 Apr 10;531:27-38.
 119. Lu Lu, He Zhang, Yudong Zhou, Jiayi Lin, Weidong Gao, Ting Yang, Jinmei Jin, Lijun Zhang, Dale G Nagle, Weidong Zhang, Ye Wu, Hongzhan Chen, Xin Luan Polymer chimera of stapled oncolytic peptide coupled with anti-PD-L1 peptide boosts immunotherapy of colorectal cancer *Theranostics* 2022 Apr 24;12(7):3456-3473.
 120. Jingjing Zhang, Xinli Zhou, Danyang Wan, Li Yu, Xu Chen, Tong Yan, Zhu Wu, Meimei Zheng, Feng Zhu, Hui Zhu TMRSS12 Functions in Meiosis and Spermiogenesis and Is Required for Male Fertility in Mice *Front Cell Dev Biol* 2022 Apr 25;10:757042.
 121. Xiaoli Huang, Hongfeng Guo, Xuejun Cheng, Jinyu Zhang, Wenzheng Qu, Qianyun Ding, Qihang Sun, Qiang Shu, Xuekun Li NAD⁺ Modulates the Proliferation and Differentiation of Adult Neural Stem/Progenitor Cells via Akt Signaling Pathway *Cells* 2022 Apr 9;11(8):1283.
 122. Weidong Li, Jiaming Zhang, Xiaoxia Yu, Fei Meng, Ju Huang, Liangran Zhang, Shunxin Wang Aristolochic acid I exposure decreases oocyte quality *Front Cell Dev Biol* 2022 Aug 11;10:838992.
 123. Bo Wang, Feng-Wan Zhang, Wei-Xiao Wang, Yan-Yan Zhao, Su-Yue Sun, Jin-Hong Yu, Michael P Vitek, George F Li, Rui Ma, Shiwei Wang, Zhiliang Hu, Wei Chen Apolipoprotein E mimetic peptide COG1410 combats pandrug-resistant *Acinetobacter baumannii* *Front Microbiol* 2022 Aug 23;13:934765.
 124. Qing He, Peng Hao, Gang He, Hantao Mai, Wenzhou Liu, Weiqiong Zhang, Kelin Zhang, Guifang Zhong, Ruilian Guo, Changzhi Yu, Yang Li, Chipiu Wong, Qian Chen, Yantao Chen IGF2BP1-regulated expression of ERRA is involved in metabolic reprogramming of chemotherapy resistant osteosarcoma cells *J Transl Med* 2022 Aug 2;20(1):348.
 125. Lei Bai, Wen-Yo Tu, Yatao Xiao, Kejing Zhang, Chengyong Shen Motoneurons innervation determines the distinct gene expressions in multinucleated myofibers *Cell Biosci* 2022 Aug 30;12(1):140.
 126. Ping Jin, Jingwen Jiang, Li Zhou, Zhao Huang, Siyuan Qin, Hai-Ning Chen, Liyuan Peng, Zhe Zhang, Bowen Li, Maochao Luo, Tingting Zhang, Hui Ming, Ning Ding, Lei Li, Na Xie, Wei Gao, Wei Zhang, Edouard C Nice, Yuquan Wei, Canhua Huang Disrupting metformin adaptation of liver cancer cells by targeting the TOMM34/ATP5B axis *EMBO Mol Med* 2022 Dec 7;14(12):e16082.
 127. Rui Chen, Tingbin Ma, Shiyue Du, Junyu Luo, Huan Zhang, Xuan Xu, Zhijian Cao, Zhangqi Yuan, Hao Sun, Mugen Liu, Bo Xiong, Qinghua Shi, Jing Yu Liu Impaired fertility in 493059J08Rik mutant male mice is associated with defective sperm energy metabolism *FASEB J* 2022 Dec;36(12):e22634.
 128. Yi Li, Yu-Chen Li, Xiao-Tian Liu, Lu Zhang, Yi-Hua Chen, Qiong Zhao, Wen Gao, Baolin Liu, Hua Yang, Ping Li Blockage of citrate export prevents TCA cycle fragmentation via Irg1 inactivation *Cell Rep* 2022 Feb 15;38(7):110391.
 129. Xuan Wang, Yuting Ji, Dekui Jin, Jingyi Qi, Xuening Hou, Wenting Zhao, Shuaishuai Zhou, Chengying Zhang, Yongting Luo, Peng An, Junjie Luo Natural Polysaccharide β -Glucan Protects against Doxorubicin-Induced Cardiotoxicity by Suppressing Oxidative Stress *Nutrients* 2022 Feb 21;14(4):906.
 130. Mi Gao, Chunxiao Liang, Wei Hong, Xiaoyuan Yu, Yumin Zhou, Ruiting Sun, Haiqing Li, Haichao Huang, Xuhong Gan, Ze Yuan, Jiahuan Zhang, Juan Chen, Qiudi Mo, Luyao Wang, Biting Lin, Bing Li, Pixian Ran Biomass-related PM2.5 induces mitochondrial fragmentation and dysfunction in human airway epithelial cells *Environ Pollut* 2022 Jan 1;292(Pt B):118464.
 131. Ruikun Wang, Jingru Chen, Wei Wang, Zhuoqian Zhao, Haoran Wang, Shiyu Liu, Fan Li, Yajuan Wan, Jie Yin, Rui Wang, Yuanke Li, Cuizhu Zhang, Hongkai Zhang, Youjia Cao CD40L-armed oncolytic herpes simplex virus suppresses pancreatic ductal adenocarcinoma by facilitating the tumor microenvironment favorable to cytotoxic T cell response in the syngeneic mouse model *J Immunother Cancer* 2022 Jan;10(1):e003809.
 132. Shun Yao, Qingru Zhou, Mingzhu Yang, Ya Li, Xiuxiu Jin, Qingge Guo, Lin Yang, Fangyuan Qin, Bo Lei Multi-mtDNA Variants May Be a Factor Contributing to Mitochondrial Function Variety in the Skin-Derived Fibroblasts of Leber's Hereditary Optic Neuropathy Patients *Front Mol Neurosci* 2022 Jul 13:15:920221.
 133. Hui Luo, Xiaohui Wang, Shuai Song, Yunhan Wang, Qinfu Dan, Hong Ge Targeting stearyl-coa desaturase enhances radiation induced ferroptosis and immunogenic cell death in esophageal squamous cell carcinoma *Oncoimmunology* 2022 Jul 15;11(1):2101769.
 134. Qian Wang, Zhenzhen Sun, Shihan Cao, Xiuli Lin, Mengying Wu, Yuanyuan Li, Jie

- Yin, Wei Zhou, Songming Huang, Aihua Zhang, Yue Zhang, Weiwei Xia, Zhanjun Jia Reduced Immunity Regulator MAVS Contributes to Non-Hypertrophic Cardiac Dysfunction by Disturbing Energy Metabolism and Mitochondrial Homeostasis *Front Immunol* 2022 Jul 1;13:919038.
135. Zhipeng Chen, Shuang Qiu, Mengting Li, Di Zhou, Shijian Ge Instant Inhibition and Subsequent Self-Adaptation of *Chlorella* sp. Toward Free Ammonia Shock in Wastewater: Physiological and Genetic Responses *Environ Sci Technol* 2022 Jul 5;56(13):9641-9650.
136. Kun Zhang, Tao Zhang, Yuhan Yang, Wenling Tu, Hongbin Huang, Yujun Wang, Yuzhuo Chen, Kejian Pan, Zhuojia Chen N6-methyladenosine-mediated LDHA induction potentiates chemoresistance of colorectal cancer cells through metabolic reprogramming *Theranostics* 2022 Jun 13;12(10):4802-4817.
137. Shun Deng, Enming Zhang, Yan Wang, Yunyang Zhao, Zezhong Yang, Bingxin Zheng, Xiaoyuan Mu, Xuangen Deng, Hai Shen, Haibo Rong, Desheng Pei In vivo toxicity assessment of four types of graphene quantum dots (GQDs) using mRNA sequencing *Toxicol Lett* 2022 Jun 15:363:55-66.
138. Zheng-Hao Zhang, Zi-Guan Zhang, Min-Wei Chen, Ying Yang, Run-Jing Li, Jia-Jia Xu, Cui Yang, Yu-Ying Li, Hong-Wei Chen, Shi-Xiao Liu, Yan-Ling Li, Ping Luo, Yi-Jiang Liu, Wen-Bo Chen, Zhong-Gui Shan, Zheng-Rong Huang Inhibition of GSDMD Activates Poly(ADP-ribosylation) and Promotes Myocardial Ischemia-Reperfusion Injury *Oxid Med Cell Longev* 2022 Jun 24:2022:1115749.
139. Yajuan Lu, Hanyu Tang, Junjie Xu, Fei Sun Toxic effects of 4-methylimidazole on the maturation and fertilization of mouse oocytes *Food Chem Toxicol* 2022 Jun:164:113051.
140. Debin Zheng, Jingfei Liu, Limin Xie, Yuhan Wang, Yinghao Ding, Rong Peng, Min Cui, Ling Wang, Yongjie Zhang, Chunqiu Zhang, Zhimou Yang Enzyme-instructed and mitochondria-targeting peptide self-assembly to efficiently induce immunogenic cell death *Acta Pharm Sin B* 2022 Jun;12(6):2740-2750.
141. Yao Liu, Yiping Lu, Bo Ning, Xiaomin Su, Binru Yang, Haiqing Dong, Bo Yin, Zhiqing Pang, Shun Shen Intravenous Delivery of Living *Listeria monocytogenes* Elicits Gasdmermin-Dependent Tumor Pyroptosis and Motivates Anti-Tumor Immune Response *ACS Nano* 2022 Mar 22;16(3):4102-4115.
142. Lianhua Sun, Dekun Gao, Junmin Chen, Shule Hou, Yue Li, Yuyu Huang, Fabio Mammano, Jianyong Chen, Jun Yang Failure Of Hearing Acquisition in Mice With Reduced Expression of Connexin 26 Correlates With the Abnormal Phasing of Apoptosis Relative to Autophagy and Defective ATP-Dependent Ca²⁺ Signaling in Köll *Front Cell Neurosci* 2022 Mar 3:16:816079.
143. Hai-Fan Yu, Zhan-Qing Yang, Ming-Yue Xu, Ji-Cheng Huang, Zhan-Peng Yue, Bin Guo Yap is essential for uterine decidualization through Rrm2/GSH/ROS pathway in response to Bmp2 *Int J Biol Sci* 2022 Mar 6;18(6):2261-2276.
144. Luyao Zhang, Kexiong Liu, Qingrui Zhuan, Zhiqiang Liu, Lin Meng, Xiangwei Fu, Gongxue Jia, Yunpeng Hou Mitochondrial Calcium Disorder Affects Early Embryonic Development in Mice through Regulating the ERK/MAPK Pathway *Oxid Med Cell Longev* 2022 May 20:2022:8221361.
145. Huili Yang, Yadong Xue, Bing Li, Yanan Lin, Haochuan Li, Zhanyong Guo, Weihua Li, Zhiyuan Fu, Dong Ding, Jihua Tang The chimeric gene *atp6c* confers cytoplasmic male sterility in maize by impairing the assembly of the mitochondrial ATP synthase complex *Mol Plant* 2022 May 2;15(5):872-886.
146. Zhumei Gao, Chuyue Zhang, Fei Peng, Qianqian Chen, Yinghua Zhao, Liangmei Chen, Xu Wang, Xiangmei Chen Hypoxic mesenchymal stem cell-derived extracellular vesicles ameliorate renal fibrosis after ischemia-reperfusion injury by restoring CPT1A mediated fatty acid oxidation *Stem Cell Res Ther* 2022 May 7;13(1):191.
147. Yajuan Lu, Hanyu Tang, Xia Wang, Junjie Xu, Fei Sun Dibutyltin dichloride exposure affects mouse oocyte quality by inducing spindle defects and mitochondria dysfunction *Chemosphere* 2022 May:295:133959.
148. Boyu Yuan, Mingming Liu, Yuhong Gong, Zifan Wang, Xinxin Jin, Gaijie Xie, Mingqiang Zhu, Xue Zhang, Siyuan Luo, Qing Qu, Yufeng Zhu, Meng Wang, Yingli Jin, Bai Li, Wei Wang Sodium butyrate exerts antioxidant stress effects and attenuates A β 25-35-induced cytotoxicity in PC12 cells *Arch Biochem Biophys* 2022 Nov 30:731:109448.
149. Yang Liu, Lei Cao, Yuting Song, Zhengwei Kang, Ting Liu, Jianhua Ding, Gang Hu, Ming Lu Mitochondrial glutamine transporter SLC1A5_var, a potential target to suppress astrocyte reactivity in Parkinson's Disease *Cell Death Dis* 2022 Nov 9;13(11):946.
150. Yi Wu, Minhua Zong, Zhenhui Zhang, Yongtai Wu, Lin Li, Xia Zhang, Hong Wu, Bing Li Selective transportation and energy homeostasis regulation of dietary advanced glycation end-products in human intestinal Caco-2 cells *Food Chem* 2022 Oct 15:391:133284.
151. Jin-Li Ding, Xiu-Hui Li, Jia-Hui Lei, Ming-Guang Feng, Sheng-Hua Ying Succinate Dehydrogenase Subunit C Contributes to Mycelial Growth and Development, Stress Response, and Virulence in the Insect Parasitic Fungus *Beauveria bassiana* *Microbiol Spectr* 2022 Oct 26;10(5):e0289122.
152. Wanjun Wang, Yan Liu, Guiying Li, Zhenni Liu, Po Keung Wong, Taicheng An Mechanism insights into bacterial sporulation at natural sphalerite interface with and without light irradiation: The suppressing role in bacterial sporulation by photocatalysis *Environ Int* 2022 Oct:168:107460.
153. Zhijie Zhu, Xuankang Wang, Zhiwen Song, Xiaoshuang Zuo, Yangguang Ma, Zhihao Zhang, Cheng Ju, Zhuowen Liang, Kun Li, Xueyu Hu, Zhe Wang Photobiomodulation promotes repair following spinal cord injury by restoring neuronal mitochondrial bioenergetics via AMPK/PGC-1 α /TFAM pathway *Front Pharmacol* 2022 Sep 12:13:991421.
154. Yuting Ji, Dekui Jin, Jingyi Qi, Xuan Wang, Chengying Zhang, Peng An, Yongting Luo, Junjie Luo Fucoidan Protects against Doxorubicin-Induced Cardiotoxicity by Reducing Oxidative Stress and Preventing Mitochondrial Function Injury *Int J Mol Sci* 2022 Sep 14;23(18):10685.
155. Shuya Gao, Qingchen Yang, Yue Peng, Weixian Kong, Zekun Liu, Zhe Li, Jiawen Chen, Mengmeng Bao, Xie Li, Yubin Zhang, Xiaohong Bian, Liang Jin, Hanwen Zhang, Yuexin Zhang, Daniel Sanchis, Fangrong Yan, Junmei Ye SIRT6 regulates obesity-induced oxidative stress via ENDOG/SOD2 signaling in the heart *Cell Biol Toxicol* 2023 Aug;39(4):1489-1507.
156. Jingwen Dong, Chenfei Zhu, Ying Huang, Quanhao Li, Jing Li, Zheng Wang, Yixin Wang, Zhanwei Zhou, Minjie Sun Reversing the PAI-1-induced fibrotic immune exclusion of solid tumor by multivalent CXCR4 antagonist nano-permeator *Acta Pharm Sin B* 2023 Jul;13(7):3106-3120.

Version 2024.03.12