NADP⁺/NADPH检测试剂盒(WST-8法)

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$$[\text{NADP}^+] = [\text{NADP}_{\text{total}}] - [\text{NADPH}]$$

$$[\text{NADP}^+]/[\text{NADPH}] = ([\text{NADP}_{\text{total}}] - [\text{NADPH}])/[\text{NADPH}]$$

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- Ma Y, Qi Y, Wang L, Zheng Z, Zhang Y, Zheng J. SIRT5-mediated SDHA desuccinylation promotes clear cell renal cell carcinoma tumorigenesis. *FREE RADICAL BIO MED.* 134:458-467.
- Xie GH, Dai HJ, Liu F, Zhang YP, Zhu L, Nie JJ, Wu JH. A Dual Role of ATM in Ischemic Preconditioning and Ischemic Injury. *Cell Mol Neurobiol.*
- Fenglan Qin, Huige Zhou, Jiayang Li, Jing Liu, Yaling Wang, Ru Bai, Shihui Liu, Manman Ma, Tao Liu, Fene Gao, Peiyao Du, Xiaoquan Lu, Chunying Chen. Hypoxia and pH co-triggered oxidative stress amplifier for tumor therapy. *Eur J Pharmacol.* doi: 10.1016/j.ejphar.2021.174187.
- Jialei Sun, Chenhao Zhou, Yue Zhao, Xiaofei Zhang, Wanyong Chen, Qiang Zhou, Bo Hu, Dongmei Gao, Lisa Raatz, Zhefang Wang, Peter J Nelson, Yuchao Jiang, Ning Ren, Christiane J Bruns, Haijun Zhou. Quiescin sulfhydryl oxidase 1 promotes sorafenib-induced ferroptosis in hepatocellular carcinoma by driving EGFR endosomal trafficking and inhibiting NRF2 activation. *Redox Biol.* doi: 10.1016/j.redox.2021.101942.
- Yi Wen, Hansen Chen, Lu Zhang, Meiling Wu, Feng Zhang, Dan Yang, Jiangang Shen, Jianping Chen. Glycyrrhetic acid induces oxidative/nitrative stress and drives ferroptosis through activating NADPH oxidases and iNOS, and depriving glutathione in triple-negative breast cancer cells. *Free Radic Biol Med.* doi: 10.1016/j.freeradbiomed.2021.07.019.
- Ting Liu, Qianqian Ma, Wenjie Li, Yan Hu, Jun Yang, Qi Yao. Ubiquitin 1 suppresses the cancer stem cell-like traits of non-small cell lung cancer cells by regulating reactive oxygen species homeostasis. *Bioengineered.* doi: 10.1080/21655979.2021.1979353.
- Yan Wang, Haojie Jin, Yafang Wang, Yang Yao, Cuixia Yang, Jihong Meng, Xiaomu Tan, Yu Nie, Lixiang Xue, Baohui Xu, Heng Zhao, Feng Wang. Sult2b1 deficiency exacerbates ischemic stroke by promoting pro-inflammatory macrophage polarization in mice. *Theranostics.* doi: 10.7150/thno.61646.
- Jian Luo, Qiang He, Jin-Zhi Xu, Chen Xu, Yin-Ze Han, Hai-Long Gao, Xian-Zhi Meng, Guo-Qing Pan, Tian Li, Ze-Yang Zhou. Microsporidia infection upregulates host energy metabolism but maintains ATP homeostasis. *J Invertebr Pathol.* doi: 10.1016/j.jip.2021.107596.
- Mingming Sun, Hao Sheng, Tingfeng Wu, Jiaqi Song, Huanran Sun, Yingzhi Wang, Jiyan Wang, Zhen Li, Huifang Zhao, Junzhen Tan, Yanping Li, Guo Chen, Qingrong Huang, Yuan Zhang, Bei Lan, Shuangping Liu, Changliang Shan, Shuai Zhang. PIKE-A promotes glioblastoma growth by driving PPP flux through increasing G6PD expression mediated by phosphorylation of STAT3. *Biochem Pharmacol.* doi: 10.1016/j.bcp.2021.114736.
- Junjie Xu, Lin Ji, Yeling Ruan, Zhe Wan, Zhongjie Lin, Shunjie Xia, Liye Tao, Junhao Zheng, Liuxin Cai, Yifan Wang, Xiao Liang, Xiujun Cai. UBQLN1 mediates sorafenib resistance through regulating mitochondrial biogenesis and ROS homeostasis by targeting PGC1β in hepatocellular carcinoma. *Signal Transduct Target Ther.* doi: 10.1038/s41392-021-00594-4.
- Wenqi Liang, Chunli Zhao, Zhongrui Chen, Zijing Yang, Ke Liu, Shusheng Gong. Sirtuin-3 Protects Cochlear Hair Cells Against Noise-Induced Damage via the Superoxide Dismutase 2/Reactive Oxygen Species Signaling Pathway. *Front Cell Dev Biol.* doi: 10.3389/fcell.2021.766512.
- Dan He, Qiangdongzi Mao, Jialin Jia, Zhiyu Wang, Yu Liu, Tingting Liu, Bangwei Luo, Zhiren Zhang. Pentose Phosphate Pathway Regulates Tolerogenic Apoptotic Cell Clearance and Immune Tolerance. *Front Immunol.* doi: 10.3389/fimmu.2021.797091.
- Yan Fang, Rongzhang Dou, Sihao Huang, Lei Han, Hang Fu, Chaogang Yang, Jialin Song, Jinsen Zheng, Xinyao Zhang, Keshu Liu, Zhenxian Xiang, Xinghua Zhang, Shuyi Wang, Bin Xiong. LAMC1-mediated preadipocytes differentiation promoted peritoneum pre-metastatic niche formation and gastric cancer metastasis. *Int J Biol Sci.* doi: 10.7150/ijbs.70524.
- Xiaoli Huang, Hongfeng Guo, Xuejun Cheng, Jinyu Zhang, Wenzheng Qu, Qianyun Ding, Qihang Sun, Qiang Shu, Xuekun Li. NAD⁺ Modulates the Proliferation and Differentiation of Adult Neural Stem/Progenitor Cells via Akt Signaling Pathway. *Cells.* doi: 10.3390/cells11081283.
- Feifei Pu, Jianxiang Liu, Doudou Jing, Fengxia Chen, Xin Huang, Deyao Shi, Wei Wu, Hui Lin, Lei Zhao, Zhenhao Zhang, Xiao Lv, Baichuan Wang, Zhicai Zhang, Zengwu Shao. LncCCAT1 interaction protein PKM2 upregulates SREBP2 phosphorylation to promote osteosarcoma tumorigenesis by enhancing the Warburg effect and lipogenesis. *Int J Oncol.* doi: 10.3892/ijo.2022.5334.
- Yixing Gao, Lan Feng, Luping Zhang, Jianhui Geng, Erlong Zhang. ATF4/MYC Regulates MTHFD2 to Promote NSCLC Progression by Mediating Redox Homeostasis. *Dis Markers.* doi: 10.1155/2022/7527996.
- Jia Zhou, Shuangming Yue, Jingjing Du, Benchu Xue, Lizhi Wang, Quanhui Peng, Huawei Zou, Rui Hu, Yahui Jiang, Zhisheng Wang, Bai Xue. Integration of transcriptomic and metabolomic analysis of the mechanism of dietary N-carbamoylglutamate in promoting follicle development in yaks. *Front Vet Sci.* doi: 10.3389/fvets.2022.946893.
- Sijun Liu, Jun Tao, Fengqi Duan, Huangjing Li, Hongmei Tan. HHcy Induces Pyroptosis and Atherosclerosis via the Lipid Raft-Mediated NOX-ROS-NLRP3 Inflammasome Pathway in apoE^{-/-} Mice. *Cells.* doi: 10.3390/cells11152438.
- Wenqiao Ding, Xueying Zhang, Xiaoyu Yin, Qing Zhang, Ying Wang, Changhong Guo, Ying Chen. Ganoderma lucidum aqueous extract inducing PHGPx to inhibit membrane lipid hydroperoxides and regulate oxidative stress based on single-cell animal transcriptome. *Sci Rep.* doi: 10.1038/s41598-022-06985-z.
- Xueyao Zhang, Yingxian Sun. Chromodomain Helicase DNA Binding Protein 1-like, a negative regulator of Forkhead box O3a, promotes the proliferation and migration of Angiotensin II-induced vascular smooth muscle cells. *Bioengineered.* doi: 10.1080/21655979.2021.2019869.
- Zhuan Feng, Yifei Qin, Fei Huo, Zhe Jian, Xiaomin Li, Jiejie Geng, Yong Li, Jiao Wu. NMN recruits GSH to enhance GPX4-mediated

- ferroptosis defense in UV irradiation induced skin injury. *Biochim Biophys Acta Mol Basis Dis*. doi: 10.1016/j.bbadis.2021.166287.
22. Na Song, Jing Wang, Zonglang Lai, Shuting Liang, Wenjuan Zou, Juan Wang, Dandan Zheng, Ying Li, Yuxi He, Jun Cheng, Yue Wu. Arisaema heterophyllum Blume Monomer Stigmasterol Targets PPAR γ and Inhibits the Viability and Tumorigenicity of Lung Adenocarcinoma Cells NCI-H1975. *Evid Based Complement Alternat Med*. doi: 10.1155/2022/5377690.
 23. Siyuan Chen, Bo Ning, Jinwen Song, Zihan Yang, Li Zhou, Zhiji Chen, Linhong Mao, Hongtao Liu, Qingliang Wang, Song He, Zhihang Zhou. Enhanced pentose phosphate pathway activity promotes pancreatic ductal adenocarcinoma progression via activating YAP/MMP1 axis under chronic acidosis. *Int J Biol Sci*. doi: 10.7150/ijbs.69526.
 24. Jing Tan, Ming Yang, Haiping Wang, Conghui Shen, Maoxiong Wu, He Xu, Yandi Wu, Yuanlong Li, Xinghui Li, Tongsheng Huang, Shijie Deng, Zhenyu Yang, Saifei Gao, Hui Li, Jianguo Zhou, Hui Chen, Nan Cao, Weibin Cai. Moderate heart rate reduction promotes cardiac regeneration through stimulation of the metabolic pattern switch. *Cell Rep*. doi: 10.1016/j.celrep.2022.110468.
 25. Xiaochun Yin, Jin Peng, Lihong Gu, Yan Liu, Xihan Li, Jinhui Wu, Bing Xu, Yuzheng Zhuge, Feng Zhang. Targeting glutamine metabolism in hepatic stellate cells alleviates liver fibrosis. *Cell Death Dis*. doi: 10.1038/s41419-022-05409-0.
 26. Yuening Sun, Wenjing Guo, Yongjian Guo, Zhangxing Lin, Dechao Wang, Qinglong Guo, Yuxin Zhou. Apoptosis induction in human prostate cancer cells related to the fatty acid metabolism by wogonin-mediated regulation of the AKT-SREBP1-FASN signaling network. *Food Chem Toxicol*. doi: 10.1016/j.fct.2022.113450.
 27. Peina Wang, Yanmei Cui, Yuanyuan Liu, Zhongda Li, Huiyuan Bai, Yashuo Zhao, Yan-Zhong Chang. Mitochondrial ferritin alleviates apoptosis by enhancing mitochondrial bioenergetics and stimulating glucose metabolism in cerebral ischemia reperfusion. *Redox Biol*. doi: 10.1016/j.redox.2022.102475.
 28. Hong Ding, Chuanfeng Tang, Wei Wang, Ying Pan, Ruiqing Jiao, Lingdong Kong. Polydatin Ameliorates High Fructose-Induced Podocyte Oxidative Stress via Suppressing HIF-1 α /NOX4 Pathway. *Pharmaceutics*. doi: 10.3390/pharmaceutics14102202.
 29. Weifan Lin, Xiangwan Lu, Hang Yang, Linxuan Huang, Wuheng Huang, Yuluan Tang, Situn Liu, Hua Wang, Yan Zhang. Metabolic heterogeneity protects metastatic mucosal melanomas cells from ferroptosis. *Int J Mol Med*. doi: 10.3892/ijmm.2022.5180.
 30. Ben Li, Chufan Wang, Peng Lu, Yumeng Ji, Xufeng Wang, Chaoyang Liu, Xiaohu Lu, Xiaohan Xu, Xiaowei Wang. IDH1 Promotes Foam Cell Formation by Aggravating Macrophage Ferroptosis. *Biology (Basel)*. doi: 10.3390/biology11101392.
 31. Zhen Zhang, Zichen Yu, Jinduo Wang, Yifa Yu, Lanxiao Li, Pengjie Sun, Xiaoguang Fan, Qingyang Xu. Metabolic engineering of *Escherichia coli* for efficient production of L-5-hydroxytryptophan from glucose. *Microb Cell Fact*. doi: 10.1186/s12934-022-01920-3.
 32. Wanzhen Li, Hongxin Jiang, Chen Bai, Shuna Yu, Yitong Pan, Chenchen Wang, Huiting Li, Ming Li, Yaxin Sheng, Fangfang Chu, Jie Wang, Yuting Chen, Jianguo Li, Jiyang Jiang. Ac2-26 attenuates hepatic ischemia-reperfusion injury in mice via regulating IL-22/IL-22R1/STAT3 signaling. *PeerJ*. doi: 10.7717/peerj.14086.
 33. Hailong Shi, Yuanyuan Cheng, Qimei Shi, Wenzhi Liu, Xue Yang, Shuang Wang, Lin Wei, Xiangming Chen, Hao Fang. Myoferlin disturbs redox equilibrium to accelerate gastric cancer migration. *Front Oncol*. doi: 10.3389/fonc.2022.905230.
 34. Yili Hu, Xiaocui Chai, Jun Men, Shen Rao, Xin Cong, Shuiyuan Cheng, Zhixian Qiao. Does Methionine Status Influence the Outcome of Selenomethionine Supplementation? A Comparative Study of Metabolic and Selenium Levels in HepG2 Cells. *Nutrients*. doi: 10.3390/nu14183705.
 35. Jingru Bai, Tong Liu, Bo Tu, Meng Yuan, Zhaoqi Shu, Minghe Fan, Sihan Huo, Yuyao Guo, Lina Wang, Hua Wang, Ying Zhao. Autophagy loss impedes cancer-associated fibroblast activation via downregulating proline biosynthesis. *Autophagy*. doi: 10.1080/15548627.2022.2093026.
 36. Luo-Yi Wang, Rui-Xue Wang, Chen Wang, Su-Fang Chen, Xiao-Jing Sun, Zhi-Ying Li, Min Chen, Mark A Little, Ming-Hui Zhao. Inhibitor of apoptosis proteins antagonist SM164 ameliorates experimental MPO-ANCA-associated vasculitis via enhancing fatty acid oxidation in neutrophils. *Rheumatology (Oxford)*. doi: 10.1093/rheumatology/keac621.

Version 2024.04.18