



碧云天生物技术/Beyotime Biotechnology  
订货热线: 400-1683301或800-8283301  
订货e-mail: order@beyotime.com  
技术咨询: info@beyotime.com  
网址: http://www.beyotime.com

## 碱性磷酸酯酶标记山羊抗兔IgG (H+L)

| 产品编号  | 产品名称                  | 包装    |
|-------|-----------------------|-------|
| A0239 | 碱性磷酸酯酶标记山羊抗兔IgG (H+L) | 0.5ml |

### 产品简介:

| 来源   | 用途             | WB         | ELISA      | IHC        |
|------|----------------|------------|------------|------------|
| Goat | WB, ELISA, IHC | 1:500-5000 | 1:500-5000 | 1:500-2000 |

WB, Western blot; IHC, Immunohistochemistry.

- 本碱性磷酸酯酶标记山羊抗兔IgG (H+L) (AP-labeled Goat Anti-Rabbit IgG (H+L))为进口分装, 用于Western、ELISA和IHC等的高纯度高特异性二抗。
- 碱性磷酸酯酶(Alkaline Phosphatase, AP/ALP/AKP/ALKP/ALPase/Alk Phos)常被称作碱性磷酸酶(EC 3.1.3.1), 是一类水解酶, 通过水解磷酸单酯将底物分子上的磷酸基团除去, 并生成磷酸根离子和自由的羟基, 其去磷酸化作用的底物包括蛋白质、核苷酸和生物碱等, 并在碱性条件下最为有效。该酶是一组同功酶的统称。常见的小牛肠碱性磷酸酶(Calf Intestinal Alkaline Phosphatase, CIAP/CIP)被广泛用于二抗等的标记最终用于蛋白和核酸等的检测, 例如在Western、ELISA、IHC、ICC、Northern、Southern等时, 催化BCIP/NBT等显色试剂显色或催化适当的化学发光试剂产生化学发光。本产品标记的就是小牛肠碱性磷酸酶。
- 本抗体为用纯化的兔IgG免疫山羊, 然后用亲和纯化柱对获得的抗血清进行纯化, 并经过人、小鼠和大鼠IgG吸附纯化的优质二抗。对人、小鼠和大鼠IgG几乎没有结合能力(Minimal Cross-reaction to Human, Mouse and Rat Serum Proteins)。
- 本碱性磷酸酯酶标记山羊抗兔IgG (H+L)用于WB、ELISA和IHC的推荐稀释比例参考上表, 实际实验操作过程中需根据抗原和抗体的具体情况适当调节二抗的稀释比例。
- 本抗体如果用于常规的Western检测, 以每次检测需10毫升1:1000稀释的二抗计, 可以检测50次。

### 包装清单:

| 产品编号  | 产品名称                  | 包装    |
|-------|-----------------------|-------|
| A0239 | 碱性磷酸酯酶标记山羊抗兔IgG (H+L) | 0.5ml |
| -     | 说明书                   | 1份    |

### 保存条件:

-20°C保存, 一年有效。

### 注意事项:

- 本产品仅限于专业人员的科学的研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

### 使用说明:

1. Western、ELISA或IHC请参考相关实验步骤进行。起始稀释浓度按照产品简介中推荐的稀释比例进行稀释。

### 使用本产品的文献:

1. Yan C, Ji Y, Huang T, Yu F, Gao Y, Gu Y, Qi Q, Du J, Dai J, Ma H, Jin G . An esophageal adenocarcinoma susceptibility locus at 9q22 also confers risk to esophageal squamous cell carcinoma by regulating the function of BARX1. Cancer Lett. 2018 May 1;421:103-111.
2. Nannan Gao, Feng Yang, Siyuan Chen, Hanxing Wan, Xiaoyan Zhao, Hui Dong . The role of TRPV1 ion channels in the suppression of gastric cancer development J EXP CLIN CANC RES. 2020 Oct 2;39(1):206.
3. 姜 艳, 金 英, 王世兴, 李亚男, 闫恩志, 齐志敏, 魏 佳. 吡格列酮对淀粉样β蛋白片段 1-42 引起的大鼠海马神经细胞凋亡相关蛋白表达的影响. 中国药理学与毒理学杂志. 2008 年 10 月 22(5):355-61.
4. Hu L, Wang Z, Du H, Huang B. . Differential accumulation of dehydrins in response to water stress for hybrid and common bermudagrass genotypes differing in drought tolerance. J Plant Physiol. 2010 Jan 15;167(2):103-9.
5. Yu Q, Guo W, Song X, Liu X, Xiang Z, He C, Burnstock G. . Expression of P2Y receptors in the rat anterior pituitary. PURINERG SIGNAL. 2011 Jun;7(2):207-19.
6. Tang DQ, Wei YQ, Yin XX, Lu Q, Hao HH, Zhai YP, Wang JY, Ren J. . In vitro suppression of quercetin on hypertrophy and extracellular matrix accumulation in rat glomerular mesangial cells cultured by high glucose. Fitoterapia. 2011 Sep;82(6):920-6.

7. Yu J, Feng L, Hu Y, Zhou Y . Effects of SAC on oxidative stress and NO availability in placenta: potential benefits to preeclampsia. *Placenta*. 2012 Jun;33(6):487-94.
8. Liu J, Luo Y, Liu Q, Zheng L, Yang Z, Wang Y, Su J, Quan F, Zhang Y . Production of cloned embryos from caprine mammary epithelial cells expressing recombinant human $\beta$ -defensin-3. *Theriogenology*. 2013 Mar 1;79(4):660-6.
9. Yang G, Cui T, Wang Y, Sun S, Ma T, Wang T, Chen Q, Li Z . Selective isolation and analysis of glycoprotein fractions and their glycomes from hepatocellular carcinoma sera. *PROTEOMICS*. 2013 May;13(9):1481-98.
10. Liu J, Luo Y, Ge H, Han C, Zhang H, Wang Y, Su J, Quan F, Gao M, Zhang Y . Anti-Bacterial Activity of Recombinant Human $\beta$ -Defensin-3 Secreted in the Milk of Transgenic Goats Produced by Somatic Cell Nuclear Transfer. *PLoS One*. 2013 Jun 14;8(6):e65379.
11. Shen XL, Zhang B, Liang R, Cheng WH, Xu W, Luo Y, Zhao C, Huang K . Central role of Nix in the autophagic response to ochratoxin A. *Food Chem Toxicol*. 2014 Jul;69:202-9.
12. He H, Yu WG, Yang JP, Ge S, Lu YH . Multiple Comparisons of Glucokinase Activation Mechanisms of Five Mulberry Bioactive Ingredients in Hepatocyte. *J AGR FOOD CHEM*. 2016 Mar 30;64(12):2475-84.
13. Ye H, Wu Q, Guo M, Wu K, Lv Y, Yu F, Liu Y, Gao X, Zhu Y, Cui L, Liang N, Yun T, Li L, Zheng X . Growth inhibition effects of ent-11 $\alpha$ -hydroxy-15-oxo-kaur-16-en-19-oic-acid on colorectal carcinoma cells and colon carcinoma-bearing mice. *Mol Med Rep*. 2016 Apr;13(4):3525-32.
14. Qian LQ, Li XQ, Ye PH, Su HY, Wang G, Liu Y, Shen GH, Gao QG . Downregulation of MACC1 inhibits the viability, invasion and migration and induces apoptosis in esophageal carcinoma cells through the phosphatase and tensin homolog/phosphoinositide 3-kinase/protein kinase B signaling pathway. *Oncol Lett*. 2017 Oct;14(4):4897-4905.
15. Luo L, Wang Y, Yin Y, Ge J, Lu X . Effects of NOB1 on the pathogenesis of osteosarcoma and its expression on the chemosensitivity to cisplatin. *Oncol Lett*. 2018 Mar;15(3):3548-3551.
16. Du H, Ma L, Chen G, Li S . The effects of oxyresveratrol abrogates inflammation and oxidative stress in rat model of spinal cord injury. *Mol Med Rep*. 2018 Mar;17(3):4067-4073.
17. Lin L, Huang S, Zhu Z, Han J, Wang Z, Huang W, Huang Z . P2X7 receptor regulates EMMPRIN and MMP-9 expression through AMPK/MAPK signaling in PMA-induced macrophages. *Mol Med Rep*. 2018 Sep;18(3):3027-3033.
18. Yang T, Cheng J, You J, Yan B, Liu H, Li F . S100B promotes chemoresistance in ovarian cancer stem cells by regulating p53. *Oncol Rep*. 2018 Sep;40(3):1574-1582.
19. Xue X, Lv Y, Leng Y, Zhang Y . Protective effects of extracellular polymeric substances from Aphanizomenon flos-aquae on neurotoxicity induced by local anesthetics. *Exp Ther Med*. 2018 Oct;16(4):3011-3019.
20. Li Z, Liu H, Zhong Q, Wu J, Tang Z . LncRNA UCA1 is necessary for TGF- $\beta$ -induced epithelial-mesenchymal transition and stemness via acting as a ceRNA for Slug in glioma cells. *FEBS Open Bio*. 2018 Oct 18;8(11):1855-1865.
21. Ma Y, Zhang C, Zhang B, Yu H, Yu Q . circRNA of AR-suppressed PABPC1 91 bp enhances the cytotoxicity of natural killer cells against hepatocellular carcinoma via upregulating UL16 binding protein 1. *Oncol Lett*. 2019 Jan 17(1):388-397.
22. Wen Q, Xu C, Zhou J, Liu NM, Cui YH, Quan MF, Cao JG, Ren KQ . 8-bromo-7-methoxychrysins suppress stemness of SMMC-7721 cells induced by co-culture of liver cancer stem-like cells with hepatic stellate cells. *BMC Cancer*. 2019 Mar 12 19(1):224.
23. Zhang L, Chen Y, Li C, Liu J, Ren H, Li L, Zheng X, Wang H, Han Z . RNA binding protein PUM2 promotes the stemness of breast cancer cells via competitively binding to neuropilin-1 (NRP-1) mRNA with miR-376a. *Biomed Pharmacother*. 2019 Jun 114:108772.
24. Niu G, Yang Y, Ren J, Song T, Hu Z, Chen L, Hong R, Xia J, Ke C, Wang X . Overexpression of CPXM2 predicts an unfavorable prognosis and promotes the proliferation and migration of gastric cancer. *Oncol Rep*. 2019 Oct 42(4):1283-1294.
25. Tian X, He Y, Han Z, Su H, Chu C . The Cytoplasmic Expression Of CLDN12 Predicts An Unfavorable Prognosis And Promotes Proliferation And Migration Of Osteosarcoma. *Cancer Manag Res*. 2019 Nov 1 11:9339-9351.
26. Wang YF, Han J . OTOR in breast carcinoma as a potent prognostic predictor correlates with cell proliferation, migration, and invasiveness. *Biochem Cell Biol*. 2019 Dec 97(6):750-757.
27. Xinyu Liu, Feiran Song, Chunna Liu, Yi Zhang . 25-OH-PPD inhibits hypertrophy on diabetic cardiomyopathy via the PI3k/Akt/GSK-3 $\beta$  signaling pathway *Exp Ther Med*. 2020 Sep;20(3):2141-2147.
28. Feng-Ying Yuan, Ming-Xing Zhang, Yi-Hua Shi, Mei-Hui Li, Jia-Yuan Ou, Wen-Fang Bai, Ming-Sheng Zhang . Bone marrow stromal cells-derived exosomes target DAB2IP to induce microglial cell autophagy, a new strategy for neural stem cell transplantation in brain injury *Exp Ther Med*. 2020 Sep;20(3):2752-2764.
29. Kepeng Ou, Youjian Li, Yiling Long, Yafei Luo, Dianyong Tang, Zhongzhu Chen . Inhibition of MORC2 Mediates HDAC4 to Promote Cellular Senescence through p53/p21 Signaling Axis Molecules. *2022 Sep 22;27(19):6247*.
30. Heng Sun, Wenjie Han, Juyi Wen, Xue Ma . IL4I1 and tryptophan metabolites enhance AHR signals to facilitate colorectal cancer progression and immunosuppression *Am J Transl Res*. 2022 Nov 15;14(11):7758-7770.

Version 2024.03.12