



碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

Hanks' Balanced Salt Solution

产品编号	产品名称	包装
C0218	Hanks' Balanced Salt Solution	500ml

产品简介:

- 碧云天生产的Hanks' Balanced Salt Solution, 简称HBSS, 含137.93mM NaCl, 5.33mM KCl, 4.17mM NaHCO₃, 0.441mM KH₂PO₄, 0.338mM Na₂HPO₄, 5.56mM D-Glucose, pH值为7.1-7.5。
- 本Hanks' Balanced Salt Solution溶液经过过滤除菌, 可以直接用于细胞培养过程中细胞的洗涤等常规用途。使用前不必再进行稀释或过滤除菌等任何处理。

包装清单:

产品编号	产品名称	包装
C0218	Hanks' Balanced Salt Solution	500ml
—	说明书	1份

保存条件:

4°C保存, 半年有效。-20°C保存一年有效。本产品可以短期室温存放。

注意事项:

- 在使用Hanks' Balanced Salt Solution的过程中要特别注意避免溶液被微生物污染。
- 本Hanks' Balanced Salt Solution溶液不含Ca²⁺与Mg²⁺, 如需使用含有Ca²⁺与Mg²⁺的Hanks' Balanced Salt Solution 溶液, 请订购C0219 Hanks' Balanced Salt Solution (with Ca²⁺ & Mg²⁺)。
- 本产品仅限于专业人员的科学研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

使用本产品的文献:

1. Li Y, Su H, Pi X, Gong Y, Xiong X, Wu G, Li Z. . The Stimulatory Activities of Baicalein and Baicalin Compounds Derived from *Scutellaria baicalensis* on insulin secretion in vitro. *ADV MATER.* 2012 Jul;554-556(10):1673-7.
2. Chai Q, Wang X, Qiang L, Zhang Y, Ge P, Lu Z, Zhong Y, Li B, Wang J, Zhang L, Zhou D, Li W, Dong W, Pang Y, Gao GF, Liu CH . A Mycobacterium tuberculosis surface protein recruits ubiquitin to trigger host xenophagy. *Nat Commun.* 2019 Apr 29 10(1):1973.
3. Meng G, Liu J, Liu S, Song Q, Liu L, Xie L, Han Y, Ji Y . Hydrogen sulfide pretreatment improves mitochondrial function in myocardial hypertrophy via a SIRT3-dependent manner. *BRIT J PHARMACOL.* 2018 Apr;175(8):1126-1145.
4. Tang Z, Xie H, Heier C, Huang J, Zheng Q, Eichmann TO, Schoiswohl G, Ni J, Zechner R, Ni S, Hao H . Enhanced monoacylglycerol lipolysis by ABHD6 promotes NSCLC pathogenesis. *EBioMedicine.* 2020 Mar 53:102696.
5. Beibei Guo, Mengwei Qi, Shuai Huang, Run Zhuo, Wenxue Zhang, Yufang Zhang, Man Xu, Mei Liu, Tuchen Guan, Yan Liu . Cadherin-12 Regulates Neurite Outgrowth Through the PKA/Rac1/Cdc42 Pathway in Cortical Neurons *Front Cell Dev Biol.* 2021 Nov 8;9:768970.
6. Zhou Y, Sun P, Wang T, Chen K, Zhu W, Wang H. . Inhibition of Calcium Influx Reduces Dysfunction and Apoptosis in Lipotoxic Pancreatic β -Cells via Regulation of Endoplasmic Reticulum Stress. *pLos one.* 2015 Jul 6;10(7):e0132411.
7. Sun WC, Liang ZD, Pei L. . Propofol-induced rno-miR-665 targets BCL2L1 and influences apoptosis in rodent developing hippocampal astrocytes. *Neurotoxicology.* 2015 Aug 5;51:87-95.
8. Sun W, Pei L. . microRNA Expression Profiling of Propofol-Treated Developing Rat Hippocampal Astrocytes. *DNA Cell Biol.* 2015 Aug;34(8):511-23.
9. Wang X, Ma J, Fu Q, Zhu L, Zhang Z, Zhang F, Lu N, Chen A. . Role of hypoxia inducible factor 1 α in autophagic cell death in microglial cells induced by hypoxia. *Mol Med Rep.* 2017 Apr;15(4):2097-2105.
10. Zhu J, Chen Z, Tian J, Meng Z, Ju M, Wu G, Tian Z. . miR-34b attenuates trauma-induced anxiety-like behavior by targeting CRHR1. *Int J Mol Med.* 2017 Jul;40(1):90-100.
11. Nie H, Wang A, He Q, Yang Q, Liu L, Zhang G, Huang Y, Ding X, Yu H, Hu S. . Phenotypic switch in lung interstitial macrophage polarization in an ovalbumin-induced mouse model of asthma. *Exp Ther Med.* 2017 Aug;14(2):1284-1292.
12. Zhao HJ, Zhao XH . Modulatory Effect of the Supplemented Copper Ion on In Vitro Activity of Bovine Lactoferrin to Murine Splenocytes and RAW264.7 Macrophages. *Biol Trace Elem Res.* 2018 Aug 16.
13. Yang T, Cui H, Wen M, Zuber J, Kogan SC, Wei G . TCEA1 regulates the proliferative potential of mouse myeloid cells. *Exp Cell Res.* 2018 Sep 15;370(2):551-560.

14. Gong H, Sheng X, Xue J, Zhu D . MicroRNA-365 regulates the occurrence and immune response of sepsis following multiple trauma via interleukin-6. *Exp Ther Med*. 2018 Oct;16(4):3745-3751.
15. Zhang Y, Sun R, Geng S, Shan Y, Li X, Fang W . Porcine Circovirus Type 2 Induces ORF3-Independent Mitochondrial Apoptosis via PERK Activation and Elevation of Cytosolic Calcium. *J Virol*. 2019 Mar 21 93(7). pii: e01784-18.
16. Wu FB, Yang YY, Wang XB, Wang Z, Zhang WW, Liu ZY, Qian YQ . A sample processing method for immunoassay of whole blood tacrolimus. *Anal Biochem*. 2019 Jul 1 576:13-19.
17. Yan X, Wang J, Zhu Y, Feng W, Zhai C, Liu L, Shi W, Wang Q, Zhang Q, Chai L, Li M . S1P induces pulmonary artery smooth muscle cell proliferation by activating calcineurin/NFAT/OPN signaling pathway. *BIOCHEM BIOPH RES CO*. 2019 Aug 27 516(3):921-927.
18. Cao T, Li X, Xu Y, Zhang S, Wang Z, Shan Y, Sun J, Fang W, Li X . Npro of Classical Swine Fever Virus Suppresses Type III Interferon Production by Inhibiting IRF1 Expression and Its Nuclear Translocation. *VIRUSES-BASEL*. 2019 Oct 31 11(11).
19. Zhang J, Xu HY, Wu YJ, Zhang X, Zhang LQ, Li YM . Neutrophil elastase inhibitory effects of pentacyclic triterpenoids from *Eriobotrya japonica* (loquat leaves). *J Ethnopharmacol*. 2019 Oct 5 242:111713.
20. Pan Y, Li W, Feng Y, Xu J, Cao H . Edaravone attenuates experimental asthma in mice through induction of HO-1 and the Keap1/Nrf2 pathway. *Exp Ther Med*. 2020 Feb 19(2):1407-1416.
21. Hui Zou, Tao Wang, Junzhao Yuan, Jian Sun, Yan Yuan, Jianhong Gu, Xuezhong Liu, Jianchun Bian, Zongping Liu . Cadmium-induced cytotoxicity in mouse liver cells is associated with the disruption of autophagic flux via inhibiting the fusion of autophagosomes and lysosomes *Toxicol Lett*. 2020 Mar 15;321:32-43.
22. Chengbao Li, Xuli Rong, Jing Qin, Sheng Wang . An improved two-step method for extraction and purification of primary cardiomyocytes from neonatal mice *J Pharmacol Toxicol Methods*. Jul-Aug 2020;104:106887.
23. Shuai Yuan, Chuanxin Zhang, Yunli Zhu, Bo Wang . Neohesperidin Ameliorates Steroid-Induced Osteonecrosis of the Femoral Head by Inhibiting the Histone Modification of lncRNA HOTAIR *DRUG DES DEV THER*. 2020 Dec 7;14:5419-5430.
24. Zhixin He, Yajing Zhang, Huijie Zhang, Chao Zhou, Qinlong Ma, Ping Deng, Muxue Lu, Zhenlin Mou, Min Lin, Lingling Yang, Yanqi Li, Yang Yue, Hui Feng Pi, Yonghui Lu, Mindi He, Lei Zhang, Chunhai Chen, Zhou Zhou, Zhengping Yu . NAC antagonizes arsenic-induced neurotoxicity through TMEM179 by inhibiting oxidative stress in Oli-neu cells *Ecotoxicol Environ Saf*. 2021 Oct 15;223:112554.
25. Chuanxin Zhang, Shuai Yuan, Yi Chen, Bo Wang . Neohesperidin promotes the osteogenic differentiation of human bone marrow stromal cells by inhibiting the histone modifications of lncRNA SNHG1 *Cell Cycle*. 2021 Oct;20(19):1953-1966.
26. Xiaojing Liang, Jiansheng Xie, Hao Liu, Rongjie Zhao, Wei Zhang, Haidong Wang, Hongming Pan, Yubin Zhou, Weidong Han . STIM1 Deficiency In Intestinal Epithelium Attenuates Colonic Inflammation and Tumorigenesis by Reducing ER Stress of Goblet Cells *Cell Mol Gastroenterol Hepatol*. 2022;14(1):193-217.
27. Qiang Yang, Feiyu Long . MiRNA-384-5p targets GABRB1 to regulate ketamine-Induced Neurotoxicity in Neurons *Turk Neurosurg*. 2022 Jan 20. doi: 10.5137/1019-5149.JTN.36367-21.2. (IF 0.963)
28. Meng Duan, Hainan Chen, Linjie Yin, Xiao Zhu, Petr Novák, Yuncheng Lv, Guojun Zhao, Kai Yin . Mitochondrial apolipoprotein A-I binding protein alleviates atherosclerosis by regulating mitophagy and macrophage polarization *Cell Commun Signal*. 2022 May 7;20(1):60.
29. Qian Li, Liuqian Wang, Di Ji, Wei Yu, Yan Zhang, Yanghong Xiang, Chao Zhou, Liting Wang, Ping Deng, Hui Feng Pi, Yonghui Lu, Qinlong Ma, Mindi He, Lei Zhang, Zhengping Yu, Anchun Deng . Metformin attenuates cadmium-induced degeneration of spiral ganglion neuron via restoring autophagic flux in primary culture *J Inorg Biochem*. 2022 Sep;234:111901.

Version 2024.03.12