



碧云天生物技术/Beyotime Biotechnology  
 订货热线: 400-1683301 或 800-8283301  
 订货 e-mail: order@beyotime.com  
 技术咨询: info@beyotime.com  
 网址: http://www.beyotime.com

## NTP set (100mM each, Nuclease free)

产品编号	产品名称	包装
D7383-1ml	NTP set (100mM each, Nuclease free)	4×250μl

### 产品简介:

- 本NTP set (100mM each, Nuclease free)为ATP、CTP、GTP和UTP四种溶液的套装，每种浓度均为100mM，并经测试不含Nuclease，即不含DNase和RNase。
- 本产品适用于常规的分子生物学和生物化学用途，例如用于SP6、T3和T7 RNA聚合酶等进行体外转录以产生RNA，或用于体外转录RNA，以及RNA的合成与扩增。
- 本NTP set溶液是用Nuclease-free水配制，并含有1mM Tris, pH7.3-7.5浓度。
- 本NTP set不含DNase、RNase、phosphatase和蛋白酶。

### 包装清单:

产品编号	产品名称	包装
D7383-1	ATP (100mM, Nuclease free)	250μl
D7383-2	CTP (100mM, Nuclease free)	250μl
D7383-3	GTP (100mM, Nuclease free)	250μl
D7383-4	UTP (100mM, Nuclease free)	250μl
—	说明书	1份

### 保存条件:

-20°C保存，2年有效。

### 注意事项:

- 可以室温溶解，溶解后宜存放于冰盒内或冰浴上，使用完毕后宜立即置于-20°C保存。
- 本产品仅限于专业人员的科学研究用，不得用于临床诊断或治疗，不得用于食品或药品，不得存放于普通住宅内。
- 为了您的安全和健康，请穿实验服并戴一次性手套操作。

### 相关产品:

产品编号	产品名称	包装
D7351	dATP(100mM)	250μl
D7353	dCTP(100mM)	250μl
D7355	dGTP(100mM)	250μl
D7357	dTTP(100mM)	250μl
D7359	dUTP(100mM)	250μl
D7366	dNTP (4管套装, 100mM)	4×50μl
D7368	dNTP (5管套装, 100mM)	5×50μl
D7371	dNTP Mixture (2.5mM each)	1ml
D7373	dNTP Mixture (25mM each)	250μl
D7376-1ml	dNTP/dUTP Mixture (2.5mM each/5mM)	1ml
D7376-5ml	dNTP/dUTP Mixture (2.5mM each/5mM)	5ml
D7378-250μl	ATP (100mM, Nuclease free)	250μl
D7378-1ml	ATP (100mM, Nuclease free)	1ml
D7378-5ml	ATP (100mM, Nuclease free)	5×1ml
D7379-250μl	CTP (100mM, Nuclease free)	250μl
D7379-1ml	CTP (100mM, Nuclease free)	1ml
D7379-5ml	CTP (100mM, Nuclease free)	5×1ml
D7380-250μl	GTP (100mM, Nuclease free)	250μl

D7380-1ml	GTP (100mM, Nuclease free)	1ml
D7380-5ml	GTP (100mM, Nuclease free)	5×1ml
D7381-250µl	UTP (100mM, Nuclease free)	250µl
D7381-1ml	UTP (100mM, Nuclease free)	1ml
D7381-5ml	UTP (100mM, Nuclease free)	5×1ml
D7383-1ml	NTP set (100mM each, Nuclease free)	4×250µl
D7385-500µl	NTP Mix (10mM each, Nuclease free)	500µl
D7385-2ml	NTP Mix (10mM each, Nuclease free)	2ml
D7387-250µl	NTP Mix (25mM each, Nuclease free)	250µl
D7387-1ml	NTP Mix (25mM each, Nuclease free)	1ml

Version 2021.05.31