



碧云天生物技术/Beyotime Biotechnology  
订货热线: 400-1683301或800-8283301  
订货e-mail: order@beyotime.com  
技术咨询: info@beyotime.com  
网址: http://www.beyotime.com

## 磷酸氢二钠(99%, Reagent grade)

产品编号	产品名称	包装
ST1655-100g	磷酸氢二钠(99%, Reagent grade)	100g
ST1655-500g	磷酸氢二钠(99%, Reagent grade)	500g
ST1655-6×500g	磷酸氢二钠(99%, Reagent grade)	6×500g

### 产品简介:

CAS Number	Chemical Formula	Molecular Weight	Purity	Grade
7558-79-4	Na <sub>2</sub> HPO <sub>4</sub>	141.96	99%	Reagent grade

### ➤ 基本信息(General Information):

Name (Chinese)	磷酸氢二钠
Name (English)	Sodium phosphate dibasic
Specifications	Reagent grade, 99%
Chemical Formula	Na <sub>2</sub> HPO <sub>4</sub>
Synonym (Chinese)	磷酸一氢钠, 磷酸氢二钠
Synonym (English)	sec-Sodium phosphate, Disodium hydrogen phosphate, Disodium phosphate, Sodium hydrogenphosphate
Beilstein Registry No.	-
EINECS Number	231-448-7
MDL Number	MFCD00003496
UNSPSC Code	-



### ➤ 产品描述(Description):

Application	Sodium phosphate dibasic has been used in the establishment of epithelial cell lines from individual rabbit lenses.
-------------	---

### ➤ 性质(Properties):

grade	reagent grade
assay	99%

### ➤ 安全信息(Safety Information):

Hazard Pictogram Codes	-
Signal Word	-
Hazard Statements	-
Precautionary Statements	-
Personal Protective Equipment	-
Hazard Codes (Europe)	-
Risk Codes (Europe)	-
Safety Codes (Europe)	-
RIDADR	NONH for all modes of transport
WGK Germany	1
RTECS	-
Flash Point (F)	-
Flash Point (C)	-

包装清单:

产品编号	产品名称	包装
ST1655-100g	磷酸氢二钠(99%, Reagent grade)	100g
ST1655-500g	磷酸氢二钠(99%, Reagent grade)	500g
ST1655-6×500g	磷酸氢二钠(99%, Reagent grade)	6×500g
—	说明书	1份

保存条件:

室温保存, 三年有效。

注意事项:

- 本产品仅限于专业人员的科学研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

Version 2021.12.14