



碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301或800-8283301
 订货e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

Triton X-100/曲拉通X-100 (Reagent grade)


产品编号	产品名称	包装
ST1723-100ml	Triton X-100/曲拉通X-100 (Reagent grade)	100ml
ST1723-500ml	Triton X-100/曲拉通X-100 (Reagent grade)	500ml
ST1723-5L	Triton X-100/曲拉通X-100 (Reagent grade)	5L
ST1723-20L	Triton X-100/曲拉通X-100 (Reagent grade)	20L

产品简介:

CAS Number	Chemical Formula	Molecular Weight	Purity	Grade
9036-19-5/9002-93-1	-	-	-	Reagent grade

➤ 基本信息(General Information):

Name (Chinese)	曲拉通X-100
Name (English)	Triton X-100
Specifications	Reagent grade
Chemical Formula	-
Synonym (Chinese)	4-(1,1,3,3-四甲基丁基)苯基-聚乙二醇, 叔辛基苯氧基聚乙烯氧基乙醇, 聚乙二醇叔辛基苯基醚
Synonym (English)	4-(1,1,3,3-Tetramethylbutyl)phenyl-polyethylene glycol; t-Octylphenoxy polyethoxyethanol; Polyethylene glycol tert-octylphenyl ether
Beilstein Registry No.	-
EINECS Number	-
MDL Number	MFCD00128254
UNSPSC Code	12161900



➤ 产品描述(Description):

Biochem/physiol Actions	广泛使用的非离子表面活性剂, 用于在温和的非变性条件下进行膜恢复。
-------------------------	-----------------------------------

➤ 性质(Properties):

HLB	13.5
CMC	0.2-0.9 mM(20-25°C)
aggregation number	100-155
mol wt	average mol wt 625
transition temp	cloud point 67°C
transition temp	pour point ~7°C
grade	laboratory grade
form	liquid
mol wt	micellar avg mol wt 80,000
description	non-ionic
solubility	water: soluble 0.1mL/mL, clear to slightly hazy, colorless to faintly yellow

➤ 安全信息(Safety Information):

Hazard Pictogram Codes	
Signal Word	Danger
Hazard Statements	H302-H315-H318-H410

Precautionary Statements	P280-P301 + P312 + P330-P305 + P351 + P338 + P310
Personal Protective Equipment	Eyeshields, Faceshields, full-face respirator (US), Gloves, multi-purpose combination respirator cartridge (US), type ABEK (EN14387) respirator filter, type N95 (US), type P1 (EN143) respirator filter
Hazard Codes (Europe)	Xn,N
Risk Codes (Europe)	22-36-51/53
Safety Codes (Europe)	26-39-61
RIDADR	UN 3082 9 / PGIII
WGK Germany	2
RTECS	MD0907700
Flash Point (F)	483.8 F
Flash Point (C)	251°C

包装清单:

产品编号	产品名称	包装
ST1723-100ml	Triton X-100/曲拉通X-100 (Reagent grade)	100ml
ST1723-500ml	Triton X-100/曲拉通X-100 (Reagent grade)	500ml
ST1723-5L	Triton X-100/曲拉通X-100 (Reagent grade)	5L
ST1723-20L	Triton X-100/曲拉通X-100 (Reagent grade)	20L
—	说明书	1份

保存条件:

室温保存。

注意事项:

- 本产品仅限于专业人员的科学研究用，不得用于临床诊断或治疗，不得用于食品或药品，不得存放于普通住宅内。
- 为了您的安全和健康，请穿实验服并戴一次性手套操作。

Version 2024.12.23